



Customer Release Notes

eRAD PACS
Version 8.0
Viewer Build 301
Release Date February 8, 2024

Table of Contents

New Features and Enhancements.....	3
1. Separatespine label state for each study	3
2. Slow referenced image searching	3
Resolved Issues.....	3
1. Deleting multiple redirected studies	3
2. Encapsulated pdf documents are not uploaded with CDUpload	3

New Features and Enhancements

1. Separate spine label state for each study

A spine label state exists for each study so when applying them, they advance sequentially even if the user interrupts labeling one study and starts labeling a second. Note that spine labels apply to studies as defined by the study UID, not the frame of reference, meaning studies whose frames of reference have been linked, e.g., through a spatial registration object, do not share spine labels. Each study must be individually labeled.

2. Slow referenced image searching

The algorithm used to find a referenced image object, such as is used in a presentation state object, has been modified to improve overall performance.

Resolved Issues

1. Deleting multiple redirected studies

When redirecting across three levels, such as launching a study from a consolidation worklist server and downloading data from its grandchild, cache entries are created to handle the redirection. When the cache is purged to make space for new data, some of the unused intermediate cache entries do not get deleted properly, causing misinformation in the cache manager.

2. Encapsulated pdf documents are not uploaded with CDUpload

When uploading data from a DICOM CD containing a DICOMDIR file, encapsulated PDF objects were skipped and not uploaded because the upload application didn't check for these object types.