



# **Customer Release Notes**

eRAD PACS  
Version 8.0  
Server Build 34  
Release Date July 16, 2021

## Table of Contents

<b>Resolved Issues.....</b>	<b>3</b>
<b>1. Utilize a Date Range filter for "last seen" field on Users page.....</b>	<b>3</b>
<b>2. Session table signoff column timeout entry uses the wrong date format.....</b>	<b>3</b>
<b>3. Study edit updates Receive Date.....</b>	<b>3</b>

## Resolved Issues

**1. Utilize a Date Range filter for "last seen" field on Users page**

A free text filter panel on the user accounts page was used for session date fields, including Last Seen, Sing On and Sign Off.

**2. Session table signoff column timeout entry uses the wrong date format**

The timed-out sign off event's date was saved in the database using a date format inconsistent with other saved dates. This has been resolved.

**3. Study edit updates Receive Date**

Editing a study with soft edit enabled always updated the Receive Date value because it stopped checking if the modified object was the MinSOPInstance as a result of a change in medsrv-8.0.32. This has been resolved.