



# **Customer Release Notes**

eRAD PACS  
Version 8.0  
Server Build 51  
Release Date April 27, 2022

## Table of Contents

<b>Resolved Issues.....</b>	<b>3</b>
<b>1. Tasks filter window allows user to go to another web page without closing.....</b>	<b>3</b>
<b>2. Unable to forward from patient folder .....</b>	<b>3</b>
<b>3. Notify action error message .....</b>	<b>3</b>
<b>4. User with empty password cannot be edited.....</b>	<b>3</b>

## Resolved Issues

**1. Tasks filter window allows user to go to another web page without closing**

If the tasks page's filter panel was open when the user called up a different web page, the filter panel was not closed and remained on the screen. This has been resolved.

**2. Unable to forward from patient folder**

Back-end support for forwarding studies from the patient folder was incomplete, resulting in no action when the button was clicked. This has been resolved.

**3. Notify action error message**

Identifying the email address offered as the default when configuring a Notify action failed for internally-defined accounts, such as the system account. In this case, the default comes up empty, requiring the user to explicitly declare the email recipient. This has been resolved.

**4. User with empty password cannot be edited**

User accounts with empty passwords, which is not a valid state, could not be corrected because the missing password was not handled and caused an error when editing. This has been resolved.