



# Customer Release Notes

eRAD PACS  
Version 8.0  
Viewer Build 62  
Release Date June 3, 2021

## Table of Contents

<b>Resolved Issues.....</b>	<b>3</b>
<b>1. Viewer crash at exit.....</b>	<b>3</b>
<b>2. Progression view can be wrong state.....</b>	<b>3</b>

## Resolved Issues

### 1. Viewer crash at exit

The viewer could crash when closing if the structure identifying a concurrent session was NULL. This has been resolved.

### 2. Progression view can be wrong state

When starting a progression with the previous view tool, progressing with the next view and then again going back using the previous view, the progression state was mishandled and the previous progression view was not displayed. This has been resolved.