



Customer Release Notes

eRAD PACS
Version 8.0
Viewer Build 61

Release Date May 20, 2021

Table of Contents

Resolved Issues.....	3
1. Viewer crash when a close command is called while creating the Patient Folder panel.....	3
2. Viewer crashed when the report was processed after the close	3
3. Viewer crash when show Patient Folder command was not finished	3

Resolved Issues

1. Viewer crash when a close command is called while creating the Patient Folder panel

If the web viewer receives a command to close the session while it is still initializing the patient folder panel, it might crash. This has been resolved.

2. Viewer crashed when the report was processed after the close

If a report is received after the viewer starts terminating the session, the viewer might crash. This has been resolved.

3. Viewer crash when show Patient Folder command was not finished

If a second web control request to open the patient folder panel is issued before the first completes, the viewer could crash. This has been resolved.