



Customer Release Notes

eRAD PACS
Version 8.0
Viewer Build 303
Release Date April 4, 2024

Table of Contents

New Features and Enhancements.....	3
1. Don't drop the context menu request if the right mouse button is down on an annotation....	3
Resolved Issues.....	3
1. Incorrect response to WMI command "addpartner"	3
2. Finding image references doesn't work in all cases.....	3

New Features and Enhancements

1. **Don't drop the context menu request if the right mouse button is down on an annotation**

If the user moved the mouse while the context (right-click) menu was displayed, the menu disappeared. Now, the menu remains visible and active until the user clicks on it or clears it.

Resolved Issues

1. **Incorrect response to WMI command "addpartner"**

The response message to WMI launch commands using the addpartner option included the required tag but left out the status value.

2. **Finding image references doesn't work in all cases**

If a multi-frame mammography (tomo) object defines a frame number "0", the CAD header would appear on the image only if the image contained findings.