



Customer Release Notes

eRAD PACS
Version 8.0
Viewer Build 63
Release Date July 8, 2021

Table of Contents

New Features and Enhancements.....	3
1. Auto select first image frame	3
2. Reduce number of points in free form ROI annotation	3
Resolved Issues.....	3
1. CDUpload write error details in log	3

New Features and Enhancements

1. Auto select first image frame

To facilitate the activation of keyboard-driven commands, the viewer automatically selects the first image frame after a hanging protocol is applied or a new view page is selected.

2. Reduce number of points in free form ROI annotation

The number of points stored in free form ROI annotations has been consolidated, resulting in less annotation data and allowing for more annotations saved per key image.

Resolved Issues

1. CDUpload write error details in log

Details about specific SSL errors are saved in the log file to aid in debug investigations. These details have been removed from the message box displayed on the user interface and replaced with the specific error code.