

Server Update

**For eRAD RIS
Version 1.0
Build 37
Update 37.013**

Table of Contents

1. Purpose	3
Defects resolved	3
2. Package Contents	4
4. Intended Audience	5
5. Who is affected.....	5
6. eRAD RIS Release Version Numbers	5
7. Installing	6

1. Purpose

This server update has been produced to fix bugs required for production setting, or to complete an enhancement requested for production.

This server update includes the following:

Defects resolved

Redmine #	Subject	Priority
1442	Save IVT with large Precert note generates String or binary data would be truncated error	Normal
2107	error on opening registration when patient has multiple orders for today	Normal
2122	selecting search result times/typing date/times in manually could cause the wrong duration override to be used	Normal
2098	Registration - When registering an exam for a patient with multiple exams for the same day	High
2081	VR2 - Messages from Merge with alternate patient ids no longer Merging	Urgent

2. Package Contents

Files included for this server update:

1. GUI ZIP Updates folder contains
 - a. CommonLib.dll
 - b. CommonLib.pdb
 - c. Ris.exe
 - d. Ris.pdb
 - e. UILib.dll
 - f. UILib.pdb
 - g. WCFClient.exe
 - h. WCFClient.pdb

2. RIS Service Updates folder contains
 - a. CDSJobManager.dll
 - b. CDSJobManager.pdb
 - c. CommonLib.dll
 - d. bommonLib.pdb
 - e. rRISServices.exe
 - f. rRISServices.pdb
 - g. ScheduledTaskManager.dll
 - h. ScheduledTaskManager.pdb
 - i. SignatureApprovalService.exe
 - j. SignatureApprovalService.pdb
 - k. WCFSservice.dll
 - l. WCFSservice.pdb
 - m. Wedge.exe
 - n. Wedge.pdb

3. SQL folder contains:
 - a. B37.013_updates.sql

4. Intended Audience

This document is created by the RIS Development team for the RadNet RIS Implementation team.

5. Who is affected

Previous versions of build 37.

Server updates may not be cumulative; therefore this server update must be applied directly to the previous server update.

6. eRAD RIS Release Version Numbers

Build	Patch	UI Version	Core Version	WS Version	DB Version	Notes
37	-	1.2.37.12843	1.2.37.12782	1.2.37.12843	1.2.37.12843	
37	1	1.2.37.13040	1.2.37.13040	1.2.37.13039	1.2.37.13003	
37	2	1.2.37.13069	1.2.37.13040	1.2.37.13069	1.2.37.13003	
37	3	1.2.37.13069	1.2.37.13040	1.2.37.13069	1.2.37.13003	Just WebAPI updated, no version changes
37	4	1.2.37.13171	1.2.37.13040	1.2.37.13165	1.2.37.13171	
						Bug was discovered with version numbers, this is why it appears we went backwards in the version number but we really did not.
37	5	1.2.37.13347	1.2.37.13040	1.2.37.12979	1.2.37.12979	
37	6	1.2.37.13665	1.2.37.13040	1.2.37.13580	1.2.37.12979	
37	7	1.2.37.13807	1.2.37.13040	1.2.37.13799	1.2.37.13807	
37	8	1.2.37.13807	1.2.37.13040	1.2.37.13827	1.2.37.13807	
37	9	1.2.37.13807	1.2.37.13040	1.2.37.13839	1.2.37.13807	
37	10	1.2.37.13895	1.2.37.13040	1.2.37.13895	1.2.37.13895	
37	11	1.2.37.14161	1.2.37.13040	1.2.37.14161	1.2.37.14161	
37	12	1.2.37.14285	1.2.37.13040	1.2.37.14300	1.2.37.14285	
37	13	1.2.37.14326	1.2.37.13040	1.2.37.14326	1.2.37.14326	

7. Installing

- 1) Copy the files in the GUI ZIP Updates folder and replace those in the rRIS_1.2.37.14285.zip folder. Rename the zip file to be **rRIS_1.2.37.14326.zip**
- 2) Stop each service listed below. Copy the files from the RIS Service Updates folder to the directory for each rRISService.
rRISService directories include:
 - a) *rRISService*
 - b) *CDS rRISService*
 - c) *Wedge rRISService*
 - d) *DDE rRISService*

Restart each service.
- 3) Run the B37.013_updates.sql file against the eRAD RIS database instance.