

Server Update

for eRAD RIS

Version 1.0

Build 37

Update 37.009

Table of Contents

1. Purpose.....	3
Defects resolved.....	3
2. Package Contents	4
3. Intended Audience	4
4. Who is affected	4
5. Installing	4

1. Purpose

This server update has been produced to fix bugs required for production setting, or to complete an enhancement requested for production.

This server update includes the following:

Defects resolved

Redmine #	Subject	Priority
1863	Signed exams fail billing – “No active billing codes for this study”	High
1866	Action data does not have all studies in the study container XML	Med-High

2. Package Contents

Files included for this server update:

1. RIS Service Update folder contains:
 - a. WCFSservice.dll
 - b. WCFSservice.pdb

3. Intended Audience

This document is created by the RIS Development team for the RadNet RIS Implementation team.

4. Who is affected

Previous versions of build 37.

Server updates may not be cumulative; therefore this server update must be applied directly to the previous server update.

5. Installing

- 1) Stop each service listed below. Copy the files from the RIS Service Updates folder to the directory for each rRISService.

rRISService directories include:

- a) *rRISService*
- b) *CDS rRISService*
- c) *Wedge rRISService*
- d) *DDE rRISService*

Restart each service