

# Server Update

For eRAD RIS  
Version 1.0  
Build 40  
Update 40.022.1

## **Table of Contents**


1. Defects Resolved .....	3
2. Code Stream.....	4
3. Package Contents .....	5
4. Intended Audience .....	5
5. Who is affected .....	5
6. eRAD RIS Release Version Numbers .....	6
7. Installing .....	7

### **1. Defects Resolved**


Redmine #	Subject	Priority
6190	Hidden text in editor counts as text when validating fields and section	High

The resolved bug listed above resolves an issue in which users had the option to hide text within the reporting window. This option was available to radiologist and transcription users. The workflow for this option has been altered. You can no longer hide text, only unhide. If there is hidden text in a report, transcription users cannot transcribe and radiologist users cannot sign. They will be presented with a message like the one displayed below.

Message presented to transcription




Message presented to radiologist



## 2. Code Stream

### eRAD RIS - Source Code Evolution



#### **Legend:**

**Light Green = Released software**

**Gray = Will not be released**

**Bright Green = Current Release**

### **3. Package Contents**

Files included for this server update:

- 1) GUI ZIP Updates folder contains
  - a) rRIS\_1.3.40.21936(3GB).zip

### **4. Intended Audience**

This document is created by the RIS Development team for the RadNet RIS Implementation team.

### **5. Who is affected**

Build 40 installs. This server update must be applied directly to the build 40.020.

## RADNET, Inc. – eRAD RIS Server Update

### 6. eRAD RIS Release Version Numbers

Build	Patch	UI Version	Core Version	WS Version	DB Version	Notes
40	16	1.1.3.40.18585 (3GB)	1.3.40.16659	1.3.40.18585	1.3.40.18646	GUI.zip, Web service, database updates
40	16.1	1.3.40.18705 (3GB)	1.3.40.16659	1.3.40.18585	1.3.40.18646	GUI.zip
40	17	1.3.40.18780 (3GB)	1.3.40.18756	1.3.40.18756	1.3.40.18646	<b>Core install required</b> for Radiologist User Group. GUI.zip, Core setup, and Web Service updates
40	18	1.3.40.19253(3GB)	1.3.40.18756	1.3.40.18756	1.3.40.19253	GUI.zip, Web Service and database updates
40	19	1.3.40.19909(3GB)	1.3.40.18756	1.3.40.19878	1.3.40.19253	GUI.zip, Web Service
40	20	1.3.40.20086(3GB)	1.3.40.18756	1.3.40.19878	1.3.40.19253	Gui.zip and Recognizer
40	21	1.3.40.20703(3GB)	1.3.40.18756	1.3.40.19878	1.3.40.19253	Gui.zip
40	22	1.3.40.20898(3GB)	1.3.40.18756	1.3.40.19878	1.3.40.19253	Gui.zip
<b>40</b>	<b>22.1</b>	<b>1.3.40.21936(3GB)</b>	<b>1.3.40.18756</b>	<b>1.3.40.19878</b>	<b>1.3.40.19253</b>	<b>Gui.zip</b>
41	-	1.3.41.19135(3GB)	1.3.41.19135	1.3.41.19135	1.3.41.19135	Full version release
41	1	1.3.41.19356(3GB)	1.3.41.19135	1.3.41.19346	1.3.41.19346	GUI.zip, Web Service and database updates
41	2	1.3.41.19620(3GB)	1.3.41.19135	1.3.41.19620	1.3.41.19346	GUI.zip and Web Service
41	2sf	1.3.41.19620SF(3G)	1.3.41.19135	1.3.41.19620	1.3.41.19346	GUI.zip
41	2.1	1.3.41.20622(3GB)	1.3.41.19135	1.3.41.20617	1.3.41.19346	GUI.zip and Web Service
41	3	1.3.41.20171(3GB)	1.3.41.19135	1.3.41.20171	1.3.41.19346	GUI.zip , Web Service and Recognizer
41	4	1.3.41.20285(3GB)	1.3.41.19135	1.3.41.20171	1.3.41.19346	GUI.zip
41	5	1.3.41.20809(3GB)	1.3.41.20708	1.3.41.20792	1.3.41.20738	GUI.zip, Core install, Web Service and database updates
41	6	1.3.41.20912(3GB)	1.3.41.20708	1.3.41.20792	1.3.41.20738	GUI.zip
41	7	1.3.41.21963(3GB)	1.3.41.20708	1.3.41.20792	1.3.41.20738	GUI.zip
41.1	-	1.3.41.20918(3GB)	1.3.41.20791	1.3.41.20918	1.3.41.20739	GUI.zip, Core install, Web Service and DB updates, Help files
41.1	1	1.3.41.20987(3GB)	1.3.41.20791	1.3.41.20918	1.3.41.20739	GUI.zip
41.1	2	1.3.41.21162(3GB)	1.3.41.20791	1.3.41.21162	1.3.41.21150	GUI.zip, Web Service and DB updates,
41.1	3	1.3.41.21264(3GB)	1.3.41.20791	1.3.41.21162	1.3.41.21150	GUI.zip
41.1	4	1.3.41.21931(3GB)	1.3.41.20791	1.3.41.21162	1.3.41.21150	GUI.zip
42	1	2.42.0.0.21577	2.42.0.0.21261	2.42.0.0.21577	2.42.0.0.21187	GUI and Web Service updates
42	2	2.42.0.1.21702	2.42.0.0.21261	2.42.0.1.21702	2.42.0.0.21187	GUI and Web Service updates
42	3	2.42.0.2.21806	2.42.0.0.21261	2.42.0.2.21806	2.42.0.0.21187	GUI and Web Service updates
42	4	2.42.0.4.21862	2.42.0.0.21261	2.42.0.2.21806	2.42.0.4.21870	GUI and Database updates
42	5	2.42.0.5.21924	2.42.0.0.21261	2.42.0.2.21806	2.42.0.4.21870	GUI

## **7. Installing**

- 1) Copy and replace the current rRIS...zip file with the rRIS\_1.3.40. 21936.zip file provided with this release.  
Be sure to make a backup of the ris.exe.config file contained within the current .zip folder.

**NOTE:** If your upgrade is coming from 40.016.1 or prior a core install will be required. Versions (builds) of 40.017 and greater do not require a core install as long as you are at least at 40.017 (Core Version 1.3.40.18756)