

HOTFIX

for eRAD RIS
Version 1.0
Build 35
Hotfix 35.009

Table of Contents

1. Purpose.....	3
2. Intended Audience	3
3. Who is affected	3
4. Installing	3

1. Purpose

This HOTFIX has been produced to fix bugs required for production setting, or to complete an enhancement requested for production.

This HOTFIX includes the following:

- 1) Outbound HL7 message queue was blocked due to person key resolution problem.
 - a) Wedge.exe, Wedge.pdb

2. Intended Audience

This document is created by the RIS Development team for the RadNet RIS Implementation team.

3. Who is affected

All previous version of build 35. Current build is 1.1.35.10752.

4. Installing

- 1) *Copy the 2 files to the directory for each rRISService. The changes to these files are only relevant to the Wedge, but it is recommended to also copy them to the other rRISService directories.*