



eRAD RIS

RELEASE ANNOUNCEMENT

Hotfix 4.2023.057.2

UPDATED AUGUST 12, 2025

TABLE OF CONTENTS

Summary	2
Release Announcement.....	2
New Features.....	2
Feature Details	3
Billing.....	3
Feature #37422 – Imagine Simplicity Integration.....	3
Version Details	4
Code Stream.....	4

PUBLICATION HISTORY

Revision	Author	Description
August 12, 2025	Michelle Mahabir/ Kevin Brooks	▪ Commercial release.

SUMMARY

Release Announcement

This hotfix of ERAD RIS 4.2023.057.2 enhances the ability of RIS to integrate with billing systems such as Imagine Simplicity.

New Features

This release introduces the following features and enhancements:

Category	Redmine #	Subject	Description
Billing	37422	Imagine Simplicity Integration	This enhancement to Billing introduces a service tool that enables the ability to submit requests through Imagine Simplicity, offering a streamlined alternative to Availity and other eligibility services.

Refer to the FEATURE DETAILS section below for configuration and usage information.

FEATURE DETAILS

Billing

Feature #37422 – Imagine Simplicity Integration

Summary

This enhancement to Billing introduces a service tool that enables the ability to submit requests through Imagine Simplicity, offering a streamlined alternative to Availity and other eligibility services.

Background

As more customers connect to a variety of insurance systems, the eligibility process needed greater flexibility. This update introduces a way to route eligibility checks based on payer, helping reduce integration complexity and improve support for different insurance platforms.

Feature Description

The new `Ansi270GenericWS` Web Service allows customers to customize how eligibility checks are handled, based on the insurance provider. It supports integration with systems like `Availity` and `Simplicity`, and automatically routes requests using the Payer ID. This makes the eligibility process more flexible and easier to manage when working with multiple insurance networks.

Service Team Upgrade Instructions

The Service Team must complete the following actions to deploy this feature:

1. Configure the Ansi270 Web Service

- Open the `Ansi270WebService.exe.config` file.
- Set the **IP address** and **port** where the web service should listen for incoming requests.

2. Update eRAD RIS Core Service Configuration

- Open the `RISServicesx64.exe.config` file in the **eRAD RIS core services** directory.
- Replace the current hostname and port for the `Ansi270` endpoint with the values you configured in Step 1.

3. Apply Imagine-Specific Configuration

- Navigate to the folder where the **Ansi270 Windows Service** is deployed.
- Open the `Ansi270WebService.dll.config` file.
- Update any configuration sections that are specific to **Imagine** (e.g., customer-specific settings or integration parameters).
 - These values should be provided in the Imagine-specific config package or documentation.


Configuration Instructions

No System Administrator actions are necessary to enable this feature.

VERSION DETAILS

Code Stream

The following source code branches have been merged into this release:



WARNING: The functionality of this release has not yet been merged forward into any other release.

Please carefully review the code stream diagram of any future releases before deploying to ensure that the functionality of this release will not be lost when upgrading.