

# **Server Update**

**For eRAD RIS  
Version 1.0  
Build 41  
Update 41.1.002**

## **Table of Contents**

1. Defects Resolved .....	3
2. Code Stream.....	4
3. Package Contents .....	5
4. Intended Audience .....	6
5. Who is affected .....	6
6. eRAD RIS Release Version Numbers .....	7
7. Installing .....	8


## **RADNET, Inc. – eRAD RIS Server Update**

### **1. Defects Resolved**

<b>Redmine #</b>	<b>Subject</b>	<b>Priority</b>
5778	Mammo follow-up orders no longer included in the list of outstanding orders when scheduling	Urgent
5776	Error after login: nullable object must have a value	High
5773	Parameter is not valid error while loading thumbnail images in document viewer	High
5772	Patient Folder - can't sort relevance or relevant flag	High
4978	User is able to schedule to a closed room	High
5775	null reference error when opening administration / management reports	Normal
5774	foreign key error when adding follow-up / contact log entries	Normal
5344	Playback with bouncing ball, preventing word selection with double click.	Normal

## 2. Code Stream

eRad RIS– source code evolution



### **Legend:**

**Green = Released software**

**Red = Will not be released**

**Blue = Future release**

### **3. Package Contents**

Files included for this server update:

- 1) GUI ZIP Updates folder contains
  - a) rRIS\_1.3.41.21162 (3GB).zip
- 2) SQL folder contains:
  - a) B41\_1\_002.sql
- 3) RIS\_Service\_Updates
  - a) ABCpdf.dll
  - b) ABCpdf9-32.dll
  - c) ABCpdf9-64.dll
  - d) CDSJobManager.dll
  - e) CDSJobManager.pdb
  - f) CdsLib.dll
  - g) CdsLib.xml
  - h) CdsLib.XmlSerializers.dll
  - i) CommonLib.dll
  - j) CommonLib.pdb
  - k) DDEngine.dll
  - l) DDEngine.pdb
  - m) Interop.FAXCOMEXLib.dll
  - n) Microsoft.Practices.EnterpriseLibrary.Caching.dll
  - o) Microsoft.Practices.EnterpriseLibrary.Caching.pdb
  - p) Microsoft.Practices.EnterpriseLibrary.Common.dll
  - q) Microsoft.Practices.EnterpriseLibrary.Common.pdb
  - r) Microsoft.Practices.EnterpriseLibrary.Logging.dll
  - s) Microsoft.Practices.EnterpriseLibrary.Logging.pdb
  - t) Microsoft.Practices.ServiceLocation.dll
  - u) Microsoft.Practices.ServiceLocation.pdb
  - v) Microsoft.Practices.Unity.dll
  - w) Microsoft.Practices.Unity.Interception.dll
  - x) Microsoft.Practices.Unity.Interception.pdb
  - y) Microsoft.Practices.Unity.pdb
  - z) Microsoft.ReportViewer.Common.dll
  - aa) Microsoft.ReportViewer.Common.pdb
  - bb) Microsoft.ReportViewer.WinForms.xml
  - cc) Newtonsoft.Json.dll
  - dd) Newtonsoft.Json.xml
  - ee) Ninject.dll
  - ff) Ninject.Web.Common.dll
  - gg) Ninject.Web.Common.xml
  - hh) Ninject.xml
  - ii) O2S.Components.PDFRender4NET.dll
  - jj) PDFToImages.exe
  - kk) rRISServices.exe
  - ll) rRISServices.pdb
  - mm) rRISServicesx64.exe
  - nn) ScheduledTaskManager.dll
  - oo) ScheduledTaskManager.pdb
  - pp) ServicesConfigWizard.exe
  - qq) ServicesConfigWizard.pdb
  - rr) SignatureApprovalService.exe
  - ss) SignatureApprovalService.pdb
  - tt) System.Net.Http.dll
  - uu) System.Net.Http.Formatting.dll
  - vv) System.Net.Http.Formatting.xml
  - ww) System.Net.Http.WebRequest.dll
  - xx) System.Net.Http.WebRequest.xml
  - yy) System.Net.Http.xml
  - zz) System.Web.dll
  - aaa) System.Web.Helpers.dll
  - bbb) System.Web.Helpers.xml
  - ccc) System.Web.Http.dll
  - ddd) System.Web.Http.SelfHost.dll
  - eee) System.Web.Http.SelfHost.xml
  - fff) System.Web.Http.WebHost.dll
  - ggg) System.Web.Http.WebHost.xml
  - hhh) System.Web.Http.xml
  - iii) System.Web.Mvc.dll
  - jjj) System.Web.Mvc.xml
  - kkk) System.Web.Optimization.dll
  - lll) System.Web.Providers.dll
  - mmm) System.Web.Razor.dll
  - nnn) System.Web.Razor.xml

## **RADNET, Inc. – eRAD RIS Server Update**

ooo)	System.Web.WebPages.Deployment.dll	www)	Tourreau.Gilles.FaxDotNet.xml
ppp)	System.Web.WebPages.Deployment.xml	xxx)	WCFSservice.dll
qqq)	System.Web.WebPages.dll	yyy)	WCFSservice.pdb
rrr)	System.Web.WebPages.Razor.dll	zzz)	WebApi.dll
sss)	System.Web.WebPages.Razor.xml	aaaa)	WebApi.pdb
ttt)	System.Web.WebPages.xml	bbbb)	Wedge.exe
uuu)	Tourreau.Gilles.FaxDotNet.dll	cccc)	Wedge.pdb
vvv)	Tourreau.Gilles.FaxDotNet.pdb		

### **4. Intended Audience**

This document is created by the RIS Development team for the RadNet RIS Implementation team.

### **5. Who is affected**

Build 41.1 installs. This server update must be applied directly to the build 41.1 core install. Server updates are not cumulative. Any previous server updates must be applied before this one.

## **RADNET, Inc. – eRAD RIS Server Update**

### **6. eRAD RIS Release Version Numbers**

<b>Build</b>	<b>Patch</b>	<b>UI Version</b>	<b>Core Version</b>	<b>WS Version</b>	<b>DB Version</b>	<b>Notes</b>
40	10	1.3.40.17435	1.3.40.16659	1.3.40.17435	1.3.40.17370	<b>64 BIT web service.</b> GUI.zip, Webservice.
40	11	1.3.40.17435	1.3.40.16659	1.3.40.17525	1.3.40.17370	Web service update
40	12	1.3.40.17685	1.3.40.16659	1.3.40.17685	1.3.40.17648	GUI.zip, Web service, database updates
40	13	1.3.40.17776	1.3.40.16659	1.3.40.17776	1.3.40.17648	GUI.zip, Web service
40	14	1.3.40.18063(3GB)	1.3.40.16659	1.3.40.18063	1.3.40.18063	GUI.zip, Web service, database updates
40	14.1	1.3.40.18063(3GB)	1.3.40.16659	1.3.40.18107	1.3.40.18063	Web service
40	15	1.3.40.18127(3GB)	1.3.40.16659	1.3.40.18127	1.3.40.18127	GUI.zip, Web service, database updates
40	16	1.1.3.40.18585 (3GB)	1.3.40.16659	1.3.40.18585	1.3.40.18646	GUI.zip, Web service, database updates
40	16.1	1.3.40.18705 (3GB)	1.3.40.16659	1.3.40.18585	1.3.40.18646	GUI.zip
40	17	1.3.40.18780 (3GB)	1.3.40.18756	1.3.40.18756	1.3.40.18646	<b>Core install required for Radiologist User Group.</b> GUI.zip, Core setup, and Web Service updates
40	18	1.3.40.19253(3GB)	1.3.40.18756	1.3.40.18756	1.3.40.19253	GUI.zip, Web Service and database updates
40	19	1.3.40.19909(3GB)	1.3.40.18756	1.3.40.19878	1.3.40.19253	GUI.zip, Web Service
40	20	1.3.40.20086(3GB)	1.3.40.18756	1.3.40.19878	1.3.40.19253	Gui.zip and Recognizer
40	21	1.3.40.20703(3GB)	1.3.40.18756	1.3.40.19878	1.3.40.19253	Gui.zip
40	22	1.3.40.20898(3GB)	1.3.40.18756	1.3.40.19878	1.3.40.19253	Gui.zip
41	-	1.3.41.19135(3GB)	1.3.41.19135	1.3.41.19135	1.3.41.19135	Full version release
41	1	1.3.41.19356(3GB)	1.3.41.19135	1.3.41.19346	1.3.41.19346	GUI.zip, Web Service and database updates
41	2	1.3.41.19620(3GB)	1.3.41.19135	1.3.41.19620	1.3.41.19346	GUI.zip and Web Service
41	2sf	1.3.41.19620SF(3G)	1.3.41.19135	1.3.41.19620	1.3.41.19346	GUI.zip
41	2.1	1.3.41.20622(3GB)	1.3.41.19135	1.3.41.20617	1.3.41.19346	GUI.zip and Web Service
41	3	1.3.41.20171(3GB)	1.3.41.19135	1.3.41.20171	1.3.41.19346	GUI.zip , Web Service and Recognizer
41	4	1.3.41.20285(3GB)	1.3.41.19135	1.3.41.20171	1.3.41.19346	GUI.zip
41	5	1.3.41.20809(3GB)	1.3.41.20708	1.3.41.20792	1.3.41.20738	GUI.zip, Core install, Web Service and database updates
41	6	1.3.41.20912(3GB)	1.3.41.20708	1.3.41.20792	1.3.41.20738	GUI.zip
41.1	-	1.3.41.20918(3GB)	1.3.41.20791	1.3.41.20918	1.3.41.20739	GUI.zip, Core install, Web Service and DB updates, Help files
41.1	1	1.3.41.20987(3GB)	1.3.41.20791	1.3.41.20918	1.3.41.20739	GUI.zip
41.1	2	1.3.41.21162(3GB)	1.3.41.20791	1.3.41.21162	1.3.41.21150	GUI.zip, Web Service and DB updates,

## **7. Installing**

- 1) Copy and replace the current rRIS...zip file with the rRIS\_1.3.41.21162.zip file provided with this release.

Be sure to make a backup of the ris.exe.config file contained within the current .zip folder.

- 2) Stop each service listed below. Copy the files from the RIS Service Updates folder to the directory for each rRISService. Be sure to make backup of all files, but especially the rRISServices.exe.config and rRISServices.vshost.exe.config files.

*rRISService directories include:*

- a) *rRISService*
- b) *CDS rRISService*
- c) *Wedge rRISService*
- d) *DDE rRISService*

Restart each service.

- 3) Run the B41\_1\_002.sql file against the eRAD RIS database instance.