



eRAD RIS

# RELEASE ANNOUNCEMENT

Build 4.2023.046

UPDATED JUNE 27, 2023

## TABLE OF CONTENTS

**Summary** ..... 3

    Release Announcement ..... 3

    New Features ..... 3

    Resolved Issues ..... 3

**Feature Details** ..... 4

    Digital Forms ..... 4

        Feature #33815 – Display most appropriate answer when portals share a Digital Forms question ..... 4

    System Config ..... 5

        Feature #31704, 33953 – Additional granularity options for Question Details by Portal type ..... 5

**Version Details** ..... 7

    Code Stream ..... 7

## PUBLICATION HISTORY

Revision	Author	Description
June 27, 2023	Kevin Brooks / Hilary Saltmarsh	<ul style="list-style-type: none"> <li>▪ Commercial release.</li> </ul>

# SUMMARY

## Release Announcement

*This release of eRAD RIS 4.2023.046 introduces enhancements to the Portal procedure picker.*

## New Features

This release introduces the following features and enhancements:

Category	Redmine #	Subject	Description
Digital Forms	33815	Display most appropriate answer when portals share a Digital Forms question	This enhancement to the Digital Form Editor adds the ability to display the most appropriate answer in a digital form when the same question appears in both the Patient and Referring portals.
System Config	31704	Additional granularity options for Question Details by Portal type	This enhancement to System Config adds additional options to the Question Details and allows for separation of configuration between Portals.
System Config	33953	Additional Granularity to Question Details by Portal Type (Using in Portal Side)	This enhancement to System Config adds additional options to the Question Details and allows for separation of configuration between Portals.

SORTED BY CATEGORY AND REDMINE

Refer to the FEATURE DETAILS section below for configuration and usage information.

## Resolved Issues

This release resolves the following issues:

Category	Redmine #	Subject
Digital Forms	34155	Resolved a Digital Forms issue where Portal Call Center only questions were not displaying the PortalAdditionalInfoToCall message.

# FEATURE DETAILS

## Digital Forms

### Feature #33815 – Display most appropriate answer when portals share a Digital Forms question

#### Summary

This enhancement to the Digital Form Editor adds the ability to display the most appropriate answer in a digital form when the same question appears in both the Patient and Referring portals.

#### Background

Previously, questions created via the Procedure Picker could be used by both the Patient and Referring Portal. A referring could create an order for an exam in the Provider Portal and answer one of these questions, and the patient viewing the order to schedule in the Patient Portal would see the answer from the referring portal and can update it. However, Digital Forms must map to either the Patient Portal or Referring Portal version of the question, and do not distinguish which is the most recent answer to the question.

The recently resolved **FEATURE #30297** and **BUG #32759** now allows these responses to be viewed and edited in digital forms.

Note that these answers may be different as changes made in one location will not be reflected in the other.

#### Feature Description

Since questions can be answered differently by the referring doctor and the patient, a new read-only **Portal Question Recent Answer** column has been added to the Digital Form Editor that prioritizes the patient response - it will always display the Patient Portal response if one exists, otherwise it will use the Provider Portal response.

The screenshot shows a web interface for editing digital forms. At the top, there are two tabs: 'Page 1' and 'Page 2'. Below the tabs is a section titled 'New Section'. Inside this section, there are three rows of data, each with a label on the left and a text input field on the right. The first row is labeled 'RP Answer' and has 'No' in the input field. The second row is labeled 'PP Answer' and has 'Yes' in the input field. The third row is labeled 'Portal Recent Answer' and has 'Yes' in the input field.

THE SAME QUESTION MAPPED TO THREE SOURCES.

#### Configuration Instructions

No System Administrator actions are necessary to enable this feature.

## System Config

### Feature #31704, 33953 – Additional granularity options for Question Details by Portal type

#### Summary

This enhancement to System Config adds additional options to the **Question Details** and allows for separation of configuration between Portals.

#### Feature #31704 – Additional granularity options for Question Details by Portal type

With this change, new sections have been created for each portal within the **Procedure Picker** » **Edit Questions** dialog. This provides the ability to independently update questions with respect to the desired portal environment.

The screenshot shows the 'Question Details' configuration window. It includes fields for Code (weight), Short text (Weight (lbs)), Type (YesNoUnsure), and Max length (0). Below these is a 'Portal Scheduling Overrides' table with columns for Practice, Selected Answer, and Padding (Hours). The table has one row with 'Yes' and '1'. The 'Referring Portal' section has a 'Required' checkbox checked, 'Referring text' (What is the patient's current weight (lbs)?), 'Order only on:' (Yes, No, Unsure), 'Call us on:' (Yes, No, Unsure), 'Subnote' (sometext), and 'Display subnote on:' (Yes, No, Unsure). The 'Patient Portal' section also has a 'Required' checkbox checked, 'Patient text' (What is your current weight (lbs)?), 'Call us on:' (Yes, No, Unsure), 'Subnote' (hj), and 'Display subnote on:' (Yes, No, Unsure).

#### Feature #33953 – Additional Granularity to Question Details by Portal Type (Using in Portal Side)

Previously, questions used on the Referring and Patient portals shared the same "question details" section, meaning there was no way to independently configure question details on both Portals.

With this change, there are new options:

- Both the new **Referring Portal** and **Patient Portal** sections can be noted as **Required**,
- Both the **Call us on:** and **Display subnote on:** questions now have been updated from a **Yes/No** radiobuttons to the checkbox states of **Yes**, **No**, and **Unsure** (which allows combinations such as **Yes/Unsure**).
- A **Subnote:** phrase can be displayed on checkbox states of **Yes**, **No**, and **Unsure**.

- A new `Order only on:` option for the Referring Portal limits the ability to schedule. Note that orders created due to the `Order only on:` question will be promoted with the `Call the Call Center` prompt in the Patient Portal.

What is the patient's current weight (lbs)?  
(lbs)

Is the patient pregnant?  
**YES** NO NOT SURE

Confirm test results.

This request will be saved as an order and cannot be scheduled

Subnote on "Yes"

Order Only on "Yes"

The screenshot shows a form with a weight input field, a pregnancy question with radio buttons (YES, NO, NOT SURE), and a confirmation message. Two callout boxes point to the YES button and the confirmation message.

## Configuration Instructions

No System Administrator actions are necessary to enable this feature.

# VERSION DETAILS

## Code Stream

The following source code branches have been merged into this release:

