


Teaching Folder

Users can have a teaching or personal folder for interesting cases.

From any worklist, select folder from the context menu.




The screenshot shows a software interface with a worklist table. At the top, there are tabs for 'Tables' and 'Pending Dictation WL (11 of 116)'. Below the tabs, there are buttons for 'Read', 'Dr. Washington', 'Mammo', and 'MR'. The table has two columns: 'Patient Name' and 'Procedures'. A context menu is open over the first row of the table, listing various actions such as 'Dictate', 'View Study', 'Patient Folder', 'Print Forms', 'Attachments', 'Image Request', 'Add to Cache', 'Assign to Chair', 'Assign to Radiologist', 'Change Assistant Radiologist', 'Add to Folder', 'Mark As Read', 'Set STAT Read', 'Set Suspended Flag', 'Audit History', and 'Copy Cell'.

Patient Name	Procedures
King, Clarice	MA Digital Mammo Screen Bilat [G0202]
St, Hilary	MA Digital Mammo Screen Bilat [G0202] - Brea
nanoff, Na...	MA Digital Mammo Screen Bilat [G0202] - Brea
opperidge,...	MA Digital Mammo
mpler, Mer...	MA Digital Mammo
olia, Julia	MA Digital Mammo
nks, Melissa	MA Digital Mammo
olia, Julia	MA Digital Mammo
olia, Julia	MA Digital Mammo
olia, Julia	MA Digital Mammo
okett, Allison	MA Digital Mammo

Create a new folder or add to an existing folder

Teaching Folder



Folders can be accessed in the File menu

