

Server Update

For eRAD RIS
Version 1.0
Build 41
Update 41.1.001

Table of Contents

| | |
|---|---|
| 1. Defects Resolved | 3 |
| 2. Code Stream..... | 4 |
| 3. Package Contents | 5 |
| 4. Intended Audience | 5 |
| 5. Who is affected | 5 |
| 6. eRAD RIS Release Version Numbers | 6 |
| 7. Installing | 7 |


RADNET, Inc. – eRAD RIS Server Update

1. Defects Resolved

| Redmine # | Subject | Priority |
|------------------|--|-----------------|
| 5781 | Signed Pending signature items won't move | High |
| 5782 | Automacros containing fields lose their field attributes | Immediate |

2. Code Stream

eRad RIS – source code evolution



Legend:

Green = Released software

Red = Will not be released

Blue = Future release

3. Package Contents

Files included for this server update:

- 1) GUI ZIP Updates folder contains
 - a) rRIS_1.3.41.20987 (3GB).zip

4. Intended Audience

This document is created by the RIS Development team for the RadNet RIS Implementation team.

5. Who is affected

Build 41.1 installs. This server update must be applied directly to the build 41.1 core install. Server updates are not cumulative. Any previous server updates must be applied before this one.

6. eRAD RIS Release Version Numbers

| Build | Patch | UI Version | Core Version | WS Version | DB Version | Notes |
|--------------|--------------|----------------------|---------------------|-------------------|-------------------|--|
| 40 | 10 | 1.3.40.17435 | 1.3.40.16659 | 1.3.40.17435 | 1.3.40.17370 | 64 BIT web service. GUI.zip, Webservice. |
| 40 | 11 | 1.3.40.17435 | 1.3.40.16659 | 1.3.40.17525 | 1.3.40.17370 | Web service update |
| 40 | 12 | 1.3.40.17685 | 1.3.40.16659 | 1.3.40.17685 | 1.3.40.17648 | GUI.zip, Web service, database updates |
| 40 | 13 | 1.3.40.17776 | 1.3.40.16659 | 1.3.40.17776 | 1.3.40.17648 | GUI.zip, Web service |
| 40 | 14 | 1.3.40.18063(3GB) | 1.3.40.16659 | 1.3.40.18063 | 1.3.40.18063 | GUI.zip, Web service, database updates |
| 40 | 14.1 | 1.3.40.18063(3GB) | 1.3.40.16659 | 1.3.40.18107 | 1.3.40.18063 | Web service |
| 40 | 15 | 1.3.40.18127(3GB) | 1.3.40.16659 | 1.3.40.18127 | 1.3.40.18127 | GUI.zip, Web service, database updates |
| 40 | 16 | 1.1.3.40.18585 (3GB) | 1.3.40.16659 | 1.3.40.18585 | 1.3.40.18646 | GUI.zip, Web service, database updates |
| 40 | 16.1 | 1.3.40.18705 (3GB) | 1.3.40.16659 | 1.3.40.18585 | 1.3.40.18646 | GUI.zip |
| 40 | 17 | 1.3.40.18780 (3GB) | 1.3.40.18756 | 1.3.40.18756 | 1.3.40.18646 | Core install required for Radiologist User Group. GUI.zip, Core setup, and Web Service updates |
| 40 | 18 | 1.3.40.19253(3GB) | 1.3.40.18756 | 1.3.40.18756 | 1.3.40.19253 | GUI.zip, Web Service and database updates |
| 40 | 19 | 1.3.40.19909(3GB) | 1.3.40.18756 | 1.3.40.19878 | 1.3.40.19253 | GUI.zip, Web Service |
| 40 | 20 | 1.3.40.20086(3GB) | 1.3.40.18756 | 1.3.40.19878 | 1.3.40.19253 | Gui.zip and Recognizer |
| 40 | 21 | 1.3.40.20703(3GB) | 1.3.40.18756 | 1.3.40.19878 | 1.3.40.19253 | Gui.zip |
| 40 | 22 | 1.3.40.20898(3GB) | 1.3.40.18756 | 1.3.40.19878 | 1.3.40.19253 | Gui.zip |
| 41 | - | 1.3.41.19135(3GB) | 1.3.41.19135 | 1.3.41.19135 | 1.3.41.19135 | Full version release |
| 41 | 1 | 1.3.41.19356(3GB) | 1.3.41.19135 | 1.3.41.19346 | 1.3.41.19346 | GUI.zip, Web Service and database updates |
| 41 | 2 | 1.3.41.19620(3GB) | 1.3.41.19135 | 1.3.41.19620 | 1.3.41.19346 | GUI.zip and Web Service |
| 41 | 2sf | 1.3.41.19620SF(3G) | 1.3.41.19135 | 1.3.41.19620 | 1.3.41.19346 | GUI.zip |
| 41 | 2.1 | 1.3.41.20622(3GB) | 1.3.41.19135 | 1.3.41.20617 | 1.3.41.19346 | GUI.zip and Web Service |
| 41 | 3 | 1.3.41.20171(3GB) | 1.3.41.19135 | 1.3.41.20171 | 1.3.41.19346 | GUI.zip , Web Service and Recognizer |
| 41 | 4 | 1.3.41.20285(3GB) | 1.3.41.19135 | 1.3.41.20171 | 1.3.41.19346 | GUI.zip |
| 41 | 5 | 1.3.41.20809(3GB) | 1.3.41.20708 | 1.3.41.20792 | 1.3.41.20738 | GUI.zip, Core install, Web Service and database updates |
| 41 | 6 | 1.3.41.20912(3GB) | 1.3.41.20708 | 1.3.41.20792 | 1.3.41.20738 | GUI.zip |
| 41.1 | - | 1.3.41.20918(3GB) | 1.3.41.20791 | 1.3.41.20918 | 1.3.41.20739 | GUI.zip, Core install, Web Service and DB updates, Help files |
| 41.1 | 1 | 1.3.41.20987(3GB) | 1.3.41.20791 | 1.3.41.20918 | 1.3.41.20739 | GUI.zip |

7. Installing

- 1) Copy and replace the current rRIS...zip file with the rRIS_1.3.41.20987.zip file provided with this release.
Be sure to make a backup of the ris.exe.config file contained within the current .zip folder.