



# eRAD RIS

## RELEASE ANNOUNCEMENT

### Rapid Release 4.2023.068.0.2

UPDATED SEPTEMBER 14, 2023

#### TABLE OF CONTENTS

**Summary** ..... 2

    Release Announcement .....2

    New Features .....2

**Feature Details** ..... 3

    Provider Images .....3

        Feature #33608 – Add Quick Launch Icon on Images Tab.....3

        Feature #34343 – Enhanced user messaging while retrieving archived studies 3

**Version Details** ..... 5

    Code Stream .....5

#### PUBLICATION HISTORY

Revision	Author	Description
September 14, 2023	Kevin Brooks, Michelle Mahabir / Hilary Saltmarsh	<ul style="list-style-type: none"> <li>Commercial release.</li> </ul>

# SUMMARY

## Release Announcement

*This rapid release of eRAD RIS 4.2023.068.0.2 introduces Provider Portal image retrieval enhancements.*

## New Features

This release introduces the following features and enhancements:

Category	Redmine #	Subject	Description
Provider Images	33608	Add Quick Launch Icon on the Images tab	This enhancement to Provider Portal adds a Quick Launch icon to directly view studies within the Portal.
Provider Images	34343	Enhanced user messaging while retrieving archived studies	This enhancement to Provider Images provides an informational message to users while attempting to retrieve archived studies.

Refer to the FEATURE DETAILS section below for configuration and usage information.

# FEATURE DETAILS

## Provider Images

### Feature #33608 – Add Quick Launch Icon on Images Tab


#### Summary

This enhancement to Provider Portal adds a Quick Launch icon to directly view studies within the Portal.

#### Feature Description

With this change, a new **Quick Launch** icon has been added to the **Images tab** within the Provider Portal. This will allow the user to launch full studies in the default viewer.

Previously, there was no way to open the default viewer until the thumbnails were loaded.



The images for this exam have been archived. We are retrieving the images and will display them below when available.



A NEW QUICK LAUNCH ICON ON THE RIGHT-HAND SIDE OF THE IMAGES TAB WILL LAUNCH FULL STUDIES.

## Configuration Instructions

No System Administrator actions are necessary to enable this feature.

### Feature #34343 – Enhanced user messaging while retrieving archived studies


#### Summary

This enhancement to Provider Images provides an informative message to users while attempting to retrieve archived studies.

#### Feature Description

With the introduction of V9 PACS, image thumbnails are no longer immediately displayed on the portal for studies that have been archived.

With this change, when a user clicks the **Images** tab of a study that has been archived, the system will present an informative message indicating that the images are not yet accessible while the thumbnails are loaded.



Text of the message is customizable via a new `PACSIImagesPreparingMessage` application configuration setting.

Note that if the study has not been archived the thumbnails will load as normal when the `Images` tab is selected.

## Known Limitations

While there are no Known Limitations for this feature, there is additional functionality planned for a future release that will replace the `PACSIImagesPreparingMessage` application configuration setting for customizable with a paragraph configuration setting to allow RIS Administrators to modify it without Service Team assistance:

- FEATURE #34372 – CREATE PARAGRAPH CONFIG FOR FEATURE# 34343

## Configuration Instructions


No System Administrator actions are necessary to enable this feature; however, optional configuration is available:

- Optionally customize the `PACSPreparingImagesMessage` message that displays to user while images are being prepared. Requires Service Team assistance (see KNOWN LIMITATIONS).

# VERSION DETAILS

## Code Stream

The following source code branches have been merged into this release:



**WARNING:** *The functionality of this release has not yet been merged forward into any other release.*

*Please carefully review the code stream diagram of any future releases before deploying to ensure that the functionality of this release will not be lost when upgrading.*