



eRAD RIS

# Commercial Management Reports

---

For v3.2022.2.14

UPDATED FEBRUARY 24, 2022

# Table of Contents

Management Reports \ Billing .....	8
Charge Audit All Status .....	8
Charge Audit Report .....	9
Charge Audit Report no header .....	10
Management Reports \ Connect (Requires Patient Portal) .....	11
Connect Patient Portal Weekly Activity .....	11
Management Reports \ Insurance Eligibility (Requires Eligibility) .....	12
Carriers Not Using Insurance Eligibility .....	12
Eligibility Failed Requests .....	13
Eligibility Time to Complete (seconds) .....	14
Insurance Eligibility Usage .....	15
Payers Not Responding .....	16
Management Reports \ Maintenance .....	17
Audit Log .....	17
Audit Log by Patient .....	18
BI-RAD to Recommendation Mapping .....	19
Confirmed Radiologist Audit .....	20
Device Availability .....	21
Document Distribution Log .....	22
EMR Subscriber Listing .....	23
FAX Error Log .....	24
Last Login by User Id .....	25
Message Queues .....	26
Modality to CPT to Procedure Mapping .....	27
New Referring Physicians by Created Date .....	28
Procedure Codes .....	29
Procedure Template Mappings .....	30
Procedures without Procedure Group .....	31
Questionnaire Data Dictionary .....	32

---

Radiologist Setup .....	33
Referring Physician List .....	34
Referring Physician Notification.....	35
Report template and macro export .....	36
Schedule Group Mappings.....	37
Site Procedure Mappings.....	38
Top RIS Errors By Date .....	39
Top RIS Errors By Resource Type .....	40
Top RIS Errors By UserID .....	41
Top RIS Errors Since Yesterday .....	42
User Group Permissions.....	43
User List.....	44
User Permissions.....	45
Users in Usergroups .....	46
XR-29 Compliance .....	47
Management Reports \ Mammography.....	48
Mammo Annual Radiologist Summary .....	48
Mammo Annual Technologist Summary.....	49
Mammo Back on Track .....	50
Mammo Document Distribution Monitoring.....	51
Mammo Ignored Reminder List .....	52
Mammo Orphaned Orders.....	53
Mammo Radiologist Summary.....	54
Mammo Technologist Summary.....	55
Mammography Missing Data.....	56
MQSA General Exam List .....	57
MQSA Layletter Audit .....	58
MQSA List of Exams Performed by Technologist.....	59
MQSA List of False Positive Exams.....	60
MQSA List of False Negative Exams .....	61
MQSA List of Positive Exam Followups (Stand Alone System).....	62
MQSA List of Positive Exam Followups .....	63

---

MQSA List of Positive Exams with Outstanding Biopsies.....	64
MQSA Location Summary Group By BI-RADS .....	65
MQSA Location Summary Group By Procedure.....	66
MQSA Outcome Report .....	67
MQSA Radiologist Reporting Summary .....	68
Management Reports \ Meaningful Use .....	69
Automated Measures by Radiologist.....	69
Automated Measures Raw Data .....	70
Automated Measures .....	71
Automated Measures Details .....	72
CQM Measures .....	73
Patient Lists.....	74
Management Reports \ Operations.....	75
Aborted Exams.....	75
Addended Mammography Exams.....	76
Appointment Schedule .....	77
Appointment Schedule by Procedure Group.....	78
Blood Work Appointment Schedule .....	79
Cancelations by Referring .....	80
Canceled Exams .....	81
Canceled Reasons .....	82
Check-in Counts .....	83
Confirmation Status Canceled or Rescheduled .....	84
Confirmed Appointments .....	85
Daily Schedule.....	86
Daily Schedule Add Ons .....	87
Exam Volume Breakdown by Payor .....	88
Finding Follow Up .....	89
Front Desk Activity .....	90
Image Request History.....	91
IVT 7 Day Trend.....	92
IVT Idle Time .....	93

---

IVT Turnaround .....	94
Modality Duration Override.....	95
MR Breast Dictated .....	96
Need Plan .....	97
NM Therapy Patient Listing.....	98
No Shows .....	99
Open Capacity .....	100
OTC QA .....	101
Patient Email and Cellphone Percentages .....	102
Patient Wait Time Trending by Site .....	103
Patient Wait Times.....	104
Patients without Insurance .....	105
Payments Received .....	106
Pre-Auth Referring Physician .....	107
Reschedules and Cancelations.....	108
Referring Physician Details.....	109
Room Closures .....	110
Room Utilization .....	111
RTAT by Radiologist.....	112
RTAT By Referring .....	113
RTAT Detail.....	114
RTAT Detail – Monthly .....	115
RTAT Detail Mammo .....	116
RTAT RAD Signing.....	117
RTAT Summary.....	118
Scheduled vs Performed .....	119
Special Accommodations .....	120
Time to Schedule (seconds) .....	121
Time to Schedule by User .....	121
Transportation Report .....	123
Unfinished Exams.....	124
Unknown Referring.....	125

---

Unread Exams .....	126
Wait Room Times .....	127
Wait Room Turnover .....	128
Management Reports \ Peer Review .....	129
Peer Review Report .....	129
Review Summary by Radiologist .....	130
Management Reports \ Sales .....	131
Referrals by Physician .....	131
Referrals by Physician Patient List .....	132
Referring Physician Volumes .....	133
Management Reports \ Statistics .....	134
Management Reports \ Statistics \ Editors .....	134
Edited Jobs Activity Report .....	134
Edited Jobs By User Detail .....	135
Editors Activity .....	136
Management Reports \ Statistics \ IVT .....	137
IVT User Activity .....	137
Management Reports \ Statistics \ Radiologists .....	138
Dictations by User .....	138
QA Radiologist Review .....	139
Radiologist Dictation Activity by Date .....	140
Radiologist Productivity by User - Mammo .....	141
Radiologist Signing Activity .....	142
Management Reports \ Statistics \ Registration .....	143
Front Desk Detailed Activity by Hour .....	143
Front Desk Detailed Activity .....	144
Management Reports \ Statistics \ Reviewers .....	145
QA Review Statistics .....	145
Management Reports \ Statistics \ Schedulers .....	146
Scheduler Activity .....	146
Scheduling Event Trend .....	147
Management Reports \ Statistics \ Technologists .....	148

Technologist Activity by Modality.....	148
Management Reports \ Utilization Management.....	149
Coverage by Medical Group.....	149
Inbound Document Activity .....	150
Medical Group Activity .....	151
Referring Physician Activity.....	152
Utilization Review User Activity .....	153
UM Daily Dashboard .....	154

## Publication History

Revision	Author	Description
October 24, 2019		<ul style="list-style-type: none"> <li>▪ Released with v3.2018.5</li> </ul>
February 24, 2022	Kevin Brooks	<ul style="list-style-type: none"> <li>▪ Republished for v3.2022.2.14. Content is unchanged.</li> </ul>

## Management Reports \ Billing

### Charge Audit All Status

Description: This billing report returning DOS, Patient Name, MRN, Accession, Stats, Location and CPT Codes for the billing period specified.

Parameters:

- Practice, Site(s), Date Range, Test Patients (to be excluded), Excluded billing codes, Excluded Patient IDs and Status Code(s) included.

### Charge Audit All Status

**Practice:** Treesdale Radiology

**Site(s):** Cedarville, Mapleton, Pine Grove, Outside Read1 (outside read)

**Date Range:** 11-04-2016 - 11-08-2016

*Excluding billing codes 99999*

*Excluded patient ids: 107761*

*Excluding patients where first or last names equal test.*

DOS	Patient Name	MRN	Accession	Status	Location	CPT
11-04-2016	Parsons, Hilary	8211	1028891	Signed1	Cedarville	71030
11-04-2016	Strickland, Rhonda	9695	1028889	Signed1	Cedarville	70460
11-04-2016	Rodriguez, Rhonda	4454	1028890	Scheduled	Cedarville	72126
11-04-2016	Hartman, Loren	1378	1028888	Scheduled	Pine Grove	74160
11-04-2016	O'Dell, Jayne	z40000004	1029890	Signed1	Cedarville	70480
11-07-2016	Hartman, Loren	1378	1029889	Scheduled	Pine Grove	74160
11-07-2016	Freddieson, Fred	y6	1027871	Scheduled	Mapleton	70450
11-07-2016	Harrison, Harry	y8	1027874	Scheduled	Pine Grove	72141
11-07-2016	Strickland, Rhonda	9695	1030889	ExamDone	Cedarville	74160
<b>Count</b>						<b>9</b>

Charge Audit All Status

Confidential and Proprietary  
Page 1 of 1

Date Run: 11-08-2016 2:59 PM  
RRIS-RSNA\FORMS\GuestUser



## Charge Audit Report

Description: This billing report returning patient and billing information for the billing period specified.

Parameters:

- Practice, Site(s), Date Range

### Charge Audit Report

Practice(s): Treesdale Radiology

Site(s): Cedarville, Mapleton, Pine Grove

Date Range: 01-31-2013 - 11-08-2016

DOS	Patient Name	MRN	SSN	Encounter #	Accession #	Status	Location	CPT	CPT Modifier
11-28-2013	Mammo, Mary	1501447RSNA		35881	1018126	Signed1	Tree\Pine	G0202	
11-28-2013	Mammo, Sally	1501455RSNA		36030	1018269	Signed1	Tree\Pine	G0202	
11-27-2015	Warren, Selma	1501585RSNA		42549	1024376	Signed1	Tree\Pine	G0202	Left

Charge Audit Report

Confidential and Proprietary  
Page 1 of 1

Date Run: 11-08-2016 3:01 PM  
RRIS-RSNA\FORMS\GuestUser

**Charge Audit Report no header**

Description: Identical to the Charge Audit Report but excluding the report header for the purpose of exporting to a csv for external billing systems.

Parameters:

- Practice, Site(s), Date Range

DOS	Patient Name	MRN	SSN	Encounter #	Accession #	Status	Location	CPT	CPT Modifier
01-29-2013	Bean, Gene	1501358RSNA		32953	1014982	Signed1	Tree\Cedar	74020	
11-28-2013	Mammo, Mary	1501447RSNA		35881	1018126	Signed1	Tree\Pine	G0202	
11-28-2013	Mammo, Sally	1501455RSNA		36030	1018269	Signed1	Tree\Pine	G0202	
11-27-2015	Warren, Selma	1501585RSNA		42549	1024376	Signed1	Tree\Pine	G0202	Left

## Management Reports \ Connect (Requires Patient Portal)

### Connect Patient Portal Weekly Activity

Description: This report displays usage counts for the Connect Patient Portal based on selected start and end dates.

Parameters:

- Site Group(s) and Date Range

### Patient Portal Weekly Activity

Site Group(s): Treesdale Radiology


Date Range: (11-01-2016 - 11-07-2016)

Description: This report displays usage counts for the Connect Patient Portal based on selected start and end dates.

Note: Last Month is the month prior to the selected start date while "12 Months" takes in account for example (Start date: July 1, 2016 presents counts for June 1, 2015 to June 30, 2016).

User Accounts	Count
Total Accounts	112
Active Accounts	14
Disabled Accounts	1
Inactive Accounts (No activity for 13 months)	97

Login Activity	Count	Last Month	12 Months
Logins	1	16	51
Unique Logins	0	2	4
New Account Logins	1	3	5
Failed Logins due to Password	0	8	14
Failed Logins due to wrong account	1	6	59
Failed Logins due to disabled account	0	3	3
Accesses via iCode	0	0	0



Help Request Activity	Count	Last Month	12 Months
Get Help Requests	0	2	2
Requested Technical Assistance	0	2	2
Standard Messages Sent	0	0	0
Exam Messages Sent	0	0	0
Requested Technical Assistance Outside of Login	0	0	0

Misc. Activity	Count	Last Month	12 Months
Authorized Access Invites Sent	0	0	0
Authorized Access Invites Validated	0	0	0

## Management Reports \ Insurance Eligibility (Requires Eligibility)

### Carriers Not Using Insurance Eligibility

Description: This report lists the carriers not using insurance eligibility and their usage (count) grouped by practice and site. This report also excludes carriers where the carrier’s value for self-pay flag or attorney flag is set to “Y”. Minimum occurrence represents the minimum number of times the carrier is referenced before it will be pulled into the report.

Parameters:

- Practice(s), Site(s), Date Range and Minimum Occurrence

### Carriers Not Using Insurance Eligibility

**Practice(s):** Treesdale Radiology

**Site(s):** Cedarville, Mapleton, Pine Grove

**Date Range:** 01-01-2016 - 11-08-2016

**Description:** This report lists the carriers not using insurance eligibility and their usage (count) grouped by practice and site. This report also excludes carriers where the carriers' value for self pay flag or the attorney flag is set to Y. Min. occurrence represents the minimum number of times the carrier is referenced before it will be pulled into the report.

Practice	Site	Carrier Code	Carrier Description	Count
Treesdale Radiology				18
	☐ Cedarville			12
		☐ 0601397	BAC INSURANCE	5
		☐ 090665	Blue Cross of CA	5
		☐ 213574	Aetna	1
		1501740DEMO	Matina, May	
		☐ 050105	PARAMOUNT INSURANCE	1
	☐ Mapleton			4
		☐ 213574	Aetna	2
		☐ 090665	Blue Cross of CA	2
		1501741DEMO	Maltolia, Ratna	
		1501739DEMO	Mangia, Gertrude	
	☐ Pine Grove			2
<b>Grand Total</b>				<b>18</b>

Carriers Not Using Insurance Eligibility

Confidential and Proprietary  
Page 1 of 1

Date Run: 11-08-2016 1:47 PM  
RRIS-RSNA\FORMS\GuestUser

### Eligibility Failed Requests

Description: This report displays insurance eligibility carrier failed insurance requests.

Parameters:

- Practice, Site(s) and Date Range

### Eligibility Failed Requests

Practice(s): North Treesdale

Site(s): Cedarville

Date: 11-08-2016

Description: This report displays insurance eligibility carrier failed insurance requests.

Carrier ↕	Failed Requests	Total
<input type="checkbox"/> BCBS FEDERAL		1
<input type="checkbox"/> BCBS HEALTHY NY		1
<input type="checkbox"/> BCBS NY		13
<input type="checkbox"/> BLUE CHOICE HMO		4
<input type="checkbox"/> MEDICARE BLUE CHOICE		10
<input type="checkbox"/> MEDICARE NY		3
<input type="checkbox"/> MVP HEALTH INS HMO		1
<input type="checkbox"/> UNIVERA HEALTHCARE		1
<b>Total</b>		<b>34</b>

Eligibility Failed Requests

Confidential and Proprietary  
Page 1 of 1

Date Run: 11-08-2016 12:00 PM  
RADNET\EcksteinR

**Eligibility Time to Complete (seconds)**

Description: This report displays the amount of time for the eligibility service to complete an insurance eligibility request.

Parameters:

- Practice, Site(s) and Date Range

**Eligibility Time To Complete (seconds)**

**Practice:** Treesdale Radiology

**Site:** Cedarville, Mapleton, Pine Grove, No site identified

**Date Range:** 01-01-2016 - 11-08-2016

**Description:** This report displays the amount of time for the eligibility service to complete an insurance eligibility request. The report is based on the practices, sites and date range (requested date) parameters.

Carrier	(seconds)			Request times (seconds)						Total
	Avg	Min	Max	1	2	3-5	6-10	11-20	20+	
NATIONWIDE INSURANCE	1	0	2	2	3	0	0	0	0	5
<b>Total</b>				2	3	0	0	0	0	5

Eligibility Time To Complete (seconds)

Confidential and Proprietary  
Page 1 of 1

Date Run: 11-08-2016 1:50 PM  
RRIS-RSNAFormsGuestUser

### Insurance Eligibility Usage

Description: This report shows the usage for the insurance eligibility based on the requesting practice, requesting sites and the created date range selected. If no site is identified in the site parameter, this indicates the request did not have a site available at the time of the request.

Parameters:

- Practice, Site(s) and Date Range

#### Insurance Eligibility Usage

Practice: Treesdale Radiology

Date Range: 01-01-2016 - 11-08-2016

Description: This report shows the usage for the insurance eligibility based on requesting practice, requesting sites, and the created date range selected.

\*No site identified in the site parameter indicates the request did not have a site available at the time of the request.

Practice	Site	Visit Count	Total Requests	Non-chargeable	Failed Requests	Error Cost	Success Cost
Treesdale Radiology		2	2	0	0	\$0.00	\$0.10
	No site identified	2	2	0	0	\$0.00	\$0.10
<b>Total</b>		<b>2</b>	<b>2</b>	<b>0</b>	<b>0</b>	<b>\$0.00</b>	<b>\$0.10</b>

Insurance Eligibility Usage

Confidential and Proprietary  
Page 1 of 1

Date Run: 11-08-2016 1:58 PM  
RRIS-RSNA\FFormsGuestUser

### Payers Not Responding

Description: This report lists any payers not responding by date and hour for the period selected

Parameters:

- Date Range

#### Payers Not Responding

Date: 11-08-2016

Description: This report lists any payers not responding by date and hr for the period selected.

Completed Date	Completed Hr	Eligibility Payer ID	Payer Name	Sample Carrier	Sample Eligible ID	Exception Details	Pass Count	Error Count	Error Rate
11-08-2016	12	00002	Aetna   AHP of New Jersey Parsippany   AHP of New Jersey Purchase NY   Ohio Valley HMO Cincinnati (Aetna)   PHP of Georgia Atlanta Aetna	600 - AETNA	8O1W8F	-	0	16	100.00 %
11-08-2016	12	AID1	New York Medicaid	706 - METROPLUS HEALTH PLAN MCD	8O1W8F	-	0	5	100.00 %
11-08-2016	11	0008	Horizon Blue Cross Blue Shield of New Jersey	693 - HORIZON NJ HEALTH	8O1V95	Net.WebException: The oper@ion has timed out @ Net.WebClient.DownloadD@ainter (Uri address, WebRequest& request) @ Net.WebClient.DownloadString (Uri address) @ Net.WebClient.DownloadString (String address) @ rRIServices.InsuranceEligibility.EligibleAPI.HttpGet(String uri) in C:\BuildAgent3\... (url:1700a600b0d4b6f6a9a0a)	59	4	6.35 %



## Management Reports \ Maintenance

### Audit Log

Description: This report displays the audit information for the period and users specified.

Parameters:

- Date Range, User, Exclude Inactive Users (Y/N) and Include Login\Logout (Y/N)

Date\Time	Patient ID	Patient Name	Description	Audit Actions	Changes
03-06-2017 9:18 AM			Login	Login	
03-06-2017 9:19 AM	9424PE	Jones, Mary	UI_OrderRetrieved_UtilizationManagement	Query	
03-06-2017 9:22 AM			Logout	Logout	
03-06-2017 10:31 AM			Login	Login	
03-06-2017 11:10 AM			Login	Login	
03-06-2017 11:45 AM			Login	Login	
03-06-2017 1:36 PM			Login	Login	
03-06-2017 1:36 PM	9424PE	Jones, Mary	Patient retrieved	Query	
03-06-2017 1:38 PM			Login	Login	
03-06-2017 1:48 PM	9424PE	Jones, Mary	Patient retrieved	Query	
03-06-2017 1:48 PM	9424PE	Jones, Mary	Patient retrieved	Query	
03-06-2017 1:48 PM	9424PE	Jones, Mary	Patient retrieved	Query	
03-06-2017 1:48 PM	9424PE	Jones, Mary	Patient retrieved	Query	
03-06-2017 1:49 PM	9424PE	Jones, Mary	Patient retrieved	Query	
03-06-2017 1:49 PM	9424PE	Jones, Mary	Patient retrieved	Query	
03-06-2017 1:49 PM	9424PE	Jones, Mary	Patient retrieved	Query	
03-06-2017 1:51 PM	9424PE	Jones, Mary	Order submitted	Order submitted	<Changes><AddedRows><c_visit /><c_visit_x_patient_insurance /></AddedRows></Changes>
03-06-2017 1:51 PM	9424PE	Jones, Mary	Order submitted	Order submitted	<Changes><AddedRows><c_order /><c_order_item procedure_code="71260" /><c_order_item procedure_code="70210" /><c_order_item_certification /><c_order_item_certification /></AddedRows></Changes>

## Audit Log by Patient

Description: This report displays the audit information for the specified patient MRN.

Parameters:

- MRN

## Audit Log by Patient

Patient: (9696) Strickland, Rhonda

Description: This report displays the audit information for the specified patient MRN.

Log Date	Log Time	User ID	Username	description	Audit Actions
11-08-2016	3:23 PM	eckstein	Eckstein, Rhonda	Accessed audit history	Accessed audit history
	2:35 PM	barnold	Arnold, Bob	Clinical Summary Retrieval	Clinical Summary Retrieval
	2:35 PM	barnold	Arnold, Bob	Patient Exam List Viewed	Patient Exam List Viewed
	2:35 PM	barnold	Arnold, Bob	Study retrieved	Query
	2:35 PM	barnold	Arnold, Bob	Confirmation SMS Sent	Confirmation SMS Sent
	2:35 PM	barnold	Arnold, Bob	Study retrieved	Query
	2:35 PM	barnold	Arnold, Bob	Study scheduled	Add

### BI-RAD to Recommendation Mapping

Description: This report displays the BI-RAD, Recommendation and the BI-RAD to Recommendation mappings configured within the system. This report is made up of 3 tables, BI-RAD table, Recommendations Table and the mappings between those called BI-RAD to Recommendation Mappings.

Parameters:

- Practice

BI-RAD Table	
Birad Code	Birad Description
0	BI-RADS 0
1	BI-RADS 1
2	BI-RADS 2
3	BI-RADS 3
4	BI-RADS 4
5	BI-RADS 5
6	BI-RADS 6

Recommendations Table			
Recommendation Code	Recommendation Description	follow up days	procedure code
P	Additional projections		
testb40	test		0 CT999
M	Magnification views		
S	Spot Compression	0	

### BI-RAD to Recommendation Mappings

Birad		Recommendation	
Practice:			
0	BI-RADS 0	P	Additional projections
		M	Magnification views
		S	Spot Compression
		V	Spot magnification view(s)
		U	Ultrasound
		O	Old film comparision
		G	Ductography
		Z	Breast MRI
		R	Technical Repeat

### Confirmed Radiologist Audit

Description: This report displays the name and UserID of users who are identified as having the Confirmed Radiologist flag enabled.

Parameters:

- Date Range

### Confirmed Radiologist Audit

The following users are identified as having the Confirmed Radiologist flag enabled.

Person Name	NPI
Eckstein (Rad), Rhonda (rhondarad)	
Saltmarsh (Rad), Hilary (hilrad)	
<b>Count: 2</b>	

Confirmed Radiologist Audit

Confidential and Proprietary  
Page 1 of 4

Date Run: 11-08-2016 1:42 PM  
RRIS-RSNA\FORMS\GuestUser



### Device Availability

Description: This report displays the document distribution log based on parameters selected.

Parameters: This report shows the device availability for the practice, sites and rooms selected.

- Practice, Sites, Rooms and whether to include exceptions or not

#### Device Availability

Practice: Treescdale Radiology

Site(s): Cedarville

Rooms: CT1CE, MA1CE, MA2CE, MR1CE, MR2CE, MR3CE, US1CE, US2CE, XR1CE

Description: This report shows the device availability for the practice, sites and rooms selected.

		Monday		Tuesday		Wednesday		Thursday		Friday					
MR1CE		Start	End	Start	End	Start	End	Start	End	Start	End	Start	End	Start	End
	Default Template	12:00 AM	1:00 PM	12:00 AM	1:00 PM	12:00 AM	1:00 PM	12:00 AM	1:00 PM	12:00 AM	1:00 PM				
		Monday		Tuesday		Wednesday		Thursday		Friday					
MR3CE		Start	End	Start	End	Start	End	Start	End	Start	End	Start	End	Start	End
	Default Template	12:00 AM	1:00 PM	12:00 AM	1:00 PM	12:00 AM	1:00 PM	12:00 AM	1:00 PM	12:00 AM	1:00 PM				

Device Availability

Confidential and Proprietary  
Page 1 of 1

Date Run: 11-08-2016 3:43 PM  
RRIS-RSNA/FormsGuestUser

### Document Distribution Log

Description: This report displays the document distribution log based on parameters selected.

Parameters:

- Date Range, and Distribution Status

#### Document Distribution Log

Date: 11-08-2016

Status(s): Error

Description: This report displays the document distribution log based on parameters selected.

Date	Status	Acc #	Patient Name	MRN	Recipient	Fax #	Retries	Job id	Del. Method	Site Code	Referring Physician
11-08-2016							12				
	Error	10045	STRICKLAND, RHONDA	3508157R	ARNOLD, BOB,	(224)123-5432	0	201d239f17c	Fax	CED	ARNOLD, MD, BOB
Activity Log: 11/8/2016 1:54:54 PM -05:00 Fax was unsuccessful. Fax Status: fjsRETRIES_EXCEEDED, Extended Status: fjsFATAL_ERROR 11/8/2016 12:56:22 PM -05:00 No Dial Tone 11/8/2016 12:55:11 PM -05:00 Job created											

**EMR Subscriber Listing**

Description: This report lists all the active subscribers for the selected EMRs.

Parameters:

- EMR(s)

**EMR Subscriber Listing**

EMR(s): Cedar129

Description: This report lists all the active subscribers for the selected EMRs..

EMR	Subscriber	NPI	Count
.Cedar129			6
	ARNOLD, BOB, MD	1234563	
	SUMMERS, SUE, MD	4567821	
	MARIO, JAMES, MD	65434543	
	NOVAK, JAMES, MD	87678685	
	WAITE, ANDREW, MD	65565454	
	JAMES, JOHN, MD	45345654	
<b>Total</b>			<b>6</b>

EMR Subscriber Listing

Confidential and Proprietary  
Page 1 of 1

Date Run: 11-08-2016 3:56 PM  
RADNET\EcksteinR

**FAX Error Log**

Description: This report displays document distribution log based on all Faxes in Error. This report is similar to the document distribution log, but with a focus on FAX errors

Parameters: (No parameters)

**FAX Error Log**

Description: This report displays document distribution log based on all Faxes in Error.

Date	Status	Acc #	Patient Name	MRN	Recipient	Fax #	Retries	Job Id	Del. Method
01-22-2014	Error	038008	JONES, JAMES	1387	ARNOLD, BOB	2421234512			Fax
Activity Log: 1/22/2014 2:34:51 PM-05:00 This fax number has been determined to have no answer from another work item.									
		038008	JONES, JAMES	1387	ARNOLD, BOB	2421234512			Fax
Activity Log: 1/22/2014 2:35:13 PM-05:00 This fax number has been determined to have no answer from another work item.									



### Last Login by User Id

Description: This report displays the audit information capturing the last login for each user with the RIS.

Parameters: (No parameters)

### Last Login by User Id

Description: This report displays the audit information capturing the last login for each user with the RIS.

User Id	User Name	Active	Last Logon
accesstest	Test, Access	Y	12-24-2015 11:20 AM
administrator	Test User, System Admin	Y	10-11-2012 8:25 AM
admintest	Test User, Administrator	Y	07-05-2012 7:23 AM
Andrew	Waite, Andrew R	Y	11-18-2015 7:51 AM
barnold	Arnold, Bob	Y	11-08-2016 2:33 PM
bob	bob, bob	Y	03-18-2016 12:08 PM
clifton	Harding, Clifton	Y	01-20-2016 8:43 AM

### Message Queues

Description: This report displays the messaging queues and current counts.

Parameters: (No parameters)

### Message Queues

Description: This report displays the messaging queues and current counts.

Message Queues	Count
DD	0
DD_Fix	0
MIRTH	0
<b>Total</b>	<b>0</b>

### Modality to CPT to Procedure Mapping

Description: This report displays the available procedures mapped to billing (CPT) codes per modality and selected scheduling group.

Parameters:

- Schedule Group(s)

### Modality to CPT to Procedure Mapping

Schedule Group(s): TreesdaleSchedGroup

Description: This report displays the available procedures mapped to billing (cpt) codes per modality and selected scheduling group.

Schedule Group	Modality	CPT Code	CPT Description	Proc Code	Procedure Description	
TreesdaleSchedGroup						2241
	(CT1CE) - CT 1 at Cedar					111
		99999	Non-Billable Code	CT50	CT Recall for Technical Repeat	
		74170	COMPUTED TOMOGRAPHY, ABDOMEN; WITHOUT CONTRAST MATERIAL, FOLLOWED BY CONTRAST MATERIAL(S) AND FURTHER SECTIONS	CT13	CT 3 Phase	
		74160	COMPUTED TOMOGRAPHY, ABDOMEN; WITH CONTRAST MATERIAL(S)	74160	CT Abdomen W	
		71275	COMPUTED TOMOGRAPHIC ANGIOGRAPHY, CHEST (NONCORONARY), WITH CONTRAST MATERIAL(S), INCLUDING NONCONTRAST IMAGES, IF PERFORMED, AND IMAGE POSTPROCESSING	CT66	CT Abdomen W & CTA Chest WWo	
		74160	COMPUTED TOMOGRAPHY, ABDOMEN; WITH CONTRAST MATERIAL(S)	CT9	CT Abdomen W & Pelvis W	
		71275	COMPUTED TOMOGRAPHIC ANGIOGRAPHY, CHEST (NONCORONARY), WITH CONTRAST MATERIAL(S), INCLUDING NONCONTRAST IMAGES, IF PERFORMED, AND IMAGE POSTPROCESSING	CT54	CT Abdomen W & Pelvis W & CTA Chest WWo	
		74170	COMPUTED TOMOGRAPHY, ABDOMEN; WITHOUT CONTRAST MATERIAL, FOLLOWED BY CONTRAST MATERIAL(S) AND FURTHER SECTIONS	74170	CT Abdomen W	

### New Referring Physicians by Created Date

Description: This report displays the information about referring physicians created between the dates specified

Parameters:

- Date Range

First Name	Middle Name /Initial	Last Name	Suffix	NPI	ID	Date Created	Active Person	Address 1	Address 2	City	State	Phone	Fax	Active Address
Ronald		Regan		23243234	1286	10-24-2013	Y	2323 North Street		Northfolk	VA	(216) 452-2461		Y
Darcy		Noye	Sr	84712547	1288	11-28-2013	Y	154 Workdown Dr.	1	Fremont	CA		4899994545	Y
Darcy		Noye	Sr	84712547	1288	11-28-2013	Y	951 Home St.		Fremont	CA		4899994545	Y
Unknown		Patient Portal	MD		1294	11-19-2015	Y	456 May St		Livermore	California	(656) 655-6566		Y
Stephen		Gaudet			1297	11-24-2015	Y	<No address on file>						
<b>Total Count</b>														<b>5</b>

New Referring Physicians by Created Date  
 Confidential and Proprietary  
 Page 1 of 1  
 Date Run: 11-07-2016 7:55 PM  
 RRIS-RSNAFormsGuestUser

### Procedure Codes

Description: This report shows all active procedures matched to billing codes (CPT), based on Practice and Sites.

Parameters:

- Practice and Site(s)

### Procedure Codes

Report Description: This report shows all active procedures matched to billing codes (CPT), based on Practice and Sites.

Procedure Code	Description	Billing Code	Body Part	Group Code	Modality
<b>Practice: Treedale Radiology</b>					
☐ <b>Site: Cedarville</b>					
Distinct Procedure Codes: 491			Distinct Billing Codes: 257		
CT13	CT 3 Phase	74170		CT	CT1CE
74160	CT Abdomen W	74160	Abdomen	CT	CT1CE
CT66	CT Abdomen W & CTA Chest w/Wo	71275	Abdomen	CT	CT1CE
CT9	CT Abdomen W & Pelvis W	74160	Abdomen	CT	CT1CE
CT54	CT Abdomen W & Pelvis W & CTA Chest W/Wo	71275	Abdomen	CT	CT1CE
74170	CT Abdomen W/Wo	74170	Abdomen	CT	CT1CE
CT23	CT Abdomen W/Wo & Pelvis W	74170	Abdomen	CT	CT1CE
CT10	CT Abdomen W/Wo & Pelvis W/Wo	74170	Abdomen	CT	CT1CE
74150	CT Abdomen Wo	74150	Abdomen	CT	CT1CE
CT8	CT Abdomen Wo & Pelvis Wo	74150	Abdomen	CT	CT1CE
CT31	CT Ankle Arthrogram	73701	Ankle joint	CT	CT1CE
75571	CT Calcium Scoring	75571	Heart	CT	CT1CE
72126	CT Cervical Spine W	72126	Cervical spine	CT	CT1CE
72127	CT Cervical Spine W/Wo	72127	Cervical spine	CT	CT1CE
72125	CT Cervical Spine Wo	72125	Cervical spine	CT	CT1CE
CT48	CT Cervical Spine Wo & Lumbar Spine Wo	72131	Cervical spine	CT	CT1CE
CT45	CT Cervical Spine Wo & Thoracic Spine Wo	72128	Cervical spine	CT	CT1CE
CT47	CT Cervical Wo, Thoracic Wo, & Lumbar Wo	72131	Cervical spine	CT	CT1CE
71260	CT Chest W	71260	Chest	CT	CT1CE

### Procedure Template Mappings

Description: This report displays the details for the procedure template mappings for the schedule group and templates selected.

Parameters:

- Schedule Group, Templates and whether or not to include exceptions.

#### Procedure Template Mappings

Schedule Group: Cedarville Facility

Template(s): All templates for the selected schedule group.

Include Exceptions: N

Description: This report displays the details for the procedure template mappings for the schedule group and templates selected.

Template Name	Date Range	Monday	Tuesday	Wednesday	Thursday	Friday	Saturday	Sunday
CTProc@Cedar	Default Template	8:00 AM 9:00 PM	8:00 AM 6:00 PM	8:00 AM 9:00 PM	8:00 AM 6:00 PM		8:00 AM 6:00 PM	
						9:00 AM 9:00 PM		9:00 AM 4:00 PM
MAProc@Cedar	Default Template	8:00 AM 6:00 PM	8:00 AM 6:00 PM	8:00 AM 6:00 PM	8:00 AM 6:00 PM		8:00 AM 6:00 PM	
						9:00 AM 2:00 PM		9:00 AM 4:00 PM
MRProc@Cedar	Default Template	8:00 AM 10:00 AM	8:00 AM 10:00 AM	8:00 AM 10:00 AM	8:00 AM 10:00 AM		8:00 AM 10:00 AM	
						9:00 AM 10:00 AM		9:00 AM 4:00 PM
		10:00 AM 11:00 AM	10:00 AM 11:00 AM	10:00 AM 11:00 AM	10:00 AM 11:00 AM	10:00 AM 11:00 AM	10:00 AM 11:00 AM	
		11:00 AM 6:00 PM	11:00 AM 6:00 PM	11:00 AM 6:00 PM	11:00 AM 6:00 PM	11:00 AM 6:00 PM	11:00 AM 6:00 PM	
USProc@Pine	Default Template	8:00 AM 6:00 PM	8:00 AM 6:00 PM	8:00 AM 6:00 PM	8:00 AM 6:00 PM		8:00 AM 6:00 PM	
						9:00 AM 2:00 PM		9:00 AM 4:00 PM
XRProc@Cedar	Default Template	8:00 AM 6:00 PM	8:00 AM 6:00 PM	8:00 AM 6:00 PM	8:00 AM 6:00 PM		8:00 AM 6:00 PM	
						9:00 AM 2:00 PM		9:00 AM 4:00 PM

Procedure Template Mappings

Confidential and Proprietary  
Page 1 of 1

Date Run: 11-08-2016 3:39 PM  
RRIS-RSNA\FORMS\GuestUser

**Procedures without Procedure Group**

Description: This report displays reports that haven't been assigned a procedure group.

Parameters: (No parameters)

**Procedures without Procedure Group**

**Description:** This report displays reports that haven't been assigned a procedure group.

Code	Description	Body Part	Laterality
15119&VRA	VR XR PELVIS		
379554757	MA SCR N MAMMO DIR DIGIT BILAT		
asdf	asf	Clavicle	
MAMDSRB	Digital Bilateral Screening Mammogram/		
MRI	KNEERT		
P22	Procedure 22		
P4	Procedure 4		
PID0004	MRI-BRAIN		
PID0752	MRI CHEST		
RAD	WRR		
US21	US Pelvic Comp and Transvaginal Comp		

### Questionnaire Data Dictionary

Description: This report lists the questions and pertinent information for the selected questionnaire(s).

Parameters:

- Questionnaire dropdown (multi-select)

## Questionnaire Data Dictionary

**Description:** This report lists the questions and pertinent information for the selected questionnaire(s).

Questionnaire	Question Short Name	Question Body	Question Guid	Page #
CT Worksheet				
	DogsName	Dog's name	d2495124-22e3-4dc8-8913-57bca11523ee	1
	question_1	Shellfish allergies?	ff3db556-91dd-479a-9900-cb1acf7f4227	1
	question_3	Types of family members	bb5953df-10c0-4ec6-a860-b756f6e52b11	1
	question_4	Region	7a55c3f9-ffe0-4868-a189-4ee9fe0ff61e	1
	question_5	Birth date	794f8de1-b834-4e27-8016-980bfedb58be	1

Questionnaire Data Dictionary

Confidential and Proprietary  
Page 2 of 14

Date Run: 07-28-2015 03:22 PM  
RMIS\dale




























## Radiologist Setup

Description: This report shows all active Radiologists in the system and indicates if they have assigned NPI and DEA values. MU OK indicator is checked if the physician has both a NPI and DEA number.

Parameters: (No parameters)

## Radiologist Setup

**Report Description:** This report shows all active Radiologists in the system and indicates if they have assigned NPI and DEA values. MU OK indicator is checked if the physician has both a NPI and DEA number.

Radiologist	NPI	DEA	MU OK
Arsenault, Janet	1234567		
Bleignier, Dave	1362		
bob, bob			
Bowles, Tony	2323		
Carter, Jeff	25896		
Connors, Jim	12121		
Connors, Jim	34		
Cunningham, Dave	24323		
Delaney, Bryan	q34		
Dibiasio, Vince	965894		
DigitalFormEditor, Hilary			
Doiron, Denis	333333333		
Eckstein (Rad), Rhonda			
Eckstein, Rhonda			
Gaudet, Stephen			
Graves, Henry	547841		
Gritschke, Kelly	23423234234		
Harding, Clifton			
Harding, Clifton L	878		
Hayes, Tim	q342		
Herman, Jack	2342s		
Hobbs, Henry	65481651		
Hults, Cliff	asdf3wfw		
John, Harry			
Kazem, Imran	5698		

## Referring Physician List

Description: This report displays a list of active referring physicians.

Parameters: (No parameters)

### Referring Physician List

Description: This report displays a list of active referring physicians.

Referring Physician	NPI	UPIN	Licence #	DEA#	Address, Phone & FAX
Adams, John	33445566				323 Little Road Utah, Salt Lake City 12413 Phone: (216) 452-2410 FAX:
Ajala, Clifton	NONPI327				9105 FRANKLIN SQUARE DR SUITE 309 CA, Fremont 94539 Phone: FAX: 4899994545
Armiger, Tom G	NONPI212				1421 S CATON AVE SUITE 203 CA, Fremont 94539 Phone: FAX: 4899994545
Arnold, Bob	NONPI180				201 E UNIVERSITY PKWY CA, Fremont 94539 Phone: (555) 123-4567 FAX: 4899994545
Arnold, Bob	NONPI180				3333 N CALVERT ST SUITE 400 CA, Fremont 94539 Phone: (555) 456-7890 FAX: 4899994545
Arnold, Bob	NONPI180				532 BALTIMORE BLVD SUITE 201 CA, Fremont 94539 Phone: (555) 789-1234 FAX: 4899994545
B16, Henry A	NONPI20166				4512 Pine Drive Apt 105 CA, Fremont 94539 Phone: FAX: 4899994545

### Referring Physician Notification

Description: This report displays new referring physicians recently added to the RIS based on the date the referring was added.

Parameters:

- Date Range based on date physician was added

### Referring Physician Notification

Date Range: 01-01-2013 - 11-08-2016

Description: This report displays new referring physicians recently added to the RIS based on the date the referring was added.

Referring Physician Name	Practice	Address Type	Address	Phone	Email
<b>Gaudet, Stephen</b> Date Added: 11/24/2015	Spec				
	Tree				
<b>Dr Noye, Darcy Sr</b> Date Added: 11/28/2013	Spec	Home	951 Home St. Fremont CA 94539		
		Work	154 Workdown Dr. 1 Fremont CA 94539		
	Tree	Home	951 Home St. Fremont CA 94539		
		Work	154 Workdown Dr. 1 Fremont CA 94539		
<b>Patient Portal, Unknown MD</b> Date Added: 11/19/2015		Work	456 May St Livermore California 94551	(656) 655-6566	
<b>Regan, Ronald</b> Date Added: 10/24/2013		Work	2323 North Street Norfolk VA 23124	(216) 452-2461	

## Report template and macro export

Description: This report is just used to export the information about the diagnostic template and macros.

Parameters: (No parameters)

name	description	active flag	last updated	last updated by user id	owner id	is template flag	macro type code	is normal	is variant	is default
MR LUMBAR SPINE	MR LUMBAR SPINE WITHOUT CONTRAST 72148	Y	11-27-2013 3:37 PM	darcyn	hardingc	Y	Template	N	N	Y
XR Elbow Min 3 Views Normal	XR Elbow	Y	11-27-2013 4:23 PM	stick		Y	Template	Y	N	N
Exam title auto macro contrast	Exam title auto macro contrast	Y	11-23-2015 12:20 AM	saltmarsh		N	ExamTitle	N	Y	N
screening normal fatty	screening normal fatty	Y	11-23-2015 12:55 AM	saltmarsh	saltmarsh	Y	Template	Y	N	N
Screening Mammo Default RSNA	Screening Mammo Default RSNA	Y	11-23-2015 5:19 PM	saltmarsh		Y	Template	N	N	Y
Patient Vitals	Patient Vitals	Y	11-27-2013 3:50 PM	stick		N	Macro	N	Y	N
Exam Title AutoMacro	Exam Title AutoMacro	Y	11-23-2015 12:21 AM	saltmarsh		N	ExamTitle	N	Y	N

### Schedule Group Mappings

Description: This report shows the schedule group mappings for the selected schedule group as it relates to the procedure, rooms and template mappings.

Parameters:

- Schedule groups, Templates and Rooms

#### Schedule Group Mappings

**Schedule Group(s):** Cedarville

**Template(s):** All templates for the selected schedule group(s).

**Room(s):** Cedar: CT1CE, Cedar: MA2CE, Cedar: MR1CE

**Description:** This report shows the schedule group mappings for the selected schedule group as it relates to the procedure, rooms and template mappings.

Site	Procedure Code	Room	Template Name	Template Key	Prep Instruction	Active
						<b># Mappings: 265</b>
Cedarville (Cedar)						
	MR TMJ (70336)	MR1CE	MRProc@Cedar	117		Y
	MR Orbit Face Neck Wo (70540)	MR1CE	MRProc@Cedar	117		Y
	MR Orbit Face Neck W (70542)	MR1CE	MRProc@Cedar	117		Y
	MR Orbit Face Neck WWo (70543)	MR1CE	MRProc@Cedar	117		Y
	MRA Brain Wo (70544)	MR1CE	MRProc@Cedar	117		Y
	MRA Neck Wo (70547)	MR1CE	MRProc@Cedar	117		Y
	MRA Neck W (70548)	MR1CE	MRProc@Cedar	117		Y
	MRA Neck WWo (70549)	MR1CE	MRProc@Cedar	117		Y
	MR Chest Wo (71550)	MR1CE	MRProc@Cedar	117		Y
	MR Chest W (71551)	MR1CE	MRProc@Cedar	117		Y
	MR Chest WWo (71552)	MR1CE	MRProc@Cedar	117		Y
	MRA Chest WWo (71555)	MR1CE	MRProc@Cedar	117		Y
	MR Cervical Spine W (72142)	MR1CE	MRProc@Cedar	117		Y
	MR Thoracic Spine Wo (72145)	MR1CE	MRProc@Cedar	117		Y
	MR Thoracic Spine W (72147)	MR1CE	MRProc@Cedar	117		Y
	MR Lumbar Spine Wo (72148)	MR1CE	MRProc@Cedar	117		Y
	MR Lumbar Spine W (72149)	MR1CE	MRProc@Cedar	117		Y
	MR Thoracic Spine WWo (72157)	MR1CE	MRProc@Cedar	117		Y
	MR Lumbar Spine WWo (72158)	MR1CE	MRProc@Cedar	117		Y
	MRA Spine WWo (72159)	MR1CE	MRProc@Cedar	117		Y
	MR Pelvis Wo (72195)	MR1CE	MRProc@Cedar	117		Y
	MR Pelvis W (72196)	MR1CE	MRProc@Cedar	117		Y
	MR Pelvis WWo (72197)	MR1CE	MRProc@Cedar	117		Y
	MRA Pelvis WWo (72198)	MR1CE	MRProc@Cedar	117		Y
	MR Upper Ext Not Joint Wo (73218)	MR1CE	MRProc@Cedar	117		Y
	MR Upper Ext Not Joint W (73219)	MR1CE	MRProc@Cedar	117		Y

Schedule Group Mappings

Confidential and Proprietary  
Page 1 of 9

Date Run: 11-08-2016 3:35 PM  
RRIS-RSNA\FORMS\QuestUser

### Site Procedure Mappings

Description: This report shows the active mappings for the site, room and procedure for the values selected.

Parameters:

- Practice, Site and Rooms

#### Site Procedure Mappings

Practice: Treeddale Radiology

Site: Cedarville

Room(s): All rooms for the selected site(s).

Description: This report shows the active mappings for the site, room and procedure for the values selected.

Site	Room Code	Room Description	Room Type	Confined Space Flag	Weight Limit	Magnet Strength	Procedure Code	Procedure Description
Cedarville (Cedar)								
	CT1CE	CT 1 at Cedar	CT		300			
							CT13	CT 3 Phase
							CT13	CT 3 Phase
							CT13	CT 3 Phase
							74160	CT Abdomen W
							74160	CT Abdomen W
							74160	CT Abdomen W
							CT66	CT Abdomen W & CTA Chest W/Wo
							CT66	CT Abdomen W & CTA Chest W/Wo
							CT66	CT Abdomen W & CTA Chest W/Wo
							CT9	CT Abdomen W & Pelvis W
							CT9	CT Abdomen W & Pelvis W
							CT9	CT Abdomen W & Pelvis W
							CT54	CT Abdomen W & Pelvis W & CTA Chest W/Wo
							CT54	CT Abdomen W & Pelvis W & CTA Chest W/Wo
							CT54	CT Abdomen W & Pelvis W & CTA Chest W/Wo
							74170	CT Abdomen W/Wo
							74170	CT Abdomen W/Wo
							74170	CT Abdomen W/Wo
							CT23	CT Abdomen W/Wo & Pelvis W
							CT23	CT Abdomen W/Wo & Pelvis W
							CT23	CT Abdomen W/Wo & Pelvis W
							CT10	CT Abdomen W/Wo & Pelvis W/Wo
							CT10	CT Abdomen W/Wo & Pelvis W/Wo
							CT10	CT Abdomen W/Wo & Pelvis W/Wo
							74150	CT Abdomen Wo

### Top RIS Errors By Date

Description: This report displays the Top RIS errors that have occurred for a specified date range.

Parameters:

- Date Range based on error date

### Top RIS Errors By Date

Description: This report displays the Top RIS errors that have occurred for a specified date range

Start date: 10/06/2014

End date: 10/06/2014

Error Description (click text to see detail)	Occurrences
Error calling DisplayData(): System.NullReferenceException: Object reference not set to an instance of an object. at MainApp.Content.Schedule.Sch...	31
Submitting insurance due to change failed: System.ArgumentException: Parent ContentBase not found at System.Windows.Forms.Control.MarshaledInvoke(...	19

### Top RIS Errors By Resource Type

Description: This report displays the Top RIS errors that have occurred for a specified date range and the resource\users selected.

Parameters:

- Date Range based on error date
- Resource types of Editor, Radiologist, Surgeon, Technologist, UM Reviewer
- User(s) filtered by the resource type selected

### Top RIS Errors By Resource Type

**Date Range:** 01-06-2015 - 05-06-2015

**Resource selected:** Radiologist

**Description:** This report displays the Top RIS errors that have occurred for a specified date range and the resource\users selected.

Error Description (click text to see detail)	Occurrences
Validation of a txtFirstName failed with: System.NullReferenceException: Object reference not set to an instance of an object. at UILib.Components.ExtendedErrorProvider.Validate(Control aCtr) in c:\svn\dev\source\ui\UILib\Components\ExtendedError...	33
Validation of a cmbNamePrefix failed with: System.NullReferenceException: Object reference not set to an instance of an object. at UILib.Components.ExtendedErrorProvider.Validate(Control aCtr) in c:\svn\dev\source\ui\UILib\Components\ExtendedError...	19



### Top RIS Errors By UserID

Description: This report displays the Top RIS errors that have occurred for the specified date range and user id.

Parameters:

- Date Range based on error date and User ID

### Top RIS Errors By UserID

**Date Range:** 01-06-2015 - 05-06-2015

**Description:** This report displays the Top RIS errors that have occurred for the specified date range and user id

Error Description (click text to see detail)	Occurrences
Validation of a txtFirstName failed with: System.NullReferenceException: Object reference not set to an instance of an object. at UI.Lib.Components.ExtendedErrorProvider.Validate(Control aCtr) in c:\svn\dev\source\ui\UI\Lib\Components\ExtendedError...	33
Validation of a cmbNamePrefix failed with: System.NullReferenceException: Object reference not set to an instance of an object. at UI.Lib.Components.ExtendedErrorProvider.Validate(Control aCtr) in c:\svn\dev\source\ui\UI\Lib\Components\ExtendedError...	19

### Top RIS Errors Since Yesterday

Description: This report displays the Top RIS errors that have occurred since yesterday.

Parameters: (No parameters)

### Top RIS Errors Since Yesterday

Description: This report displays the Top RIS errors that have occurred since yesterday

Error Description (click text to see detail)	Occurrences
Error calling DisplayData(): System.NullReferenceException: Object reference not set to an instance of an object. at MainApp.Content.Schedule.Sch...	35
Submitting insurance due to change failed: System.ArgumentException: Parent ContentBase not found at System.Windows.Forms.Control.MarshaledInvoke(...	20

### User Group Permissions

Description: This report displays the associated permissions and access levels assigned to the selected user group.

Parameters:

- User Group, and Access Level(s)

### User Group Permissions

User Group(s): Add Insurance

Access Levels(s): All

Description: This report displays the associated permissions and access levels assigned to the selected user group.

User Group	Access Code	Access Description	Access Level
Add Insurance	Clinical.ChangeSupervisingRadiologist	Access to WL context menu for 'Change Supervising Radiologist'	Full
	Clinical.MarkForAddendum	Access to WL context menu for 'Request Addendum' and 'Cancel Addendum Request'	Full
	Clinical.PatientMergeMixer	Patient merge screen access	Full
	Config.Configure	Controls Access to the Configure menu option from Administration menu	Full
	Config.LookupEditor.Carrier	lookup table Carrier for Carrier	Full
	DataNugget.Attachments	Access to see the nugget for Documents	Full
	Menu.Administration	Controls access to the Administration Menu	Full
	WL.UnknownReconciliation	WL access for Unknown Referring & Carrier	Full
<b>Add Insurance Permissions Count</b>			<b>8</b>
<b>Permissions Count</b>			<b>8</b>

### User List

Description: This report lists the users based on active or inactive parameter as well as the practices those users have been assigned.

Parameters:

- Active\Inactive (Users)

### User List

**Description:** This report lists the users based on active or inactive parameter. Including indicators to identify which practices and schedule groups assigned.

User Id	Last name	First name	Practices				Schedule Groups						
					Rainbow Radiology	Spectrum Radiology	Treesdale Radiology			<Not associated to a schedule group>	Greenville Facility	SpectrumSchedGroup	TreesdaleSchedGroup
darcy	Aiken	Darcy					X						X
barnold	Arnold	Bob								X			
janet	Arsenault	Janet					X						X
marc	Arsenault	Marc					X						X
bob	bob	bob	X		X	X	X	X					X
jconners	Connors	Jim				X	X						X
dcunningham	Cunningham	Dave				X	X						X
vince	Dibiasio	Vince				X	X				X		X

### User Permissions

Description: Displays the user groups, and permissions associated to the specified user.

Parameters:

- User ID

### User Permissions

**Description:** Displays the user groups, and permissions associated to the specified user.

User	Groups
(eckstein): Eckstein, Rhonda	Demo Super User
<b>Count: 1</b>	
User Specific Access	
Ability to add/modify delete billing codes (Clinical.BillingCodes.Modify)	
Ability to override utilization review and proceed to schedule. (Clinical.UtilizationManagement.ScheduleOverride)	
Ability to release a tentatively signed report that was signed by somebody else (Reporting.ReleaseOthersTentativeSignedReports)	
Access ePrescription (eRx) Administrator Role (MU.ePrescribeAdministrator)	
Access ePrescription (eRx) Radiologist Role (MU.ePrescribeRadiologist)	
Access ePrescription (eRx) Staff Role (MU.ePrescribeStaff)	
Access for the Activity WL under the Administration menu (WL.Activity)	
Access in the report editing window to export the current reports audio and text to an external location. (Reporting.ExportAudioAndReport)	
Access in the report editing window to import external audio and reports into current report window (Reporting.ImportAudioAndReport)	
Access MU controls on a form (radio buttons and grids) (MU.MUEnabledControls)	
Access Quick Add Edit Dialog (on mu tab) (MU.QuickAddEddDialog)	
Access to 2D Barcode workflow for Reception (Clinical.ReceptionBarcodeWorkflow)	
Access to add CTRN notification (clinical.CTRNNotification)	
Access to 'All Edits Pending QA' menu option from Editor main menu (WL.MT.EditsPendingQA)	
Access to 'All Problem' menu option from Editor main menu (WL.MT.ReportDraftedProblem)	
Access to 'All Problem' menu option from Radiologist main menu (WL.Rad.TentativelySigned)	

### Users in Usergroups

Description: This report displays the list of users assigned to the selected user group(s).

Parameters:

- User group(s), active user flag

### Users in Usergroups

User Group(s): Scheduling (Scheduling)

Description: This report displays the list of users assigned to the selected user group(s).

User id	Name	Active User
<b>Group: Scheduling (Scheduling)</b>		
darcy	Aiken, Darcy j	Y
marc	Arsenault, Marc	Y
jconners	Connors, Jim	Y
dcunningham	Cunningham, Dave	Y
sgaudet	Gaudet, Stephen	Y
chults	Hults, Cliff	Y
skoepfel	Koeppel, Seth	Y
darcyn	Noye, Darcy Thane	Y
darcynrad	Noye, Darcy, Sr	Y
ian	Power, Ian	Y
Fpresent	Present, Florence	Y
reyesr	Reys, Rey	Y
dtaylor	Taylor, Dan	Y
schedtest	Test, Scheduler	Y
schedtest1	Test, Scheduler 1	Y
jtollack	Tollack, James	Y
andrew	Waite, Andrew R	Y
rziegler	Ziegler, Randy	Y
<b>Total Count</b>		<b>18</b>

## XR-29 Compliance

Description: This report displays the modalities for the specified practices, sites and modality types with the compliance date for each

Parameters:

- Site(s) and Practice(s), Modality Type(s), Active Flag Value(s)

### XR-29 Compliance

**Practice(s):** Treesdale Radiology

**Site(s):** All sites for the selected practice(s).

**Modality Type(s):** CT

**Active Flag values:** Y

**Description:** This report display all the active modalities for the specified practices, sites and modality types, returning the list and XR-29 Compliance Date for each.

Practice \ Site	Modality Code	Modality Description	AE Title	Active	XR-29 Compliance Date
<b>Treesdale Radiology (Tree)</b>					<b>3</b>
<b>Cedarville (Cedar)</b>					<b>1</b>
	CT1CE	CT 1 at Cedar	CT1CE	Y	08-08-2016
<b>Mapleton (Maple)</b>					<b>1</b>
	CT1MP	CT 1 at Maple	CT1MP	Y	
<b>Pine Grove (Pine)</b>					<b>1</b>
	CT1PG	CT 1 at Pine	CT1PG	Y	11-01-2016

XR-29 Compliance

Confidential and Proprietary  
Page 1 of 1

Date Run: 11-08-2016 02:12 PM  
RRIS-RSNA\FORMS\GuestUser



## Management Reports \ Mammography

### Mammo Annual Radiologist Summary

Description: This report displays annual Radiologist exams read summary. All Exams are Signed exams.

Parameters:

- Practice, Site(s), Signing Radiologists, Procedure Codes and Reporting Year

### Mammo Annual Radiologist Summary

**Practice:** Treesdale Radiology

**Site(s):** All Sites

**Reporting Year:** 2015

**Signing Radiologist:** Noye, Darcy (darcyn), Saltmarsh, Hilary (saltmarsh)

**Procedure Codes:** 76641, 76641B, 76642, 76642B, G0202, G0202IMP, G0202UIMP, G0202UNI, G0202WC, G0204, G0204/US, G0204CB, G0204CB/US, G0204IMP, G0206, G0206/US, G0206CB, G0206CB/US, G0206IMP, MATO1, MATO2

**Description:** This report displays annual Radiologist exams read summary. All Exams are Signed exams.

Radiologist Reporting Totals for 2015	Accession	Patient Name	Signed Date	BI-RADS®
<input type="checkbox"/> Noye, Darcy Thane	1			
<input type="checkbox"/> November	1			
<input type="checkbox"/> MA Digital Mammo Screen Bilat [G0202] - Breast - Bilateral	1			
	1024217	Test Portal, middletest	11-26-2015	0
<input type="checkbox"/> Saltmarsh, Hilary	11			
<b>Report Total:</b>				<b>12</b>

Mammo Annual Radiologist Summary

Confidential and Proprietary  
Page 1 of 1

Date Run: 11-08-2016 5:31 PM  
RRIS-RSNAVFormsGuestUser



### Mammo Annual Technologist Summary

Description: This report displays annual Technologist exams performed summary.

Parameters:

- Practice, Site(s), Signing Radiologists, Procedure Codes and Reporting Year

### Mammo Annual Technologist Summary

**Practice:** Treesdale Radiology

**Site(s):** Cedarville, Mapleton, Pine Grove, Outside Read1 (outside read)

**Performed Year:** 2015

**Performed by Technologist:** Aiken, Darcy (darcy), Arsenault, Janet (janet), bob, bob (bob), Conners, Jim (jconners), Cunningham, Dave (dcunningham), Dibiasio, Vince (vince), DigitalFormEditor, Hilary (dfeditor), Doiron, Denis (denis), Eckstein, Rhonda (eckstein), Gaudet, Stephen (sgaudet), Harding, Clifton (clifton), Hults, Cliff (chults), Koepfel, Seth (skoepfel), Mills, Terry (terry), Noye, Darcy (darcyn), Noye, Darcy (darcynrad), Powell, Stick (stick), Present, Florence (fpresent), Reys, Rey (reysr), Saltmarsh, Hilary (saltmarsh), stickgod, stickgod (stickgod), Taylor, Dan (dtaylor), Test User, System Admin (administrator), Test, TechRad (techrad), Tollack, James (jtollack), Waite, Andrew (andrew), Ziegler, Randy (rzeigler)

**Procedure Codes:** 76641, 76641B, 76642, 76642B, G0202, G0202IMP, G0202UIMP, G0202UNI, G0202VVC, G0204, G0204/US, G0204CB, G0204CB/US, G0204IMP, G0206, G0206/US, G0206CB, G0206CB/US, G0206IMP, MATO1, MATO2

**Description:** This report displays annual Technologist exams performed summary.

Technologist Performed Totals for 2015	Accession	Patient Name	Exam Status	BI-RADS®	Performed Date
<input type="checkbox"/> Arsenault, Janet	1				
<input type="checkbox"/> Conners, Jim	7				
<input type="checkbox"/> April	7				
<input type="checkbox"/> MA Digital Mammo Diag Bilat [G0204] - Breast - Bilateral	7				
	1023211	Zulu Average,	Signed1	4	04-06-2015
	1027770	Adams Average, Demo	Signed1	4	04-06-2015
	1027780	O'Brian Hillary, Demo	Signed1	4	04-06-2015
	1027790	O'Conner Average, Demo	Signed1	4	04-06-2015
	1027800	O'Dell Jayne, Demo	Signed1	4	04-06-2015
	1027810	O'Neil Average, Demo	Signed1	4	04-06-2015
	1027820	ד'אור דוד, Demo	Signed1	4	04-06-2015
<input type="checkbox"/> Present, Florence	1				
<input type="checkbox"/> Saltmarsh, Hilary	2				
<input type="checkbox"/> Test User, System Admin	1				
<b>Report Total:</b>					<b>12</b>



### Mammo Back on Track

Description: The intention of this report is to identify any patient in the system that has order type of "CallBack" or "ShortTermFollowup" where the Mammo orders have come for an original study with a diagnosis of "BI-RADS® 1 or 2".

Parameters:

- Practice, Site(s)

### Mammo Back on Track

Practice: Tressdale

Site(s) CEDARVILLE

Date Range: 10-01-2016 - 11-08-2016

Description: The intention of this report is to identify any patient in the system that has order type of "CallBack" or "ShortTermFollowup" where the mammo orders have come for an original study with a diagnosis of "BIRAD 1 or 2". This is based on the Practice, Site(s) and date range (Target Date) selected.

Patient	MRN	Study Status	Requested Date	Target Date	Order Status	Original Accession	Recommendation	Birad
CANTON, SUE	01036	Signed1	07-01-2016	10-29-2016	ShortTermFollowUp	100327	4 months	2-Benign finding
RODRIGUEZ, RHONDA	010782	Signed1	04-25-2016	10-22-2016	ShortTermFollowUp	10028	6 months	2-Benign finding
JAMES, JOHN	01082	Signed1	05-05-2016	11-01-2016	ShortTermFollowUp	100315	6 months	2-Benign finding
JONES, SARA	1144	Signed1	04-06-2016	10-03-2016	ShortTermFollowUp	10030	6 months	2-Benign finding
O'DELL, JAYNE	13306	Signed1	07-25-2016	10-23-2016	ShortTermFollowUp	10036	3 months	2-Benign finding
O'CONNOR, WILLIAM	30621	Signed1	04-28-2016	10-25-2016	ShortTermFollowUp	100289	6 months	2-Benign finding
MAYA, JUAN	3338	Signed1	04-29-2016	10-26-2016	ShortTermFollowUp	10031	6 months	2-Benign finding
<b>Total</b>								<b>7</b>

Mammo Back on Track

Confidential and Proprietary  
Page 1 of 1

Date Run: 11-08-2016 11:03 PM  
RADNETEcksteinR

### Mammo Document Distribution Monitoring

Description: This report displays mammography documents queued in the Document Distribution Worklist. Items older than 2 days will be red. Layletters have to be distributed within 30 days of exam.

Parameters:

- Practice, Site(s)

### Mammo Document Distribution Monitoring

**Practice:** Treesdale Radiology

**Site(s):** Cedarville, Mapleton, Pine Grove, Outside Read1 (outside read)

**Description:** This report displays mammography documents queued in the Document Distribution Worklist. Items older than 2 days will be red. Layletters have to be distributed within 30 days of exam.

Document Type / Status	ACC#	Recipient:	BI-RADS®	Date Created	Time in DDE Queue
<input type="checkbox"/> LAYLETTER		3			
<input type="checkbox"/> Hold		3			
	1012647	Chan, Anne	0	10-25-2012	1475
	1014612	sally, sally	4	07-05-2013	1222
	1014833	Olivens, Jenny	0	07-05-2013	1222
<b>Report Total:</b>					<b>3</b>

Mammo Document Distribution Monitoring

Confidential and Proprietary  
Page 1 of 1

Date Run: 11-08-2016 5:25 PM  
RRIS-RSNA\FORMS\GuestUser

### Mammo Ignored Reminder List

Description: This report displays the list of patients that have Ignored Mammography Reminder letters.

Parameters:

- Practice, Site(s), and Referring Physician

### Mammo Ignored Reminder List

Practice: Treesdale

Site(s): Cedarville

Referring Physician: BOB ARNOLD, MD, FRED FREDERICKSON, MD

Description: This report displays the list of patients that have Ignored Mammography Reminder letters.

Referring Name / Reminder Ignored	Target Date	Patient ID	Patient Name	Procedures	Last Letter Sent	Last Letter Date
<input type="checkbox"/> BOB ARNOLD, MD	1					
<input type="checkbox"/> Invitation						
<input type="checkbox"/> FRED FREDERICKSON, MD	1					
<input type="checkbox"/> Invitation						
	06-10-2016	00902678BI		Breast Imaging to be Determined [BID]	Reminder Letter3	09-09-2016
<b>Report Total:</b>						<b>2</b>

Mammo Ignored Reminder List

Confidential and Proprietary  
Page 1 of 1

Date Run: 11-10-2016 4:49 PM  
RADNETE:ksteinR

### Mammo Orphaned Orders

Description: This report will identify patients that have possible orphaned mammography orders. The report looks for any studies in the system that are past the outstanding Mammography Order date.

Parameters:

- Practice and Site(s)

### Mammo Orphaned Orders

Practice: Treesdale

Description: This report will identify patients that have possible orphaned mammography orders. The report looks for any studies in the system that are past the outstanding Mammography Order date.

Practice / Site	Patient ID	Orphan Filler Order Number	Orphan Order Request Date	Orphan Order Procedure
<input type="checkbox"/> Treesdale		1		
<input type="checkbox"/> Cedarville		1		
	7793	2004	05-27-2015	US Pelvic TA and TV [US500] - Pelvis   MA Digital Screening Mammo [G0202] - Breast - B
<b>Total Patients:</b>		<b>1</b>		

Mammo Orphaned Orders

Confidential and Proprietary  
Page 1 of 1

Date Run: 11-10-2016 10:35 AM  
RADNET\EcksteinR

I

### Mammo Radiologist Summary

Description: This report displays Radiologist exams read summary for a specified period. All Exams are Signed exams.

Parameters:

- Practice, Site(s), Signing Radiologists, Procedure Codes and Date Range

### Mammo Radiologist Summary

**Practice:** Treesdale Radiology

**Site(s):** Cedarville, Mapleton, Pine Grove, Outside Read1 (outside read)

**Date Range:** 01-01-2015 - 11-08-2016

**Procedure Codes:** 76641, 76641B, 76642, 76642B, G0202, G0202IMP, G0202UIIMP, G0202UNI, G0202WC, G0204, G0204/US, G0204CB, G0204CB/US, G0204IMP, G0206, G0206/US, G0206CB, G0206CB/US, G0206IMP, MATO1, MATO2

**Description:** This report displays Radiologist exams read summary for a specified period. All Exams are Signed exams.

Radiologist Reporting Totals:	Accession	Patient Name	Signed Date	BI-RADS®
<input type="checkbox"/> Noye, Darcy Thane	1			
<input checked="" type="checkbox"/> November 2015	1			
<input checked="" type="checkbox"/> Saltmarsh (Rad), Hilary	1			
<input type="checkbox"/> Saltmarsh, Hilary	15			
<input type="checkbox"/> November 2015	9			
<input type="checkbox"/> MA Digital Mammo Diag Bilat [G0204] - Breast - Bilateral	8			
	1023211	Zulu Average,	11-25-2015	4
	1027770	Adams Average, Demo	11-25-2015	4
	1027780	O'Brian Hillary, Demo	11-25-2015	4
	1027790	O'Conner Average, Demo	11-25-2015	4
	1027800	O'Dell Jayne, Demo	11-25-2015	4
	1027810	O'Neil Average, Demo	11-25-2015	4
	1027820	דוד אורן, Demo	11-25-2015	4
	1024375	Flotz Gloria,	11-27-2015	3
<input checked="" type="checkbox"/> MA Digital Mammo Screen Unilat [G0202UNI] - Breast - Left	1			
<input checked="" type="checkbox"/> December 2015	2			
<input checked="" type="checkbox"/> February 2016	1			
<input checked="" type="checkbox"/> March 2016	1			
<input checked="" type="checkbox"/> August 2016	1			
<input checked="" type="checkbox"/> September 2016	1			
<b>Report Total:</b>	<b>17</b>			

### Mammo Technologist Summary

Description: This report displays Technologist exams performed summary for a specified period.

Parameters:

- Practice, Site(s), Performed by Technologists, Procedure Codes and Date Range

### Mammo Technologist Summary

**Practice:** Treesdale Radiology

**Site(s):** All sites for the selected practice(s).

**Date Range:** 01-01-2015 - 11-08-2016

**Procedure Codes:** All active BI-RADS® enabled procedures.

**Description:** This report displays Technologist exams performed summary for a specified period.

Technologist Performed Totals:	Accession	Patient Name	Exam Status	BI-RADS®	Performed Date
<input type="checkbox"/> Saltmarsh, Hilary		2			
<input type="checkbox"/> December 2015		2			
<input type="checkbox"/> MA Digital Mammo Screen Bilat [G0202] - Breast - Bilateral		2			
	1025281	Olivens Jenny,	Signed1	0	12-16-2015
	1025409	Olivens Jenny,	Signed1	3	12-17-2015
<b>Report Total:</b>					<b>2</b>

Mammo Technologist Summary

Confidential and Proprietary  
Page 1 of 1

Date Run: 11-08-2016 5:21 PM  
RRIS-RSNA\FORMS\GuestUser

### Mammography Missing Data

Description: This report displays the Mammo studies (based on procedure codes that have the BI-RADS® flag enabled) for the practice(s), site(s) and date selected (signed date).

Three conditions will get the exam on this report, empty BI-RADS®, empty Recommendation or if the BI-RADS® density code is empty but the procedure is configured to be required.

Parameters:

- Practice(s), Site(s) and Date Range

### Mammography Missing Data

Practice: Treedale

Site(s): Cedarville

Date Range: 10-01-2014 - 11-10-2016

Description: This report displays the mammo studies (based on procedure codes that have the BI-RADS® flag enabled) for the practice(s), site(s) and date selected (signed date).

Three conditions will get the exam on this report, empty BI-RADS®, empty Recommendation or if the BI-RADS® density code is empty but the procedure is configured to be required.

Acc #	Patient	BI-RADS®	Recom. Code	Density Code	Invalid Recom.	DOS	Dictated	Dictated By
18779	STRICKLAND, RHONDA				N	03-25-2014	10-06-2014	Eckstein MD, Rhonda
MRN: 4335	MAM, DIAG BIL & CAD & BREAST US [BORGMA200]							
20155	DEMOTON, JAYNE				N	08-21-2014	11-06-2014	Saltmarsh MD, Hilary
MRN: 01212	MAM, SCREEN & CAD [BORGMA156]							
100063	BLOOM, APRIL	6	B		Y	10-13-2014	10-14-2014	Powell MD, Stephen
MRN: 3752	MR Breast W and or WO Bilateral W Cad [77059] - Breast - B							

Mammography Missing Data

Confidential and Proprietary  
Page 1 of 1

Date Run: 11-10-2016 5:02 PM  
RADNETEcksteinR



### MQSA General Exam List

Description: This report displays the general mammography exam list for selected sites/practice. All Exams are Signed exams.

Parameters:

- Practice, Site(s), Date Range, Signing Radiologists, Procedure Codes, and Biopsy Inclusion

### MQSA General Exam List

**Practice:** Treesdale Radiology

**Site(s):** Cedarville, Mapleton, Pine Grove, Outside Read1 (outside read)

**Date Range:** 02-01-2015 - 11-30-2015

**Signing Radiologist:** Saltmarsh, Hilary (saltmarsh)

**Procedure Codes:** 76641, 76641B, 76642, 76642B, G0202, G0202IMP, G0202UIMP, G0202UNI, G0202WC, G0204, G0204/US, G0204CB, G0204CB/US, G0204IMP, G0206, G0206/US, G0206CB, G0206CB/US, G0206IMP, MATO1, MATO2

**Description:** This report displays the general mammography exam list for selected sites/practice. All Exams are Signed exams.

BI-RADS® / Recommendation	Accession	Patient Name	Signed Date	Signing Radiologist	Procedure	Biopsy	Biopsy Date	Biopsy Result	
<input checked="" type="checkbox"/> BI-RAD® 3	2								
<input type="checkbox"/> BI-RAD® 4	7								
<input type="checkbox"/> Biopsy should be considered	7								
	1023211	Zulu Average,	11-25-2015	Saltmarsh, Hilary	MA Digital Mammo Diag Bilat [G0204] - Breast - Bilateral	N			
	1027770	Adams Average, Demo	11-25-2015	Saltmarsh, Hilary	MA Digital Mammo Diag Bilat [G0204] - Breast - Bilateral	N			
	1027780	O'Brian Hillary, Demo	11-25-2015	Saltmarsh, Hilary	MA Digital Mammo Diag Bilat [G0204] - Breast - Bilateral	N			
	1027790	O'Conner Average, Demo	11-25-2015	Saltmarsh, Hilary	MA Digital Mammo Diag Bilat [G0204] - Breast - Bilateral	N			
	1027800	O'Dell Jayne, Demo	11-25-2015	Saltmarsh, Hilary	MA Digital Mammo Diag Bilat [G0204] - Breast - Bilateral	N			
	1027810	O'Neil Average, Demo	11-25-2015	Saltmarsh, Hilary	MA Digital Mammo Diag Bilat [G0204] - Breast - Bilateral	N			
<b>Report Total:</b>									<b>9</b>

### MQSA Layletter Audit

Description: This report displays the list of mammography exams that LAYLETTERS were not distributed to via Document Distribution Engine within 30 days, by date, location. All Exams are Signed exams.

Parameters:

- Practice, Site(s), and Date Range

### MQSA Layletter Audit

Practice: Treesdale

Site(s): Cedarville

Date Range: 07-01-2015 - 11-10-2016

Description: This report displays the list of mammography exams that LAYLETTERS were not distributed to via Document Distribution Engine within 30 days, by date, location. All Exams are Signed exams.

BI-RADS®	Patient ID	Patient Name	Accession	Signed Date	DDE Status	DDE Last Updated Date	DDE Last Updated User
<input checked="" type="checkbox"/> BI-RAD® 0							
<input checked="" type="checkbox"/> BI-RAD® 1							
<input checked="" type="checkbox"/> BI-RAD® 2							
<input checked="" type="checkbox"/> BI-RAD® 3							
<input checked="" type="checkbox"/> BI-RAD® 4							
<input checked="" type="checkbox"/> BI-RAD® 5	32145	Madden, Rhonda	5432333	07-13-2015	Completed	07-13-2015	eckstein
<b>Total:</b>							<b>344</b>

MQSA Layletter Audit

Confidential and Proprietary  
Page 1 of 1

Date Run: 11-10-2016 2:24 PM  
RADNET\EcksteinR

### MQSA List of Exams Performed by Technologist

Description: This report displays the list of mammography exams performed by Technologist.

Parameters:

- Practice, Site(s), Date Range, Performing Technologist and Procedure Codes

### MQSA List of Exams Performed by Technologist

**Practice:** Treesdale Radiology

**Site(s):** Cedarville, Mapleton, Pine Grove, Outside Read1 (outside read)

**Date Range:** 01-01-2014 - 11-08-2016

**Performed by Technologist:** Saltmarsh, Hilary (saltmarsh)

**Procedure Codes:** 76641, 76641B, 76642, 76642B, G0202, G0202IMP, G0202UIMP, G0202UNI, G0202WC, G0204, G0204/US, G0204CB, G0204CB/US, G0204IMP, G0206, G0206/US, G0206CB, G0206CB/US, G0206IMP, MATO1, MATO2

**Description:** This report displays the list of mammography exams performed by Technologist.

Technologist:	Procedure	Accession	Patient Name	Exam Status	BI-RADS®	Performed Date
☐ Saltmarsh, Hilary		2				
	☐ MA Digital Mammo Screen Bilat [G0202] - Breast - Bilateral	2				
		1025281	Olivens Jenny,	Signed1	0	12-16-2015
		1025409	Olivens Jenny,	Signed1	3	12-17-2015
<b>Report Total:</b>						<b>2</b>

MQSA List of Exams Performed by Technologist

Confidential and Proprietary  
Page 1 of 1

Date Run: 11-08-2016 4:57 PM  
RRIS-RSNA\FORMS\GuestUser

### MQSA List of False Positive Exams

Description: This report displays the general mammography exam list for selected sites/practice. All Exams are Signed exams.

Parameters:

- Practice, Site(s), Date Range, Signing Radiologists, Procedure Codes and Biopsy Inclusion

### MQSA List of False Positive Exams

Practice: Treesdale Radiology

Site(s): Cedarville, Mapleton, Pine Grove, Outside Read1 (outside read)

Date Range: 01-01-2015 - 11-01-2016

Signing Radiologist: Saltmarsh, Hilary (saltmarsh)

Procedure Codes: 76641, 76641B, 76642, 76642B, G0202, G0202IMP, G0202UIMP, G0202UNI, G0202WC, G0204, G0204/US, G0204CB, G0204CB/US, G0204IMP, G0206, G0206/US, G0206CB, G0206CB/US, G0206IMP, MATO1, MATO2

Description: This report displays the general mammography exam list for selected sites/practice. All Exams are Signed exams.

BI-RADS® / Recommendation	Accession	Patient Name	Signed Date	Signing Radiologist	Procedure	Biopsy	Biopsy Date	Biopsy Result	Biopsy Technique	Pathology Result
<input type="checkbox"/> BI-RAD® 4		1								
<input type="checkbox"/> Biopsy should be considered		1								
	1025446	Murcia Hilary,	03-11-2016	Saltmarsh, Hilary	MA Digital Mammo Diag Bilat [G0204] - Breast - Bilateral	Y	02-28-2016	Benign	Stereotactic Core Bi	Ductal carcinoma in-situ (DCIS), Axillary nodal metastases
<b>Report Total:</b>										<b>1</b>

MQSA List of False Positive Exams

Confidential and Proprietary  
Page 1 of 1

Date Run: 11-08-2016 4:50 PM  
RRIS-RSNAVFormsGuestUser

### MQSA List of False Negative Exams

Description: This report displays the general mammography exam list for selected sites/practice. All Exams are Signed exams.

Parameters:

- Practice, Site(s), Date Range, Signing Radiologists, Procedure Codes and Biopsy Inclusion

### MQSA List of False Negative Exams

Practice: Treesdale Radiology

Site(s): Cedarville, Mapleton, Pine Grove, Outside Read1 (outside read)

Date Range: 01-01-2015 - 11-08-2016

Signing Radiologist: Saltmarsh, Hilary (saltmarsh)

Procedure Codes: 76641, 76641B, 76642, 76642B, G0202, G0202IMP, G0202UIMP, G0202UNI, G0202WVC, G0204, G0204/US, G0204CB, G0204CB/US, G0204IMP, G0206, G0206/US, G0206CB, G0206CB/US, G0206IMP, MATO1, MATO2

Description: This report displays the general mammography exam list for selected sites/practice. All Exams are Signed exams.

BI-RADS® / Recommendation	Accession	Patient Name	Signed Date	Signing Radiologist	Procedure	Biopsy	Biopsy Date	Biopsy Result
<input checked="" type="checkbox"/> BI-RAD® 0	2							
<input checked="" type="checkbox"/> BI-RAD® 1	1							
<input checked="" type="checkbox"/> BI-RAD® 2	1							
<input checked="" type="checkbox"/> BI-RAD® 3	3							
<input type="checkbox"/> BI-RAD® 4	8							
<input type="checkbox"/> Biopsy should be considered	8							
	1023211	Zulu Average,	11-25-2015	Saltmarsh, Hilary	MA Digital Mammo Diag Bilat [G0204] - Breast - Bilateral	N		
	1027770	Adams Average, Demo	11-25-2015	Saltmarsh, Hilary	MA Digital Mammo Diag Bilat [G0204] - Breast - Bilateral	N		
	1027780	O'Brian Hillary, Demo	11-25-2015	Saltmarsh, Hilary	MA Digital Mammo Diag Bilat [G0204] - Breast - Bilateral	N		
	1027790	O'Conner Average, Demo	11-25-2015	Saltmarsh, Hilary	MA Digital Mammo Diag Bilat [G0204] - Breast - Bilateral	N		
	1027800	O'Dell Jayne, Demo	11-25-2015	Saltmarsh, Hilary	MA Digital Mammo Diag Bilat [G0204] - Breast - Bilateral	N		
	1027810	O'Neil Average, Demo	11-25-2015	Saltmarsh, Hilary	MA Digital Mammo Diag Bilat [G0204] - Breast - Bilateral	N		
	1027820	ד"ר אור דוד, Demo	11-25-2015	Saltmarsh, Hilary	MA Digital Mammo Diag Bilat [G0204] - Breast - Bilateral	N		
	1025446	Murcia Hilary,	03-11-2016	Saltmarsh, Hilary	MA Digital Mammo Diag Bilat [G0204] - Breast - Bilateral	Y	02-28- 2016	Benign
<b>Report Total:</b>	<b>15</b>							

**MQSA List of Positive Exam Followups (Stand Alone System)**

Description: This report displays the list of Follow-up (Orders / Exams) for Positive Exams, by date, location. All Exams are Signed exams. This report is intended for use with RIS Mammography Stand Alone module.

Parameters:

- Practice, Site(s) and Date Range

**MQSA List of Positive Exam Followups (Stand Alone System)**

Practice: Treesdale Radiology

Site(s): Cedarville, Mapleton, Pine Grove, Outside Read1 (outside read)

Date Range: 01-01-2015 - 11-08-2016

Description: This report displays the list of Follow-up (Orders / Exams) for Positive Exams, by date, location. All Exams are Signed exams. This report is intended for use with RIS Mammography Stand Alone module.

Original ACC# - Patient Name	Next Accession	Next Accession Status	Next Procedure performed
<input checked="" type="checkbox"/> BI-RAD@ 0	3		
<input type="checkbox"/> BI-RAD@ 4	8		
<input checked="" type="checkbox"/> 1023211 - Zulu Average,			
<input checked="" type="checkbox"/> 1025446 - Murcia Hilary,			
<input type="checkbox"/> 1027770 - Adams Average, Demo			
	1027769	Signed1	MR Cervical Spine WWVo [72156] - Cervical spine
<input type="checkbox"/> 1027780 - O'Brian Hillary, Demo			
	1027779	Signed1	MR Cervical Spine WWVo [72156] - Cervical spine
<input type="checkbox"/> 1027790 - O'Conner Average, Demo			
	1027789	Signed1	MR Cervical Spine WWVo [72156] - Cervical spine
<input checked="" type="checkbox"/> 1027800 - O'Dell Jayne, Demo			
<input checked="" type="checkbox"/> 1027810 - O'Neil Average, Demo			
<input checked="" type="checkbox"/> 1027820 - ד"ר אורן דוד, Demo			
<b>Report Total:</b>			<b>11</b>

MQSA List of Positive Exam Followups (Stand Alone System)

Confidential and Proprietary  
Page 1 of 1

Date Run: 11-08-2016 4:48 PM  
RRIS-RSNA\FORMS\GuestUser



### MQSA List of Positive Exam Followups

Description: This report displays the list of Follow-up (Orders / Exams) for Positive Exams, by date, location. All Exams are Signed exams.

Parameters:

- Practice, Site(s) and Date Range

### MQSA List of Positive Exam Followups

Practice: Treesdale Radiology

Site(s): Cedarville, Mapleton, Pine Grove, Outside Read1 (outside read)

Date Range: 01-01-2015 - 11-08-2016

Description: This report displays the list of Follow-up (Orders / Exams) for Positive Exams, by date, location. All Exams are Signed exams.

Original ACC# - Patient Name	Requested Date	Next Order Status	Requested by	Next Accession	Next Accession Status	Next Procedure performed
<input checked="" type="checkbox"/> BI-RAD@ 0	3					
<input type="checkbox"/> BI-RAD@ 4	8					
<input checked="" type="checkbox"/> 1023211 - Zulu Average,						
<input checked="" type="checkbox"/> 1027770 - Adams Average, Demo						
<input type="checkbox"/> 1027780 - O'Brian Hillary, Demo						
		No Followup				
<input type="checkbox"/> 1027790 - O'Conner Average, Demo						
		No Followup				
<input type="checkbox"/> 1027800 - O'Dell Jayne, Demo						
		No Followup				
<input type="checkbox"/> 1027810 - O'Neil Average, Demo						
		No Followup				
<input checked="" type="checkbox"/> 1027820 - מורציה הילרי, Demo						
<input checked="" type="checkbox"/> 1025446 - Murcia Hilary,						
<b>Report Total:</b>						<b>11</b>

MQSA List of Positive Exam Followups

Confidential and Proprietary  
Page 1 of 1

Date Run: 11-08-2016 4:46 PM  
RRIS-RSNA\FORMS\GuestUser

### MQSA List of Positive Exams with Outstanding Biopsies

Description: This report displays the list of Positive Mammography Exams (for selected time period, location, radiologist) that do not have pathology result recorded, nor does a later exam exist with pathology results. All Exams are Signed exams.

Parameters:

- Practice, Site(s), Date Range and Signing Radiologists

### MQSA List of Positive Exams with Outstanding Biopsies


Practice: Treesdale Radiology

Site(s): Cedarville, Mapleton, Pine Grove, Outside Read1 (outside read)

Date Range: 01-01-2015 - 11-08-2016

Signing Radiologist: All radiologists for selected practice

Description: This report displays the list of Positive Mammography Exams ( for selected time period, location, radiologist) that do not have pathology result recorded, nor does a later exam exist with pathology results. All Exams are Signed exams.

BI-RADS®	Accession	Patient Name	Procedure	Signed Date	Signing Rad	Biopsy Performed
<input type="checkbox"/> BI-RAD® 0	3					
	1018485	James Sally,	MA Digital Mammo Screen Bilat [G0202] - Breast - Bilateral	12-28-2015	Saltmarsh, Hilary	N
	1025281	Olivens Jenny,	MA Digital Mammo Screen Bilat [G0202] - Breast - Bilateral	12-16-2015	Saltmarsh, Hilary	N
	1024217	Test Portal, middletest	MA Digital Mammo Screen Bilat [G0202] - Breast - Bilateral	11-26-2015	Noye, Darcy Thane	N
<input type="checkbox"/> BI-RAD® 4	7					
	1023211	Zulu Average,	MA Digital Mammo Diag Bilat [G0204] - Breast - Bilateral	11-25-2015	Saltmarsh, Hilary	N
	1027770	Adams Average, Demo	MA Digital Mammo Diag Bilat [G0204] - Breast - Bilateral	11-25-2015	Saltmarsh, Hilary	N
	1027780	O'Brian Hillary, Demo	MA Digital Mammo Diag Bilat [G0204] - Breast - Bilateral	11-25-2015	Saltmarsh, Hilary	N
	1027790	O'Conner Average, Demo	MA Digital Mammo Diag Bilat [G0204] - Breast - Bilateral	11-25-2015	Saltmarsh, Hilary	N
	1027800	O'Dell Jayne, Demo	MA Digital Mammo Diag Bilat [G0204] - Breast - Bilateral	11-25-2015	Saltmarsh, Hilary	N
	1027810	O'Neil Average, Demo	MA Digital Mammo Diag Bilat [G0204] - Breast - Bilateral	11-25-2015	Saltmarsh, Hilary	N
	1027820	ד'אור דוד, Demo	MA Digital Mammo Diag Bilat [G0204] - Breast - Bilateral	11-25-2015	Saltmarsh, Hilary	N
<b>Report Total:</b>						<b>10</b>



### MQSA Location Summary Group By BI-RADS

Description: This report displays the Mammography Site Reporting Summary for selected sites/practice. All Exams are Signed exams.

**Sensitivity Definition:** Site/Practice TP Total for selected period / (System Wide TP + System Wide FN) within the same time period.

**Specificity Definition:** Site/Practice TN Total for selected period / (System Wide TN + System Wide FP) within the same time period.

**PPV3 Definition:** Total (Screening or Diagnostic) TP exams for selected period / (System Wide Total of Biopsies) within the same time period.

**Recall Rate:** Total (Screening or Diagnostic) Positive exams for selected period / (System Wide Total Mammo Exams) within the same time period.

Parameters:

- Practice, Site(s), Date Range, Procedure Codes and Signing Radiologists

### MQSA Location Summary Group By BI-RADS

Practice: Treesdale

Site(s): Cedarville

Date Range: 11-08-2016 - 11-10-2016

Signing Radiologist: All radiologists for selected practice

Procedure Codes: MA999, G0204, G0204I, G0204

Description: This report displays the Mammography Site Reporting Summary for selected sites/practice. All Exams are Signed exams.

**Sensitivity Definition:** Site/Practice TP Total for selected period / ( System Wide TP + System Wide FN) w ithin the same time period.

**Specificity Definition:** Site/Practice TN Total for selected period / ( System Wide TN + System Wide FP) w ithin the same time period.

**PPV3 Definition:** Total (Screening or Diagnostic) TP exams for selected period / ( System Wide Total of Biopsies) w ithin the same time period.

**Recall Rate:** Total (Screening or Diagnostic) Positive exams for selected period / ( System Wide Total Mammo Exams) w ithin the same time period.

Site / Location:	Total	Positive	TP	TN	FP1	FN	Sensitivity	Specificity	PPV3	Recall Rate
CEDARVILLE	146									
☑ DIAG	18	3	0	14	0	0	0.000	1.000	0.00	
☑ SCREEN	122	10	0	109	0	0	0.000	1.000	0.00	8.20%
☐ US_DIAG	6		0	6	0	0	0.000	1.000	0.00	
☑ BI-RADS 1	3			3						
☐ BI-RADS 2	2			2						
US Breast Unilateral Or Bilateral	1			1						
US Breast Bilateral	1			1						
☐ BI-RADS 3	1			1						
US Breast Unilateral	1			1						
<b>Report Total:</b>	<b>146</b>	<b>13</b>		<b>129</b>			<b>0.000</b>	<b>1.000</b>	<b>0.00</b>	

### MQSA Location Summary Group By Procedure

Description: This report displays the Mammography Site Reporting Summary for selected sites/practice. All Exams are Signed exams.

**Sensitivity Definition:** Site/Practice TP Total for selected period / (System Wide TP + System Wide FN) within the same time period.

**Specificity Definition:** Site/Practice TN Total for selected period / (System Wide TN + System Wide FP) within the same time period.

**PPV3 Definition:** Total (Screening or Diagnostic) TP exams for selected period / (System Wide Total of Biopsies) within the same time period.

**Recall Rate:** Total (Screening or Diagnostic) Positive exams for selected period / (System Wide Total Mammo Exams) within the same time period.

Parameters:

- Practice, Site(s), Date Range, Procedure Codes and Signing Radiologists

**Signing Radiologist:** Saltmarsh, Hilary (saltmarsh)

**Procedure Codes:** 76641, 76641B, 76642, 76642B, G0202, G0202IMP, G0202UIMP, G0202UNI, G0202WC, G0204, G0204/US, G0204CB, G0204CB/US, G0204IMP, G0206, G0206/US, G0206CB, G0206CB/US, G0206IMP, MATO1, MATO2

**Description:** This report displays the Mammography Site Reporting Summary for selected sites/practice. All Exams are Signed exams.

**Sensitivity Definition:** Site/Practice TP Total for selected period / (System Wide TP + System Wide FN) within the same time period.

**Specificity Definition:** Site/Practice TN Total for selected period / (System Wide TN + System Wide FP) within the same time period.

**PPV3 Definition:** Total (Screening or Diagnostic) TP exams for selected period / (System Wide Total of Biopsies) within the same time period.

**Recall Rate:** Total (Screening or Diagnostic) Positive exams for selected period / (System Wide Total Mammo Exams) within the same time period.

Site / Location:	Accession	Positive	TP	TN	FP1	FN	Sensitivity	Specificity	PPV3	Recall Rate	Biopsy Count
Cedarville		12									
<input type="checkbox"/> DIAG		10	8	2	1		0	0.667	0.00		
<input type="checkbox"/> MA Digital Mammo Diag Bilat		9	8	1	1						
<input type="checkbox"/> BI-RADS®1		1		1							
	1025483	N		1							
<input type="checkbox"/> BI-RADS®4		8	8			1					
	1023211	Y									
	1027770	Y									
	1027780	Y									
	1027790	Y									
	1027800	Y									
	1027810	Y									
	1027820	Y									
	1025446	Y				1					
<input type="checkbox"/> MA Digital Mammo Diag Unilat		1		1							
<input type="checkbox"/> BI-RADS®2		1		1							
	1025595	N		1							
<input type="checkbox"/> SCREEN		2	1				0		0.00	0.50	
Mapleton		1									
<input type="checkbox"/> SCREEN		1	1				0		0.00	1.00	
Pine Grove		2									
<input type="checkbox"/> DIAG		1		1			0	1.000	0.00		
<input type="checkbox"/> SCREEN		1					0		0.00		
<b>Report Total:</b>		<b>15</b>	<b>10</b>	<b>0</b>	<b>3</b>	<b>1</b>	<b>0</b>	<b>0.750</b>	<b>0.00</b>	<b>0.67</b>	

### MQSA Outcome Report

Description: This report displays outcome results for BI-RADS® 4 and BI-RADS® 5 exams. All Exams are Signed exams.

Parameters:

- Practice, Site(s), Date Range, Signing Radiologists and Procedure Codes

### MQSA Outcome Report

Practice: Treesdale Radiology

Site(s): Cedarville, Mapleton, Pine Grove, Outside Read1 (outside read)

Date Range: 01-01-2015 - 11-08-2016

Signing Radiologist: Saltmarsh, Hilary (saltmarsh)

Procedure Codes: 76641, 76641B, 76642, 76642B, G0202, G0202IMP, G0202UIIMP, G0202UNI, G0202WC, G0204, G0204/US, G0204CB, G0204CB/US, G0204IMP, G0206, G0206/US, G0206CB, G0206CB/US, G0206IMP, MATO1, MATO2

Description: This report displays outcome results for BI-RADS® 4 and BI-RADS® 5 exams. All Exams are Signed exams.

Signing Radiologist	# Studies	# Benign ( - )	# Malignant ( + )	No Biopsy on file	% ( - )	% ( + )	% NA
<input type="checkbox"/> Saltmarsh, Hilary	8	1	0	7	12.50 %	0.00 %	87.50 %
	MRN: MASTER_CLONE_1	Zulu Average,	ACC:1023211	NA			
	MRN:z10000001	Adams Average, Demo	ACC:1027770	NA			
	MRN:z20000002	O'Brian Hilary, Demo	ACC:1027780	NA			
	MRN:z30000003	O'Conner Average, Demo	ACC:1027790	NA			
	MRN:z40000004	O'Dell Jayne, Demo	ACC:1027800	NA			
	MRN:z50000005	O'Neil Average, Demo	ACC:1027810	NA			
	MRN: 1501633RSNA	Murcia Hilary,	ACC:1025446	Benign			
<b>Total:</b>	<b>8</b>	<b>1</b>	<b>0</b>	<b>7</b>	<b>12.50 %</b>	<b>0.00 %</b>	<b>87.50 %</b>

### MQSA Radiologist Reporting Summary

Description: This report displays the Radiologist Reporting Summary for selected sites/practice. All Exams are Signed exams.

Parameters:

- Practice, Site(s), Date Range, Signing Radiologists and Procedure Codes

### MQSA Radiologist Reporting Summary

**Practice:** Treesdale Radiology

**Site(s):** Cedarville, Mapleton, Pine Grove, Outside Read1 (outside read)

**Date Range:** 01-01-2015 - 11-08-2016

**Signing Radiologist:** Noye, Darcy (darcyn), Saltmarsh, Hilary (saltmarsh)

**Procedure Codes:** 76641, 76641B, 76642, 76642B, G0202, G0202IMP, G0202UIIMP, G0202UNI, G0202WC, G0204, G0204/US, G0204CB, G0204CB/US, G0204IMP, G0206, G0206/US, G0206CB, G0206CB/US, G0206IMP, MATO1, MATO2

**Description:** This report displays the Radiologist Reporting Summary for selected sites/practice. All Exams are Signed exams.

Site	Signing Radiologist	Accession	Patient Name	TP	TN	FP1	FN
Cedarville	<input type="checkbox"/> Saltmarsh, Hilary	12					
	<input type="checkbox"/> DIAG	10		0	2	1	0
	<input type="checkbox"/> Procedure: MA Digital Mammo Diag Bilat [G0204] - Breast - Bilateral	9					
	<input type="checkbox"/> BI-RADS®1	1					
		1025483	Test Hilary,	N	Y	N	N
	<input type="checkbox"/> BI-RADS®4	8					
		1023211	Zulu Average,	N	N	N	N
		1027770	Adams Average, Demo	N	N	N	N
		1027780	O'Brian Hilary, Demo	N	N	N	N
		1027790	O'Conner Average, Demo	N	N	N	N
		1027800	O'Dell Jayne, Demo	N	N	N	N
		1027810	O'Neil Average, Demo	N	N	N	N
		1027820	ד'אור דוד, Demo	N	N	N	N
		1025446	Murcia Hilary,	N	N	Y	N
	<input type="checkbox"/> Procedure: MA Digital Mammo Diag Unilat [G0206] - Breast - Left	1					
<input type="checkbox"/> SCREEN	2			0	0	0	0
Mapleton	<input type="checkbox"/> Saltmarsh, Hilary	1					
Pine Grove	<input type="checkbox"/> Noye, Darcy	1					
	<input type="checkbox"/> Saltmarsh, Hilary	2					
<b>Report Total:</b>		<b>16</b>					

## Management Reports \ Meaningful Use

### Automated Measures by Radiologist

Description: Identical to the report below called Automated Measures, but returns a page for each radiologist associated with the practices selected. Each measure contains a hyperlink to a sub-report that will provide the underlying patient data.

Parameters:


- Date Range, MU Stage (1, 2 or 3) and Practice

#### Automated Measures by Radiologist Stage 2

**Radiologist:** Eckstein, Rhonda

**Practice(s):** Rainbow Radiology, Spectrum Radiology, Treesdale Radiology

**Reporting Period:** 01-01-2016 - 12-31-2016



Meaningful Use Measure	Numerator	Denominator	%	Goal	Pass/Fail
Demographics	0	1	0.00%	80%	Fail
Patient reminders	0	0	0.00%	10%	Fail
Medication reconciliation	0	0	0.00%	50%	Fail
Vitals and BMI	0	1	0.00%	80%	Fail
Vitals - BP Out of Scope	1	1	100.00%	80%	Pass
Vitals - Ht/Wt Out of Scope	0	1	0.00%	80%	Fail
Smoking status	0	1	0.00%	80%	Fail
Family health history	0	1	0.00%	20%	Fail
Electronic notes	0	1	0.00%	30%	Fail
CPOE med	0	0	0.00%	60%	Fail
CPOE lab	0	0	0.00%	30%	Fail
CPOE rad	0	0	0.00%	30%	Fail
eRx	0	0	0.00%	50%	Fail
Lab results	0	0	0.00%	55%	Fail
Imaging	0	0	0.00%	10%	Fail
Patient education	0	0	0.00%	10%	Fail
Clinical summary	0	0	0.00%	50%	Fail

### Automated Measures Raw Data

Description: This report displays the percentage of radiologists where the mu measures were recorded.

Parameters:

- Date Range, MU Stage (1 or 2) and Practice(s)

### Automated Measures Raw Data - Stage 2

**Date Range:** 01-01-2016 - 12-31-2016

**Practice(s):** Treesdale Radiology

**Description:** This report displays the percentage of radiologists where the mu measures were recorded.

Description	Numerator	Denominator	%	MU rule	Pass	Signed by user id	User Name
Demographics	0	1	0.00%	80%	N	eckstein	Eckstein, Rhonda
Patient reminders	0	0	0.00%	10%	N	eckstein	Eckstein, Rhonda
Medication reconciliation	0	0	0.00%	50%	N	eckstein	Eckstein, Rhonda
Vitals and BMI	0	1	0.00%	80%	N	eckstein	Eckstein, Rhonda
Vitals - BP Out of Scope	1	1	100.00%	80%	Y	eckstein	Eckstein, Rhonda
Vitals - Ht/Wt Out of Scope	0	1	0.00%	80%	N	eckstein	Eckstein, Rhonda
Smoking status	0	1	0.00%	80%	N	eckstein	Eckstein, Rhonda
Family health history	0	1	0.00%	20%	N	eckstein	Eckstein, Rhonda
Electronic notes	0	1	0.00%	30%	N	eckstein	Eckstein, Rhonda

### Automated Measures

Description: This report displays the percentage of patients where the mu requirement was recorded. Each measure contains a hyperlink to a sub-report that will provide the underlying patient data.

Parameters:

- Date Range, Practice(s), Radiologist and MU Stage (1, 2 or 3)

#### Automated Measures Stage 2


**Radiologist:** Eckstein, Rhonda

**Practice(s):** Rainbow Radiology, Spectrum Radiology, Treesdale Radiology

**Reporting Period:** 01-01-2016 - 12-31-2016

**Description:** This report displays the percentage of patients where the mu requirement was recorded.

\*Click [here](#) to get the details.



Meaningful Use Measure	Numerator	Denominator	%	Goal	Pass/Fail
Demographics	0	1	0.00%	80%	Fail
Patient reminders	0	0	0.00%	10%	Fail
Medication reconciliation	0	0	0.00%	50%	Fail
Vitals and BMI	0	1	0.00%	80%	Fail
Vitals - BP Out of Scope	1	1	100.00%	80%	Pass
Vitals - Ht/Wt Out of Scope	0	1	0.00%	80%	Fail
Smoking status	0	1	0.00%	80%	Fail
Family health history	0	1	0.00%	20%	Fail
Electronic notes	0	1	0.00%	30%	Fail
CPOE med	0	0	0.00%	60%	Fail
CPOE lab	0	0	0.00%	30%	Fail
CPOE rad	0	0	0.00%	30%	Fail
eRx	0	0	0.00%	50%	Fail
Lab results	0	0	0.00%	55%	Fail
Imaging	0	0	0.00%	10%	Fail
Patient education	0	0	0.00%	10%	Fail
Clinical summary	0	0	0.00%	50%	Fail

### Automated Measures Details

Description: This report displays the underlying patient data for each of the measures on the “Automated Measures” and “Automated Measures by Radiologist” reports. It can be accessed by clicking the blue hyperlinks in these reports.

Patient Name	Patient ID	Accession #	Practice Code	Scheduled Date	Signed Date	Office Visit	Status Code	Patient Education
Carter, Jeff William	4	1	ADV	12-04-2014	12-05-2014	N	Signed1	N
		3398	ADV	07-20-2015	07-20-2015	N	Signed1	N
		4878	ADV	10-06-2015	11-02-2015	N	Signed1	N
		5105	ADV	10-22-2015	11-02-2015	N	Signed1	N
		5571	ADV	11-02-2015	11-02-2015	N	Signed1	N
Hextall, Ron	12	107	ADV	03-27-2015	06-29-2015	N	Signed1	N
		120	ADV	04-17-2015	06-29-2015	N	Signed1	N
Judson, Jenny	3392	1896	ADV	05-14-2015	10-28-2015	N	Signed1	Y
Obermann, Mary	654	4413	ADV	08-13-2015	09-02-2015	N	Signed1	N
		4510	ADV	08-27-2015	08-27-2015	N	Signed1	N
test, a12501 Ben	893247	5816	ADV	12-04-2015	12-04-2015	N	Signed1	N
Test, IVT	50	5108	ADV	10-22-2015	10-22-2015	N	Signed1	N



**CQM Measures**

Description: This report displays the patient specific details used to generate the CMS measures report.

Note: If the measure is not displayed in the list there are no patients for the CMS measure.

Parameters:

- Date Range, Practice(s), Radiologist and MU Stage

**CQM Measures Stage 2**

**Radiologist:** Eckstein, Rhonda (eckstein)

**Date Range:** 01-01-2016 - 12-31-2016

**Practice(s):** Rainbow Radiology, Spectrum Radiology, Treesdale Radiology

**Description:** This report displays the patient specific details used to generate the CMS measures report. Note: If the measure is not displayed in the list there are no patients for the CMS measure.

CQM	Code	Description	Patient Population	Denominator	Denominator Exclusions	Numerator	Denominator Exceptions
CMS50		Closing the referral loop	0	0	0	0	0
CMS65		Hypertension: Improvement in blood pressure	0	0	0	0	0
CMS69	0421a	Preventive Care and Screening: Body Mass Index (BMI) Screening and Follow-Up	0	0	0	0	0
CMS69	0421b	Preventive Care and Screening: Body Mass Index (BMI) Screening and Follow-Up	0	0	0	0	0
CMS117	0038	Childhood Immunization Status	0	0	0	0	0
CMS125	0031	Breast Cancer Screening	0	0	0	0	0
CMS127	0043	Pneumonia Vaccination Status for Older Adults	0	0	0	0	0
CMS138	0028	Preventive Care and Screening: Tobacco Use: Screening and Cessation Intervention	0	0	0	0	0
CMS155	0024a	Weight Assessment and Counseling for Nutrition and Physical Activity for Children and Adolescents	0	0	0	0	0
CMS155	0024b	Weight Assessment and Counseling for Nutrition and Physical Activity for Children and Adolescents	0	0	0	0	0
CMS155	0024c	Weight Assessment and Counseling for Nutrition and Physical Activity for Children and Adolescents	0	0	0	0	0
CMS155	0024d	Weight Assessment and Counseling for Nutrition and Physical Activity for Children and Adolescents	0	0	0	0	0
CMS155	0024e	Weight Assessment and Counseling for Nutrition and Physical Activity for Children and Adolescents	0	0	0	0	0
CMS155	0024f	Weight Assessment and Counseling for Nutrition and Physical Activity for Children and Adolescents	0	0	0	0	0
CMS155	0024g	Weight Assessment and Counseling for Nutrition and Physical Activity for Children and Adolescents	0	0	0	0	0
CMS155	0024h	Weight Assessment and Counseling for Nutrition and Physical Activity for Children and Adolescents	0	0	0	0	0
CMS155	0024i	Weight Assessment and Counseling for Nutrition and Physical Activity for Children and Adolescents	0	0	0	0	0
CMS165	0018	Controlling High Blood Pressure	0	0	0	0	0



### Patient Lists

Description: This report displays the patients who meet the entered parameters. Note: The date(s) fields indicate the starting date the problem, medication, allergy or lab occurred.

Parameters:

- From and To Age, From and To Date (DOS), Primary Language, Sex, Race, Ethnicity, Problem, Medication, Medication Allergy, Exclude Deceased, Lab Test, Lab results, Preferred Communication and Max Results.

### Patient Lists

Description: This report displays the patients who meet the entered parameters. Note: The date(s) fields indicate the starting date the problem, medication, allergy or lab occurred.

**Parameters**

MRN	Last Name	First Name	Age	Sex	Primary Language	Race	Ethnicity	Preferred Communication
101003	Dean	Theodore		63 M	English	Declined to Specify	Declined to Specify	Not assigned
104782	Tester	Test		6 M	English	White	Hispanic	Not assigned
10531	B Goode	Johnny		32 M	English	No race assigned		Not assigned
10542	Adams	Average		22 M	English	No race assigned		Phone
10556	Zulu	Average		27 M	English	No race assigned		Not assigned
14194	Moran	Keisha		24 F	English	White	Not Hispanic	Not assigned
1500007	Luoto	Demetra		50 M	English	American Indian or Alaska Native	Hispanic	Not assigned
1500015	Janicki	Jerrie		26 M	English	American Indian or Alaska Native	Hispanic	Not assigned
1500045	Henrikson	Francesca		9 M	English	American Indian or Alaska Native	Hispanic	Not assigned
1500463	Hapke	Towanda		17 M	English	White	Not Hispanic	Not assigned
1500638	Gremo	Audrea		32 M	English	American Indian or Alaska Native	Hispanic	Not assigned

## Management Reports \ Operations

### Aborted Exams

Description: This report shows the counts for all aborted exams for a selected practice, grouped by site, and reason with a drill down to the accession

Parameters:

- Practice, Site(s), Date Range, and reasons

### Aborted Exams

Practice: Treesdale

Site(s): CEDARVILLE

Date Range: 01-01-2016 - 11-08-2016

Description: This report shows the counts for all aborted exams for a selected practice, grouped by site, and reason with a drill down to the accession

Accession #	BRE	CT	DIA	MR	US	Total
CEDARVILLE	2	1	2	1	4	10
☐ Other	2	1	1	1	2	7
☐ Radiologist Changed/Cancelled	0	0	0	0	1	1
☐ Scheduled Incorrectly	0	0	1	0	1	2
<a href="#">123456</a>						
<a href="#">213123</a>						
<b>Total</b>	<b>2</b>	<b>1</b>	<b>2</b>	<b>1</b>	<b>4</b>	<b>10</b>

Aborted Exams

Confidential and Proprietary  
Page 1 of 4

Date Run: 11-08-2016 10:48 PM  
RADNET\E:ksteinR

### Addended Mammography Exams

Description: This report displays a list of mammography exams that have been addended and dictated within the selected Date Range. Identification of a mammography exam is based on that procedures' mammography type value being configured with RIS (Diagnostic or Screening).

Parameters:

- Practice(s), Site(s), Date Range

### Addended Mammography Exams

**Practice(s):** Treesdale Radiology

**Site(s):** All sites for the selected practice(s).

**Date Range:** 01-01-2014 - 11-08-2016

**Description:** This report displays a list of mammography exams that have been addended and dictated within the selected date range. Identification of a mammography exam is based on the procedures' mammo type value being configured.

Acc #	MRN	Patient Name	Status	Dictated by	Signed by
1024376	1501585RSN A	Selma Warren	Signed1	saltmarsh	saltmarsh
	<b>Procedure(s):</b> MA Digital Mammo Screen Unilat [G0202UNI] - Breast - Left			11-27-2015 8:47 AM	11-27-2015 8:48 AM

Addended Mammography Exams

Confidential and Proprietary  
Page 1 of 1

Date Run: 11-08-2016 7:09 PM  
RRIS-RSNA\FORMS\GuestUser

## Appointment Schedule

Description: This report displays the details for exams scheduled on the specified date for the practice, site(s) and scanners selected.

Parameters:

- Practice, Site(s), Modality Type(s), Scanner(s), Schedule Date, Display Scheduler Notes, Display Reason for Exam

## Appointment Schedule

**Practice:** Treesdale Radiology

**Site(s):** Cedarville, Mapleton, Pine Grove

**Date:** 11-07-2016

**Modality Type(s):** AN, BD, CT, EN, FL, MA, MR, NM, PT, RT, SC, CH, SP, US, XR

**Scanner(s):** CT1CE, CT1MP, CT1PG, DE1MP, DE1PG, MA1CE, MA1MP, MA1PG, MA2CE, MA2PG, MA3PG, MR1CE, MR1MP, MR1PG, MR2CE, MR2PG, MR3CE, NM1PG, NM2PG, PT1MP, PT1PG, US1CE, US1MP, US1PG, US2CE, US2MP, US2PG, US3PG, XR1CE, XR1MP, XR1PG, XR2PG

**\*Scheduler notes suppressed**

**Description:** This report displays the details for exams scheduled on the specified date for the practice, site(s) and scanners selected.

Scanner	Time Sched	Patient Name	MRN#	DOB	Exam(s)
CT1CE	1:05 PM	Strickland, Rhonda	9695	07-17-1983	CT Abdomen W [74160] - Abdomen <b>Body Part\Lat:</b> Abdomen/
		<b>Home Phone:</b> (256) 224-5160	<b>Age:</b> 33 y	<b>Referring:</b> Arnold, Bob	<b>Copay:</b>
		<b>Cell Phone:</b> (256) 224-5160			
		Current Exam Status: <b>ExamDone</b>		<b>Enc #:</b> 48469	<b>Acc #:</b> 1030889
		Insurance: NATIONWIDE INSURANCE		Insurance:	
CT1MP	10:50 AM	Freddieson, Fred	y6	05-09-1973	CT Head Wo [70450] - Head <b>Body Part\Lat:</b> Head/
		<b>Home Phone:</b> (902) 374-5967	<b>Age:</b> 43 y	<b>Referring:</b> Flowers, Tim	<b>Copay:</b>
		<b>Cell Phone:</b>			
		Current Exam Status: <b>Scheduled</b>		<b>Enc #:</b> 46452	<b>Acc #:</b> 1027871
		Insurance:		Insurance:	
CT1PG	9:00 AM	Hartman, Loren	1378	08-02-1995	CT Abdomen W [74160] - Abdomen <b>Body Part\Lat:</b> Abdomen/
		<b>Home Phone:</b>	<b>Age:</b> 21 y	<b>Referring:</b> Arnold, Bob	<b>Copay:</b>
		<b>Cell Phone:</b>			
		Current Exam Status: <b>Scheduled</b>		<b>Enc #:</b> 47469	<b>Acc #:</b> 1029889
		Insurance: PARAMOUNT INSURANCE		Insurance:	
MR1PG	11:00 AM	Harrison, Harry	y8	03-07-1982	MR Cervical Spine Wo [72141] - Cervical spine <b>Body Part\Lat:</b> Cervical spine/
		<b>Home Phone:</b> (902) 675-6867	<b>Age:</b> 34 y	<b>Referring:</b> Howard, Tom	<b>Copay:</b>
		<b>Cell Phone:</b>			
		Current Exam Status: <b>Scheduled</b>		<b>Enc #:</b> 46454	<b>Acc #:</b> 1027874
		Insurance:		Insurance:	

## Appointment Schedule by Procedure Group

Description: This report displays the details for exams scheduled on the specified date for the practice, site(s) and procedure groups selected.

Parameters:

- Practice, Site(s), Schedule Date, Display Scheduler Notes, Display Reason for Exam and Procedure Group(s).

### Appointment Schedule by Procedure Group

**Practice:** Treesdale Radiology

**Site(s):** Cedarville, Mapleton, Pine Grove

**Date:** 11-07-2016

**Procedure Group(s):** DEXA, Diagnostic CT, Diagnostic Imagt, EPET Procedures, Evaluation Mgmt, Fluoroscopy, Intervntnl Rad, Intervntnl Surg, Mammography, MRI and MRA, Nuclear Med, OTHER, PET, Plain Film/Xray/CR, Radiation Oncol, Ultrasound Gen, Ultrasound Vasc

**Description:** This report displays the details for exams scheduled on the specified date for the practice, site(s) and procedure groups selected.

Scanner	Time Sched	Patient Name	MRN#	DOB	Exam(s)
CT1CE	1:05 PM	Strickland, Rhonda	9695	07-17-1983	CT Abdomen W [74160] - Abdomen <b>Body Part/Lat:</b> Abdomen/
		<b>Home Phone:</b> (256) 224-5160	<b>Age:</b> 33 y	<b>Referring:</b> Arnold, Bob	<b>Copay:</b>
		<b>Cell Phone:</b> (256) 224-5160			
		<b>Current Exam Status:</b> ExamDone		<b>Enc #:</b> 48469	<b>Acc #:</b> 1030889
		Insurance: NATIONWIDE INSURANCE		Insurance:	
		Scheduler Notes: 11-07-16 eckstein - Informed patient about estimated financial responsibility. 11-07-16 eckstein - send patient back to office			
		Reason for Exam: Pain			
CT1MP	10:50 AM	Freddieson, Fred	y6	05-09-1973	CT Head Wo [70450] - Head <b>Body Part/Lat:</b> Head/
		<b>Home Phone:</b> (902) 374-5967	<b>Age:</b> 43 y	<b>Referring:</b> Flowers, Tim	<b>Copay:</b>
		<b>Cell Phone:</b>			
		<b>Current Exam Status:</b> Scheduled		<b>Enc #:</b> 46452	<b>Acc #:</b> 1027871
		Insurance:		Insurance:	
		Scheduler Notes:			
		Reason for Exam: chest pain			
CT1PG	9:00 AM	Hartman, Loren	1378	08-02-1995	CT Abdomen W [74160] - Abdomen <b>Body Part/Lat:</b> Abdomen/
		<b>Home Phone:</b>	<b>Age:</b> 21 y	<b>Referring:</b> Arnold, Bob	<b>Copay:</b>
		<b>Cell Phone:</b>			
		<b>Current Exam Status:</b> Scheduled		<b>Enc #:</b> 47469	<b>Acc #:</b> 1029889
		Insurance: PARAMOUNT INSURANCE		Insurance:	
		Scheduler Notes:			
		Reason for Exam: Pain			
MR1PG	11:00 AM	Harrison, Harry	y8	03-07-1982	MR Cervical Spine Wo [72141] - Cervical spine <b>Body Part/Lat:</b> Cervical spine/

### Blood Work Appointment Schedule

Description: This report displays the details for exams scheduled on the specified date for the practice, site(s) and scanners selected. \*To run as designed this report is looking for a sedation type with a code of BLDWK.

Parameters:

- Practice, Site(s), Modality Type(s), Scanner(s), Schedule Date, Display Scheduler Notes, Display Reason for Exam, Sedation Type Code

### Blood Work Appointment Schedule

**Practice:** Treesdale Radiology

**Site(s):** Cedarville, Mapleton, Pine Grove

**Date:** 11-10-2016

**Modality Type(s):** AN, BD, CT, EN, FL, MA, MR, NM, PT, RT, SC, CH, SP, US, XR

**Scanner(s):** CT1CE, CT1MP, CT1PG, DE1MP, DE1PG, MA1CE, MA1MP, MA1PG, MA2CE, MA2PG, MA3PG, MR1CE, MR1MP, MR1PG, MR2CE, MR2PG, MR3CE, NM1PG, NM2PG, PT1MP, PT1PG, US1CE, US1MP, US1PG, US2CE, US2MP, US2PG, US3PG, XR1CE, XR1MP, XR1PG, XR2PG

**\*Scheduler notes suppressed**

**Description:** This report displays the details for exams scheduled on the specified date for the practice, site(s) and scanners selected. \*To run as designed this report is looking for a sedation type with a code of BLDWK.

Scanner	Time Sched	Patient Name	MRN#	DOB	Exam(s)
CT1CE	10:30 AM	Hall, Josh	3233	12-11-1971	CT Abdomen W [74160] - Abdomen <b>Body Part/Lat:</b> Abdomen/
	<b>Home Phone:</b>		<b>Age:</b> 44 y	<b>Referring:</b> Arnold, Bob	
	<b>Cell Phone:</b>			<b>Copay:</b>	
	Current Exam Status: <b>Scheduled</b>			<b>Enc #:</b> 48483	<b>Acc #:</b> 1030910
	Insurance: Private Pay			Insurance:	
	Sedation Type: Bloodwork				
					<b>Total Appointments: 1</b>

### Cancelations by Referring

Description: This report displays the cancelations with details made by referring physicians selected by practice, site(s), date range, physician(s) and modality.

Parameters:

- Practice, Site(s), Date Range, Referring Physician(s) and Modality(s)

#### Cancelations by Referring

Practice: Treeddale Radiology

Site(s): Cedarville, Mapleton, Pine Grove

Date Range: 09-25-2015 - 11-08-2016

Description: This report displays the cancelations with details made by referring physicians selected by practice, site(s), date range, physician(s) and modality.

Referring Physician	Site	Accession #	Patient	Procedure & Cancellation Reason	
Edwards, John				Total for Edwards, John:	2
	<input type="checkbox"/> Mapleton			Total Canceled at Mapleton:	2
		1027777	O'Brian, Hillary	PT/CT Sodium Fluoride PET Bone Scan [PT4] Patient not available	PT
		1027797	O'Dell, Jayne	PT/CT Sodium Fluoride PET Bone Scan [PT4] No Show	PT
<b>Total Cancelations:</b>					<b>2</b>

Cancelations by Referring

Confidential and Proprietary  
Page 1 of 1

Date Run: 11-08-2016 12:43 PM  
RRIS-RSNA/FormsGuestUser



### Canceled Exams

Description: This report shows the cancelation counts by reason and procedure group for the selected Practice, Site(s), and Date Range. Ability to link to patient details via blue hyperlinks.

Parameters:

- Practice, Site(s), and Date Range

#### Canceled Exams

Practice: Treedale Radiology

Site(s): All sites for the selected practice(s).

Date Range: 07-01-2016 - 11-08-2016

Description: This report shows the cancelation counts by reason and procedure group for the selected Practice, Site(s), and Date Range (based on DOS).

Site	Cancel Date	CT	MR	PT	UG	XR	Total
<b>Cedarville</b>		2	0	0	0	1	3
<input type="checkbox"/> No Show						1	1
<a href="#">1025605</a>	09-21-2016					1	1
<input type="checkbox"/> Patient not available		1					1
<input type="checkbox"/> Patient Sick		1					1
<b>Mapleton</b>		1	0	1	0	0	2
<input type="checkbox"/> Patient not available				1			1
<input type="checkbox"/> Room Issue		1					1
<b>Pine Grove</b>		0	1	0	0	1	2
<input type="checkbox"/> Unable to finance			1			1	2
<a href="#">1025563</a>	06-30-2016					1	1
<a href="#">1027873</a>	11-07-2012		1				1
		3	1	1	1	1	7

Canceled Exams

Confidential and Proprietary  
Page 1 of 1

Date Run: 11-08-2016 12:51 PM  
RRIS-RSNA\FORMS\GuestUser

### Accession Milestone Times

Patient: 1501558RSNA - Ziegler, Randy

Current Exam Status: **Cancelled by patient\_key:170783 6/30/2016 12:34 PM**

Site: Pine Grove

Accession #	Procedure Category	Procedure Description
1025563	XR	XR Orbits for Foreign Body [70030] - Head

Exam Milestone Times					
Order Date	Scheduled Date	Checking In Date	Arrival Date	Started Date	Performed End Date
06-30-2016 10:28 AM	07-01-2016 9:00 AM				
patient_key:170783	patient_key:170783				

Report Milestone Times				
Version	Dictated date	Drafted Date	Edited Date	Signed Date

Turnaround Times (in Hrs)				
<i>**calculated at 8:00-6:00 Mon-Fri</i>				
Done To Dictated	Dictated To Drafted	Drafted To Edited	Edited To Signed	Done To Signed
0.00	0.00	0.00	0.00	0.00

Turnaround Times (in Hrs)				
<i>**calculated at 24 / 7</i>				
Done To Dictated	Dictated To Drafted	Drafted To Edited	Edited To Signed	Done To Signed

### Canceled Reasons

Description: This report shows the cancellation counts by reason and procedure group for the selected Practice, Site(s), and Date Range (based on DOS).

Parameters:

- Practice, Site(s), Date Range based on date of service

#### Canceled Reasons

**Practice(s):** Treesdale Radiology

**Site(s):** Cedarville, Mapleton, Pine Grove

**Date Range:** 01-01-2016 - 11-08-2016

**Reason(s):** Lack of Transport,Late,No Show,Patient not available,Patient Sick,Rescheduled,Room Issue,Unable to finance,Unknown

**Room(s):** All rooms for selected sites.

**Rescheduled:** Yes,No

**Description:** Based on scheduled date of service, this report shows the exams that were canceled for the reason(s) selected, and displays which of those exams have been rescheduled.

Scheduled date	Procedure(s)	Room	First Name	Last Name	Acc. #	Rescheduled?
Reason: No Show						1/2
05-13-2016 8:09 PM	CT Abdomen W & Pelvis W & CTA Chest W/Wo [CT54] - Abdomen	CT 1 at Cedar	Hilary	Test	<a href="#">1025550</a>	N
09-29-2016 11:01 AM	US Abdominal Complete [76700] - Abdomen	US 1 at Cedar	Hilary	Test	<a href="#">1025605</a>	Y
Reason: Patient not available						1/7
08-25-2016 11:20 AM	CT Abdomen W & Pelvis W [CT9] - Abdomen	CT 1 at Cedar	AAAAlpha	BBBBeta	<a href="#">1025598</a>	N
11-01-2016 11:53 AM	PT/CT Sodium Fluoride PET Bone Scan [PT4]	PET 1 at Maple	Hillary	O'Brian	<a href="#">1027777</a>	N
03-11-2016 3:00 PM	MR Cervical Spine W [72142] - Cervical spine	MR1 at Maple	Jane	O'Dell	<a href="#">1025448</a>	N
03-11-2016 3:30 PM	MA Digital Mammo Screen Bilat [G0202] - Breast - Bilateral	Mammo 1 at Maple	Jane	O'Dell	<a href="#">1025447</a>	N
01-12-2016 8:00 AM	MR Chest W [71551] - Chest	MR 1 at Cedar	Hilary	Saltmarsh	<a href="#">1025418</a>	N
01-12-2016 10:30 AM	MR Chest W [71551] - Chest	MR 1 at Cedar	Hilary	Saltmarsh	<a href="#">1025419</a>	N
04-30-2016 1:00 PM	XR Chest Min 4 Views [71030] - Chest	XR 1 at Maple	Hilary	Test	<a href="#">1025537</a>	Y
Reason: Patient Sick						1/1
07-04-2016 11:30 AM	CT Abdomen W [74160] - Abdomen	CT 1 at Cedar	Stick	Stick	<a href="#">1025573</a>	Y
Reason: Room Issue						1/1
11-01-2016 10:50 AM	CT Head Wo [70450] - Head	CT 1 at Maple	Fred	Freddieson	<a href="#">1027870</a>	Y
Reason: Unable to finance						2/2
11-01-2016 11:00 AM	MR Cervical Spine Wo [72141] - Cervical spine	MR 1 at Pine	Harry	Harrison	<a href="#">1027873</a>	Y
07-01-2016 9:00 AM	XR Orbits for Foreign Body [70030] - Head	XR 1 at Pine	Randy	Ziegler	<a href="#">1025563</a>	Y
<b>Total Count</b>						<b>13</b>

Canceled Reasons

Confidential and Proprietary  
Page 1 of 1

Date Run: 11-08-2016 7:11 PM  
RRIS-RSNA\Forms\GuestUser

**Check-in Counts**

Description: This report displays the counts of exam check-ins for the practice(s), site(s) and Date Range selected.

Parameters:

- Practice(s), Site(s) and Date Range based on date patient was checked in.

**Check-in Counts**

**Practice:** Treesdale Radiology

**Site(s):** Cedarville, Mapleton, Pine Grove

**Date Range:** 11-01-2016 - 11-08-2016

**Description:** This report displays the counts of exam check-ins for the practice(s), site(s) and date range selected (Check-in date) order by check-in date descending.

Date	Checked in by	Patient Name	Accession #	Procedures	Status
11-07-2016					1
	<input type="checkbox"/> Eckstein, Rhonda(eckstein)				1
11-04-2016					3
	<input type="checkbox"/> Eckstein, Rhonda(eckstein)				1
	<input type="checkbox"/> Powell, Stick(stick)				1
	<input type="checkbox"/> Saltmarsh, Hilary(saltmarsh)				1
11-01-2016					2
	<input type="checkbox"/> Powell, Stick(stick)				2
		O'Brian, Hillary Demo	1027888	PT/CT Sodium Fluoride PET Bone Scan [PT4]	CheckingIn
		O'Dell, Jayne Demo	1027889	PT/CT Sodium Fluoride PET Bone Scan [PT4]	CheckingIn
<b>Total Count</b>					<b>6</b>

Check-in Counts

Confidential and Proprietary  
Page 1 of 1

Date Run: 11-08-2016 12:57 PM  
RRIS-RSNA\FORMS\GuestUser

### Confirmation Status Canceled or Rescheduled

Description: This report displays all scheduled exams where confirmation status is canceled or rescheduled.

Parameters:

- Practice, Site(s)

### Confirmation Status Canceled or Rescheduled

**Practice(s):** Treesdale Radiology

**Site(s):** All sites for the selected practice(s).

**Date Range:** 11-01-2016 - 11-10-2016

**Description:** This report displays all scheduled exams where confirmation status is canceled or rescheduled

Practice Site	Patient Name	MRN	Scheduled Date	Accession	Procedures	Exam Status	Confirmation Status	Total
<b>Total</b>								0

Confirmation Status Canceled or Rescheduled

Confidential and Proprietary  
Page 1 of 1

Date Run: 11-10-2016 2:45 PM  
RRIS-RSNA\FFormsGuestUser

### Confirmed Appointments

Description: This report displays the confirmed/not confirmed based on practice, and site(s), date range (DOS), confirmation Flag (confirmed\not confirmed), modality type(s), and Confirmation Status(s).

Parameters:

- Practice, Site(s), Date Range (based on DOS), Confirmation Flag (Confirmed vs Not Confirmed appointments), modality types and the confirmation statuses.

#### Confirmed Appointments

Practice: Treesdale

Site(s): CEDARVILLE

Date: 10-31-2016

Modality Type(s): CR, DX, US

Confirmation Status(s): All confirmation statuses.

Description: This report displays the confirmed/not confirmed based on practice, and site(s), date range (DOS), confirmation Flag (confirmed/not confirmed), modality type (s), and Confirmation Status(s).

Acc #	MRN	Patient	Procedures	Procedure Group	Confirmed Status
<b>Total for 10-31-2016:</b>					<b>4</b>
<b>CEDARVILLE</b>					<b>Site Count: 4</b>
1004425556	3446	JAMES, STEVEN	US Abdomen Limited [76705] - Abdomen	US	MessageLeft
1004422734	3444	ROSSI, FRANK	US Abdomen Complete [76700] - Abdomen	US	Confirmed
1004412789	2974	VELLUX, SUE	DX DEXA AXIAL [77080] - DEXA	DEXA	MessageLeft
1004403503	3431	RICHARDS, JANIE	US Pelvic TA and TV [US500] - Pelvis	US	MessageLeft
<b>Total for report</b>					<b>4</b>

Confirmed Appointments

Confidential and Proprietary  
Page 1 of 1

Date Run: 11-08-2016 10:36 PM  
RADNETE:ksteinR

### Daily Schedule

Description: This report displays the details for exams scheduled on the specified date for the practice, site(s) and scanners selected.

Parameters:

- Practice, Site(s), Scanners, and Date (based on scheduled DOS)

### Daily Schedule

**Practice:** Treesdale Radiology

**Site(s):** Cedarville, Mapleton, Pine Grove

**Date:** 11-07-2016

**Scanner(s):** CT1CE, MA1CE, MA2CE, MR1CE, MR2CE, MR3CE, US1CE, US2CE, XR1CE, CT1MP, DE1MP, MA1MP, MR1MP, PT1MP, US2MP, US1MP, XR1MP, CT1PG, DE1PG, MA1PG, MA2PG, MA3PG, MR1PG, MR2PG, NM1PG, NM2PG, PT1PG, US1PG, US2PG, US3PG, XR1PG, XR2PG

**Description:** This report displays the details for exams scheduled on the specified date for the practice, site(s) and scanners selected.

Scanner	Time Sched	Patient Name	MRN#	DOB	Exam(s)
CT1CE	1:05 PM	Strickland, Rhonda	9695	07-17-1983	CT Abdomen W [74160] - Abdomen <b>Body Part\Lat:</b> Abdomen/
<b>Home Phone:</b> (256) 224-5160			<b>Age:</b> 33 y	<b>Referring:</b> Arnold, Bob	
<b>Cell Phone:</b> (256) 224-5160				<b>Copay:</b>	
Current Exam Status: <b>ExamDone</b>				<b>Enc #:</b> 48469	<b>Acc #:</b> 1030889
Insurance: NATIONWIDE INSURANCE				Insurance:	
Scheduler Notes: 11-07-16 eckstein - Informed patient about estimated financial responsibility. 11-07-16 eckstein - send patient back to office					
Chief Complaint: Pain					
Alternate MRN's:					
CT1MP	10:50 AM	Freddieson, Fred	y6	05-09-1973	CT Head Wo [70450] - Head <b>Body Part\Lat:</b> Head/
<b>Home Phone:</b> (902) 374-5967			<b>Age:</b> 43 y	<b>Referring:</b> Flowers, Tim	
<b>Cell Phone:</b>				<b>Copay:</b>	
Current Exam Status: <b>Scheduled</b>				<b>Enc #:</b> 46452	<b>Acc #:</b> 1027871
Insurance:				Insurance:	
Scheduler Notes:					
Chief Complaint: chest pain					
Alternate MRN's:					
CT1PG	9:00 AM	Hartman, Loren	1378	08-02-1995	CT Abdomen W [74160] - Abdomen <b>Body Part\Lat:</b> Abdomen/
<b>Home Phone:</b>			<b>Age:</b> 21 y	<b>Referring:</b> Arnold, Bob	
<b>Cell Phone:</b>				<b>Copay:</b>	
Current Exam Status: <b>Scheduled</b>				<b>Enc #:</b> 47469	<b>Acc #:</b> 1029889
Insurance: PARAMOUNT INSURANCE				Insurance:	
Scheduler Notes:					
Chief Complaint: Pain					
Alternate MRN's:					

### Daily Schedule Add Ons

Description: This report displays the details for exams scheduled on the specified date for the practice, site(s) and scanners selected. As well it filters to pull exams that were added after or including to the specified "added after" date and time.

Parameters:

- Practice, Site(s), Scanners, Date (based on scheduled DOS) and Added After (date and time to see what exams were added after the scheduled date for the same DOS).

### Daily Schedule Add Ons

Practice: Treesdale Radiology

Site(s): Cedarville, Mapleton, Pine Grove

Date: 11-07-2016

Scanner(s): CT1CE, CT1MP, CT1PG, DE1MP, DE1PG, MA1CE, MA1MP, MA1PG, MA2CE, MA2PG, MA3PG, MR1CE, MR1MP, MR1PG, MR2CE, MR2PG, MR3CE, NM1PG, NM2PG, PT1MP, PT1PG, US1CE, US1MP, US1PG, US2CE, US2MP, US2PG, US3PG, XR1CE, XR1MP, XR1PG, XR2PG

Description: This report displays the details for exams scheduled on the specified date for the practice, site(s) and scanners selected. As well it filters to pull exams that were added after or including to the specified "added after" date and time.

Scanner	Time Sched	Patient Name	MRN#	DOB	Exam(s)
CT1CE	1:05 PM	Strickland, Rhonda	9695	07-17-1983	CT Abdomen W [74160] - Abdomen <b>Body Part/Lat:</b> Abdomen/
		<b>Home Phone:</b> (256) 224-5160	<b>Age:</b> 33 years	<b>Referring:</b> Arnold, Bob	
		<b>Cell Phone:</b> (256) 224-5160			<b>Copay:</b>
<b>Current Exam Status:</b> ExamDone			<b>Enc #:</b> 48469	<b>Acc #:</b> 1030889	
Insurance: NATIONWIDE INSURANCE Insurance:					
Scheduler Notes: 11-07-16 eckstein - Informed patient about estimated financial responsibility. 11-07-16 eckstein - send patient back to office					
Chief Complaint: Pain					
Alternate MRN's:					
<i>Added to schedule: 11-07-2016 10:33 AM by eckstein</i>					
					<b>Total Appointments: 1</b>

### Exam Volume Breakdown by Payor

Description: Study count by insurance payor and modality type based on Practice, Site(s), Date Range, and Payor Type.

Parameters:

- Practice, Site(s), Date Range, Payor Type

#### Exam Volume Breakdown by Payor

Practice: Treesdale Radiology

Site(s): Cedarville, Mapleton, Pine Grove

Date Range: 01-01-2012 - 11-08-2016

Payor Types: All payor types selected.

Description: This report displays accession count by insurance payor and modality type based on practice, site(s), date range, and payor type.

Payor Type	CT	MA	MR	US	XR	Accession Count	
<input type="checkbox"/> Blue Cross (P)	1					1	1.49 %
<input type="checkbox"/> HMO (I)	1					1	1.49 %
<input type="checkbox"/> Other (Z)	37	4	7	15	2	65	97.01 %
<b>Total Counts</b>	<b>39</b>	<b>4</b>	<b>7</b>	<b>15</b>	<b>2</b>	<b>67</b>	
	<b>58.21 %</b>	<b>5.97 %</b>	<b>10.45 %</b>	<b>22.39 %</b>	<b>2.99 %</b>		

Exam Volume Breakdown by Payor

Confidential and Proprietary  
Page 1 of 1

Date Run: 11-08-2016 1:12 PM  
RRIS-RSNA\FORMS\GuestUser

#### Exam Volume Breakdown by Payor

Practice: Treesdale Radiology

Site(s): Cedarville, Mapleton, Pine Grove

Date Range: 01-01-2012 - 11-08-2016

Payor Types: All payor types selected.

Description: This report displays accession count by insurance payor and modality type based on practice, site(s), date range, and payor type.

Payor Type	CT	MA	MR	US	XR	Accession Count	
<input type="checkbox"/> Blue Cross (P)	1					1	1.49 %
Blue Cross of CA (090665)	1					1	
<input type="checkbox"/> HMO (I)	1					1	1.49 %
Aetna (213574)	1					1	
<input type="checkbox"/> Other (Z)	37	4	7	15	2	65	97.01 %
AIG INSURANCE (5001270)					1	1	
BROADSPIRE INSURANCE (500606)		1			1	2	
NATIONWIDE INSURANCE (0501224)	1	1				2	
PENN NATIONAL INSURANCE (0501219)	1					1	
PHARMA PLUS MUTUAL INSURANCE (5001117)	7		7		1	15	
STATE FARM INSURANCE (0501185)		1				1	
TRANSGUARD GENERAL INSURANCE (050138)		1				1	
WEST AMERICAN INSURANCE (0501103)	28			14		42	
<b>Total Counts</b>	<b>39</b>	<b>4</b>	<b>7</b>	<b>15</b>	<b>2</b>	<b>67</b>	
	<b>58.21 %</b>	<b>5.97 %</b>	<b>10.45 %</b>	<b>22.39 %</b>	<b>2.99 %</b>		

Exam Volume Breakdown by Payor

Confidential and Proprietary  
Page 1 of 1

Date Run: 11-08-2016 1:12 PM  
RRIS-RSNA\FORMS\GuestUser



## Finding Follow Up

Description: This report displays information on the finding followups by referring and patient for the site(s), period (based on original sign date), findings and follow up reason(s) specified.

Parameters:

- Practice, Site(s), Date Range, Finding(s), Follow Up and Resolved Reason

### Finding Follow Up

**Practice:** Treesdale Radiology

**Site(s):** Cedarville, Mapleton, Pine Grove

**Date Range:** 01-01-2016 - 11-08-2016

**Findings:** Suspicious with other significant finding

**Follow Up Resolved Reason(s):** <Unresolved>

**Description:** This report displays information on the findings followups by referring and patient for the sites, period (based on original sign date), findings and follow up reasons specified.

Referring	Patient	Accession	Finding	Recommendation	Follow Up Reason	Follow Up Type	Follow Up Notes	Date Time
Test, Henry								
Sloth, Gloria MRN: 1501644RSNA								
		1025456	Suspicious with other significant finding	Tissue sampling		Called Patient-1st Attempt	Left message for patient to schedule bx.	03-16-2016 09:51 PM
		1025456	Suspicious with other significant finding	Tissue sampling		Called Referring	Spoke with Nurse Lauren about Bx rec. Will fax order.	06-28-2016 12:28 PM

Finding Follow Up

Confidential and Proprietary  
Page 1 of 1

Date Run: 11-08-2016 01:08 PM  
RRIS-RSNA\FORMS\GuestUser

### Front Desk Activity

Description: This report displays the counts of checkins, scheduled events, # walkins, # verifications for users within the selected sites/practice.

Parameters:

- Practice, Site(s), Scanners, and Date Range (based on created date, check-in date and insurance verification date).

### Front Desk Activity

**Practice:** Treesdale Radiology

**Site(s):** All sites for the selected practice(s).

**Date Range:** 11-01-2016 - 11-08-2016

**Description:** This report displays the counts of checkins, scheduled events, # walkins, # verifications for users within the selected sites/practice.

\*Outside reads are excluded from this report

	Check in	Scheduled	Walk in	Verified
11-01-2016	2	2	2	0
▣ Mapleton	2	2	2	0
11-04-2016	3	4	4	0
▣ Cedarville	3	4	4	0
Eckstein, Rhonda (eckstein)	1	2	2	
Powell, Stick (stick)	1	1	1	
Saltmarsh, Hilary (saltmarsh)	1	1	1	
11-07-2016	1	1	1	0
▣ Cedarville	1	1	1	0
Eckstein, Rhonda (eckstein)	1	1	1	
<b>Total</b>	<b>6</b>	<b>7</b>	<b>7</b>	<b>0</b>

### Image Request History

Description: This report displays the history of Image Requests where requests fall between the specified request dates for the selected practices and locations.

Parameters:

- Practice(s), Site(s), Scanners, and Date (based on image requested date)

### Image Request History

Practice(s): Treesdale Radiology

Site(s): All sites for the selected practice(s).

Date Range: 01-01-2015 - 11-08-2016

Description: This report displays the history of Image Requests where requests fall between the specified request dates for the selected practices and locations.

Requested Date	Completed by	Completed Date	Request Type	Requested By	Image Medium	Delivery Method	Status	MRN	Accession #
<b>Treesdale Radiology</b>									<b>19</b>
<input type="checkbox"/> Cedarville									<b>12</b>
<input type="checkbox"/> 06-30-2016									<b>1</b>
<input type="checkbox"/> 04-19-2016									<b>1</b>
<input type="checkbox"/> 04-06-2016									<b>1</b>
<input type="checkbox"/> 03-24-2016									<b>1</b>
<input type="checkbox"/> 03-23-2016									<b>3</b>
<input type="checkbox"/> 03-11-2016									<b>2</b>
<input type="checkbox"/> 12-22-2015									<b>1</b>
<input type="checkbox"/> 12-17-2015									<b>1</b>
			Referring	George Washington	CD	Courier	Pending	1501325RSNA	1025409
<input type="checkbox"/> 12-16-2015									<b>1</b>
			Referring	George Washington	CD	Courier	Pending	1501325RSNA	1025281
<input type="checkbox"/> Mapleton									<b>7</b>
<input type="checkbox"/> 11-25-2015									<b>7</b>
			Referring	John Adams	CD	Courier	Pending	MASTER_CLONE_1	1023209
			Referring	John Adams	CD	Courier	Pending	z10000001	1027769
			Referring	John Adams	CD	Courier	Pending	z20000002	1027779
			Referring	John Adams	CD	Courier	Pending	z30000003	1027789
			Referring	John Adams	CD	Courier	Pending	z40000004	1027799
			Referring	John Adams	CD	Courier	Pending	z50000005	1027809
			Referring	John Adams	CD	Courier	Pending	z50000006	1027819
<b>Total Requests:</b>									<b>19</b>

### IVT 7 Day Trend

Description: This report displays the counts of IVT Approvals over the 7 days preceding the selected parameter based on chosen practice/sites.

Parameters:

- Practice, Site(s) and Date


### IVT 7 Day Trend

Practice: Treesdale

Site(s): CEDARVILLE

Date Range: 11-01-2016 - 11-08-2016

Description: This report displays the counts of IVT Approvals over the 7 days preceding the selected parameter based on chosen practice/sites.



IVT 7 Day Trend

Confidential and Proprietary  
Page 1 of 1

Date Run: 11-08-2016 10:29 PM  
RADNETecksteinR

## IVT Idle Time

Description: This report displays the time an exam sits in the IVT queue before worked based on the difference between the study created and IVT requested date.

Parameters:

- Practice and Site(s)

## IVT Idle Time

**Practice:** Treesdale Radiology

**Site(s):** Cedarville

**Description:** This report displays the time an exam sits in the IVT queue before worked based on the difference between the study created and IVT requested date.

Patient Name	IV status	Pre Cert Status	Minutes in Queue	Total
Acosta,Hannah	Required		295647	
Adams,Judith	N		351149	
	Required		505362	
			1585156	
			2085153	
			2612003	
Adamson,Adam	Required		2116701	
Aguilar,Annie	Required		152979	
			2116872	
Alberg,Tijuana	Required		2610967	
Ali,Amanda	Required		354845	
Alvarez,Bridgette	Required		297274	
Anthony,Anna	Required		2116701	
Ashley,Beth	Required		281722	
Austin,Geoffrey	Required		191425	
Bean,Gene	Required		1988382	
Beatton,Betsy	N		505538	
Bell,Bart	Required		2128381	
Bonilla,Jim	Required		152289	
Borne,Bernie	Required		2018774	
			2018951	
Brilhante,Dora	Required		2611909	
Brones,Fantasia	Required		343387	
Browning,Fernando	Required		152287	
Bueche,Brittany	Required		152996	

## IVT Turnaround

Description: This report displays the report turnaround times starting at requested date to insurance verification approved date for the selected Practice/Site(s) and date range (based on IVT approval date).

Parameters:

- Practice, Site(s) and Date

## IVT Turnaround

**Practice:** Treesdale Radiology

**Site(s):** All sites for the selected practice(s).

**Date Range:** 01-01-2015 - 11-08-2016

**Description:** This report displays the report turnaround times starting at requested date to insurance verification approved date for the selected Practice/Site(s) and date range (based on date IVT approval was made).

Site	Patient Name	Requested Date	Approved Date	Total Turnaround (hrs:mi)	IV Total
Cedarville	Ibarra, Ira	11-01-2012 10:19 AM	12-15-2015 3:56 PM	27342:37	1
	Igorson, Ivan	11-01-2016 10:19 AM	12-15-2015 3:56 PM	-7716:23	1
	Matina, May	05-03-2016 7:26 PM	05-03-2016 7:30 PM	00:04	1
	Saltmarsh, Hilary	06-30-2016 8:46 AM	06-30-2016 1:38 PM	04:52	1
	Test, Hilary	05-14-2016 12:09 AM	05-14-2016 12:15 AM	00:06	1
Mapleton	Maltolia, Ratna	05-03-2016 7:32 PM	05-03-2016 7:34 PM	00:02	1
	Mangia, Gertrude	05-03-2016 7:14 PM	05-03-2016 7:26 PM	00:12	1
<b>Total</b>					<b>7</b>

IVT Turnaround

Confidential and Proprietary  
Page 1 of 1

Date Run: 11-08-2016 7:18 PM  
RRIS-RSNA\FORMS\GuestUser

### Modality Duration Override

Description: This report lists modality duration overrides for the selected location(s).

Parameters:

- Practice(s), and site(s)

### Modality Duration Override

Practice(s): Tressdale

Site(s): CEDARVILLE

Description: This report lists modality duration overrides for the selected location(s).

Modality	Procedure Code	Override Duration	Regular Duration
<input checked="" type="checkbox"/> MRI1			
	MR703	60	30
<input checked="" type="checkbox"/> US1			
	93930	45	30
	93970	45	30
	US400	45	30
	US526	45	30
	US527	45	30
	US913	45	30
<input checked="" type="checkbox"/> US2			
	93930	45	30
	93970	45	30
	US400	45	30
	US526	45	30
	US527	45	30
	US913	45	30
<input checked="" type="checkbox"/> US3			
	93930	45	30
	93970	45	30
	US400	45	30
	US526	45	30
	US527	45	30
	US913	45	30

**MR Breast Dictated**

Description: This report displays the studies which include procedures descriptions starting with MR Breast within the selected sites/practice and date range (based on original dictation date).

Parameters:

- Practice, Site(s), Date Range and Procedure Code(s)

**MR Breast Dictated**

**Practice(s):** Treesdale Radiology

**Site(s):** All sites for the selected practice(s).

**Date Range:** 01-01-2015 - 11-08-2016

**Procedure Codes included:** 77058, 77059, MR68, MR69, MR77, MR78, MR79, MR80

**Description:** This report displays the studies which include procedures descriptions starting with MR Breast within the selected sites/practice and date range (based on original dictation date).

Site	Patient Name	DOS	Dictated	MRN	Accession
<b>Cedarville</b>					<b>1</b>
	Strickland, Rhonda	11-26-2015	10-11-2016	1501529RSNA	1023939
	Procedure(s): MR Breast Biopsy Bilateral [MR78]				
<b>Total</b>					<b>1</b>

MR Breast Dictated

Confidential and Proprietary  
Page 1 of 1

Date Run: 11-08-2016 7:19 PM  
RRIS-RSNA\FORMS\GuestUser



### Need Plan

Description: This report displays exams needing an insurance plan assigned based on selected practice and site(s), excluding the selected ignored status(s). \*\*Need plan is identified as selected carrier code of 99999 (description of 'Need-Plan') configured in the carrier table of RIS.

Parameters:

- Practice, Site(s) and Ignored status(s)

### Need Plan

Practice: Treedale  
 Site(s): CEDARVILLE

Description: This report displays exams needing an insurance plan assigned based on selected practice and site(s). Excluding the selected ignored status(s).  
 \*\*Need plan is identified as selected carrier code of 99999 (description of 'Need Plan') configured in the carrier table of RIS.

Patient Name	MRN	Procedures	Accession	Status
<b>CEDARVILLE</b>				<b>Site Total :</b>
<b>10-01-2014</b>			<b>Date Total :</b>	<b>1</b>
	52901	XR Hand Minimum 3 Views [73130] - Hand - R	10001234	Signed1
<b>10-06-2014</b>			<b>Date Total :</b>	<b>1</b>
	54061	XR Lumbar Spine 2 or 3 Views [72100] - Lumbar spine	10003456	Signed1
<b>10-08-2014</b>			<b>Date Total :</b>	<b>1</b>
	5470	XR Osseous Survey Limited [77074] - Body	10005412	Signed1
<b>10-09-2014</b>			<b>Date Total :</b>	<b>1</b>
	55012	XR Tibia and Fibula 2 Views [73590] - Tib-Fib - L	10006523	Signed1
<b>10-21-2014</b>			<b>Date Total :</b>	<b>1</b>
	57934	XR Hand Minimum 3 Views [73130] - Hand - R	10008523	Signed1
<b>09-25-2015</b>			<b>Date Total :</b>	<b>1</b>
	17496	XR Hand Minimum 3 Views [73130] - Hand - R		Signed1
<b>Grand Total</b>				<b>6</b>

Need Plan

Confidential and Proprietary  
 Page 1 of 1

Date Run: 11-10-2016 5:24 PM  
 RADNETE:ksteinR

**NM Therapy Patient Listing**

Description: This report shows all exams within the specified practices/sites and date range where the procedure description contains I-131 (based on scheduled DOS).

Parameters:

- Practice(s), Site(s) and Date Range based on scheduled DOS

**NM Therapy Patient Listing**

**Practice:** Treesdale Radiology

**Site(s):** All sites for the selected practice(s).

**Date Range:** 11-09-2016 - 11-10-2016

**Description:** This report shows all exams within the specified practices/sites and date range where the procedure description contains I-131 (based on scheduled DOS).

Scheduled	Patient Name	MRN #	DOB	Encounter ID
11-10-2016	Becker, Julia	5523	06-09-1967	48479
<b>Exam:</b> NM Thyroid Mult Uptake and Scan_I-131 [78007] - Neck			<b>Location:</b> Pine Grove	
<b>Status:</b> Signed1				
<b>Referring:</b> Dr Arnold Bob				
11-10-2016	Becker, Julia	5523	06-09-1967	48480
<b>Exam:</b> NM Thyroid Whole Body w/i-131 [78018B]			<b>Location:</b> Pine Grove	
<b>Status:</b> Signed1				
<b>Referring:</b> Dr Arnold Bob				

NM Therapy Patient Listing

Confidential and Proprietary  
Page 1 of 1

Date Run: 11-10-2016 3:11 PM  
RRIS-RSNA\FORMS\GuestUser

### No Shows

Description: This report displays exams that were never arrived for the selected Practice\Site(s), date range (based on DOS) and procedure groups.

Parameters:

- Practice, Site(s), Date Range based on DOS

### No Shows

**Practice:** Treesdale Radiology

**Site(s):** All sites for the selected practice(s).

**Date Range:** 01-01-2016 - 11-08-2016

**Procedure Group(s):** CT, DE, EM, FL, MA, MR, NM, OT, PE, PT, RD, RO, RS, SR, UG, UV, XR

**Description:** This report displays exams that were never arrived for the selected Practice\Site(s), date range (based on DOS) and procedure groups.

Site	CT	MA	MR	PT	UG	XR	No Show/Total	%
Cedarville	13	4	8	0	28	5	<b>58/121</b>	47.9%
Mapleton	3	2	3	2	4	2	<b>16/29</b>	55.2%
Outside Read1	2	0	0	0	0	0	<b>2/2</b>	100.0%
Pine Grove	3	0	2	0	1	0	<b>6/9</b>	66.7%
<b>No Show/Total Exams</b>	<b>21/55</b>	<b>6/23</b>	<b>13/24</b>	<b>2/2</b>	<b>33/44</b>	<b>7/13</b>	<b>82/161</b>	
<b>% No Show</b>	<b>38.2 %</b>	<b>26.1 %</b>	<b>54.2 %</b>	<b>100.0 %</b>	<b>75.0 %</b>	<b>53.8 %</b>	<b>50.9 %</b>	<b>50.9 %</b>

No Shows

Confidential and Proprietary  
Page 1 of 1

Date Run: 11-08-2016 7:42 PM  
RRIS-RSNA\FORMS\GuestUser

### Open Capacity

Description: This report displays the open capacity on rooms (open, scheduling and available minutes) filtered by practice, site, and date.

Parameters:

- Practice, Sites and date

### Open Capacity

**Practice:** Treesdale Radiology

**Site(s):** All sites for the selected practice(s).

**Date:** 11-07-2016

**Description:** This report displays the open capacity on rooms (open, scheduling and available minutes) filtered by practice, site, and date. The room must have a template associated in order to display on this report.

Site	Room	Room Times	Minutes Open	Minutes Scheduled	Minutes Closed	Minutes Available
Cedarville	Total		960	0	0	960
	MR1CE	Total	480	0	0	480
		8:00 AM - 12:00 PM	240	0	0	480
	MR3CE	Total	480	0	0	480
		8:00 AM - 12:00 PM	240	0	0	480
	<b>Total</b>			<b>960</b>	<b>0</b>	<b>0</b>

Open Capacity

Confidential and Proprietary  
Page 1 of 1

Date Run: 11-08-2016 7:46 PM  
RRIS-RSNA\FORMS\GuestUser

### OTC QA

Description: This report shows all payment details for the selected date by patient within a selected Practice, grouped by site and patient.

Parameters:

- Practice, Sites and date

### OTC QA

Practice: Treesdale Radiology

Site(s): Cedarville

Date: 11-10-2016

Report Description: This report shows all payment details for the selected date by patient within a selected Practice, grouped by site and patient.

\*Note: ? indicates a copay was entered as a note, but not entered into the copay field within RIS.

Site	Patient Name	Potential OTC	Amount Collected	Copay Owning	Receipt #	Payment Entered by	Payment Notes
Cedarville	Site Total	\$ 185.00	\$ 185.00				
	Hall, Josh MRN: 3233 Arrived by:	\$ 15.00	\$ 15.00				
			\$ 15.00		Cedar-10	eckstein	Cash
	James, Sally MRN: 1184 Arrived by:	\$ 50.00	\$ 50.00				
			\$ 50.00		Cedar-8	eckstein	Check#1232
	Rodriguez, Rhonda MRN: 4454 Arrived by:	\$ 120.00	\$ 120.00				
			\$ 120.00		Cedar-9	eckstein	CreditCard #4442 exp:12/16 cardholder:Rodriguez, Rhonda
	<b>Total</b>	<b>\$ 185.00</b>	<b>\$ 185.00</b>				

Payment Method	Amount Collected
Cash	\$ 15.00
Check	\$ 50.00
CreditCard	\$ 120.00
<b>Total</b>	<b>\$ 185.00</b>

OTC QA

Confidential and Proprietary  
Page 1 of 1

Date Run: 11-10-2016 10:26 AM  
RRIS-RSNA\FORMS\GuestUser

### Patient Email and Cellphone Percentages

Description: This report shows the % of collected cell phone and email addresses for patients who've had studies within the specified practices, sites and date range. Dates are based on when the study was created rather than date of service.

Note: These numbers could change if a user were to update the patient cellphone or email addresses between the time this report was run and current date.

Parameters:

- Practice, Sites and Date Range based on created date of the study

#### Patient Email and Cellphone Percentages

Practice: Treedale Radiology

Site(s): All sites for the selected practice(s).

Date Range: 06-01-2016 - 11-08-2016

Description: This report shows the % of collected cell phone and email addresses for patients who've had studies within the specified practices, sites and date range. Dates are based on when the study was created rather than date of service.

Note: These numbers could change if a user were to update the patient cellphone or email addresses between the time this report was run and current date.

Practice/Site	Jun 2016			Jul 2016			Aug 2016			Sep 2016			Oct 2016			Nov 2016			Total		
	% with Email	% with Cell	Patient Count	% with Email	% with Cell	Patient Count	% with Email	% with Cell	Patient Count	% with Email	% with Cell	Patient Count	% with Email	% with Cell	Patient Count	% with Email	% with Cell	Patient Count	% with Email	% with Cell	Patient Count
Treedale Radiology	14%	29%	7	9%	64%	11	0%	50%	4	0%	33%	3	0%	67%	3	38%	63%	8	14%	53%	36
<b>Total</b>	<b>14%</b>	<b>29%</b>	<b>7</b>	<b>9%</b>	<b>64%</b>	<b>11</b>	<b>0%</b>	<b>50%</b>	<b>4</b>	<b>0%</b>	<b>33%</b>	<b>3</b>	<b>0%</b>	<b>67%</b>	<b>3</b>	<b>38%</b>	<b>63%</b>	<b>8</b>	<b>14%</b>	<b>53%</b>	<b>36</b>

Patient Email and Cellphone Percentages

Confidential and Proprietary  
Page 1 of 1

Date Run: 11-08-2016 7:49 PM  
RRIS-RSNA FormsGuestUser

### Patient Wait Time Trending by Site

Description: This report displays the patient wait times for sites within the selected practice and date range. Broken down into percentages less than 5 minutes, between 5 & 20, and finally greater than a 20 minute wait.

Parameters:

- Practice, Sites and date

### Patient Wait Time Trending by Site

**Practice:** Treesdale Radiology

**Site(s):** All sites for the selected practice(s).

**Date Range:** 09-01-2016 - 11-08-2016

**Description:** This report displays the patient wait times for sites within the selected practice and date range. Broken down into percentages less than 5 minutes, between 5 & 20, and finally greater than a 20 minute wait.

	Cedarville				Mapleton			
	<5	between	>20	# Exams	<5	between	>20	# Exams
<b>Nov 2016</b>	100.00 %	0.00 %	0.00 %	4	100.00 %	0.00 %	0.00 %	7
<b>Oct 2016</b>	100.00 %	0.00 %	0.00 %	1	0.00 %	0.00 %	0.00 %	0
<b>Sep 2016</b>	0.00 %	0.00 %	100.00 %	2	0.00 %	0.00 %	0.00 %	0
<b>Total</b>	71.43 %	0.00 %	28.57 %	7	100.00 %	0.00 %	0.00 %	7

Exam Totals By Room		Exam Counts
<input type="checkbox"/>	<b>Cedarville</b>	6
	CT1CE	3
	MR1CE	1
	US1CE	1
	XR1CE	1
<input type="checkbox"/>	<b>Mapleton</b>	7
	CT1MP	7
<b>Total</b>		13

Patient Wait Time Trending by Site

Confidential and Proprietary  
Page 1 of 1

Date Run: 11-08-2016 7:50 PM  
RRIS-RSNA\FORMS\GuestUser

### Patient Wait Times

Description: This report details the average patient wait times by room type for the selected practice, site(s). Negative numbers indicating the patient was early are resolved to 0 for averaging.

Parameters:

- Practice, Sites and date

### Patient Wait Times

Practice: Treesdale Radiology

Site(s): All sites for the selected practice(s).

Date Range: 09-01-2016 - 11-08-2016

Description: This report details the average patient wait times by room type for the selected practice,site(s). Negative numbers indicating the patient was early are resolved to 0 for averaging.

Site	Accession	CT	MA	MR	US	XR	Average
<input type="checkbox"/> Cedarville		49.33	11355.00	298.00	184.00	181.00	1738.00
	<a href="#">1025596</a>		11355.00				
	<a href="#">1025601</a>			298.00			
	<a href="#">1027404</a>				184.00		
	<a href="#">1028889</a>	143.00					
	<a href="#">1028891</a>					181.00	
	<a href="#">1029890</a>	5.00					
	<a href="#">1030889</a>	0.00					
<input type="checkbox"/> Mapleton		164.00					164.00
<b>Totals</b>		<b>129.60</b>	<b>11355.00</b>	<b>298.00</b>	<b>184.00</b>	<b>181.00</b>	<b>951.00</b>

Patient Wait Times

Confidential and Proprietary  
Page 1 of 1

Date Run: 11-08-2016 7:52 PM  
RRIS-RSNA\FORMS\GuestUser

### Accession Milestone Times

Patient: 1501615RSNA - Test, Hilary

Current Exam Status: **Started**

Site: Cedarville

Accession #	Procedure Category	Procedure Description
1027404	UG	US Abdominal Complete [76700] - Abdomen

Exam Milestone Times					
Order Date	Scheduled Date	Checking In Date	Arrival Date	Started Date	Performed End Date
09-13-2016 3:02 PM	10-28-2016 3:49 PM	10-28-2016 6:49 PM	10-28-2016 6:49 PM	10-28-2016 6:53 PM	
saltmarsh	saltmarsh	saltmarsh	saltmarsh	janet	

Report Milestone Times				
Version	Dictated date	Drafted Date	Edited Date	Signed Date

Turnaround Times (in Hrs)				
**calculated at 8:00-6:00 Mon-Fri				
Done To Dictated	Dictated To Drafted	Drafted To Edited	Edited To Signed	Done To Signed
0.00	0.00	0.00	0.00	0.00

Turnaround Times (in Hrs)				
**calculated at 24 / 7				
Done To Dictated	Dictated To Drafted	Drafted To Edited	Edited To Signed	Done To Signed



**Patients without Insurance**

Description: This report lists the exams where no insurance was identified for the selected practice(s), site(s), and date range.

Parameters:

- Practice(s), Site(s), Date Range based on DOS

**Patients Without Insurance**

Practice: Treeddale Radiology

Site(s): Cedarville

Date Range: 10-01-2016 - 11-08-2016

Description: This report lists the exams where no insurance was identified for the selected practice(s), site(s), and date range.

Appointment Date	Patient	MRN	Procedure Code	Procedure Description	Appt Time
<b>10-19-2016</b>					<b>Count:</b> 1
	Rhonda Rodriguez	4454	71260	CT Chest W [71260] - Chest	10:40 AM
	Site: Cedarville	Room:	CT 1 at Cedar	Accession: 1026702	
	Chief Complaint: pain Order Notes:				
<b>10-28-2016</b>					<b>Count:</b> 2
	HILARY TESTEROSA	1501663RSNA	CT22	CT Head W & Sinus W [CT22] - Head	3:35 PM
	Site: Cedarville	Room:	CT 1 at Cedar	Accession: 1027403	
	Chief Complaint: 10-28-16 saltmarsh - good news? uhoh Order Notes:				
	Hilary Test	1501615RSNA	76700	US Abdominal Complete [76700] - Abdomen	3:49 PM
	Site: Cedarville	Room:	US 1 at Cedar	Accession: 1027404	
	Chief Complaint: text gone from when this was log controlled Order Notes:				
<b>11-04-2016</b>					<b>Count:</b> 3
	Hilary Parsons	8211	71030	XR Chest Min 4 Views [71030] - Chest	8:47 AM
	Site: Cedarville	Room:	XR 1 at Cedar	Accession: 1028891	
	Chief Complaint: Tightness in chest Order Notes:				
	Rhonda Rodriguez	4454	72126	CT Cervical Spine W [72126] - Cervical spine	11:00 AM
	Site: Cedarville	Room:	CT 1 at Cedar	Accession: 1028890	
	Chief Complaint: Pain in neck Order Notes:				
	Jayne Demo O'Dell	z40000004	70480	CT Ear Orbit Sella Wo [70480] - Head	2:20 PM
	Site: Cedarville	Room:	CT 1 at Cedar	Accession: 1029890	
	Chief Complaint: test Order Notes:				

### Payments Received

Description: This report shows the payments collected by Date and Site(s) selected. Grouped by Site and payment type.

Parameters:

- Practice(s), Site(s), and payment date

#### Payments Received

Practice: Treeddale Radiology

Site(s): Cedarville, Mapleton, Pine Grove

Payment Date: 11-10-2016

Report Description: This report shows the payments collected by Date and Site(s) selected. Grouped by Site, payment type.

Encounter	Patient	MRN	Note	CHK/CC #	Rcpt #	Payment	Void	Net Payment
<b>Cedarville</b>						\$ 185.00		\$ 185.00
<b>Cash</b>						\$ 15.00		\$ 15.00
48478	Josh Hall Carrier: PHARMA PLUS MUTUAL INSURANCE	3233			Cedar-10	\$ 15.00		\$ 15.00
<b>Check</b>						\$ 50.00		\$ 50.00
48473	Sally James Carrier: PHARMA PLUS MUTUAL INSURANCE	1184		Chk#: 1232	Cedar-8	\$ 50.00		\$ 50.00
<b>Credit Card</b>						\$ 120.00		\$ 120.00
48471	Rhonda Rodriguez Carrier: PHARMA PLUS MUTUAL INSURANCE	4454		CC# 4442	Cedar-9	\$ 120.00		\$ 120.00
<b>Grand Total</b>						\$ 185.00		\$ 185.00

Payments Received

Confidential and Proprietary  
Page 1 of 1

Date Run: 11-10-2016 4:09 PM  
RRIS-RSNAVFormsGuestUser

**Pre-Auth Referring Physician**

Description: This report displays the number of orders placed by pre-auth referring physicians broken down by modality type for the selected Practice/Site(s) and date range.

Parameters:

- Practice(s), Site(s), Date Range and Referring Physicians

**Pre-Auth Referring Physician**

Practice: Treesdale

Site(s): Cedarville

Date Range: 09-01-2016 - 11-10-2016

Description: This report displays the number of orders placed by pre-auth referring physicians broken down by modality type for the selected Practice/Site(s) and date range.

Referring Physician	Procedure Group	Total
EARNOLD, MD, BOB	Total	2
	CT	1
	MR	1
<b>Grand Total</b>		<b>2</b>

Pre-Auth Referring Physician

Confidential and Proprietary  
Page 1 of 1

Date Run: 11-10-2016 9:57 AM  
RADNETEcksteinR

### Reschedules and Cancelations

Description: This report displays the reschedules, cancelations and discontinues by users selected by practice, site

Parameters:

- Practice, Site(s) and Date Range

### Reschedules and Cancelations

**Practice:** Treesdale Radiology

**Site(s):** All sites for the selected practice(s).

**Date Range:** 08-01-2016 - 11-08-2016

**Description:** This report displays the reschedules, cancelations and discontinues by users selected by practice, site.

Site	User	Cancelations	Reschedules	Total
<b>Cedarville</b>		<b>2</b>		<b>2</b>
	(stick) Powell, Stick	1		1
	(saltmarsh) Saltmarsh, Hilary	1		1
<b>Mapleton</b>		<b>1</b>	<b>1</b>	<b>2</b>
	(stick) Powell, Stick	1	1	2
<b>Pine Grove</b>			<b>1</b>	<b>1</b>
	(stick) Powell, Stick		1	1
<b>Grand Total</b>		<b>3</b>	<b>2</b>	<b>5</b>

Reschedules and Cancelations

Confidential and Proprietary  
Page 1 of 1

Date Run: 11-08-2016 8:01 PM  
RRIS-RSNA\FORMS\GuestUser

## Referring Physician Details

Description: This report displays the exams requested by a physician within the dates selected based on date of service. The report is grouped by date the physicians' requested the study.

Parameters:

- Date Range (DOS), last name, first name, middle name (optional) and physician selected
- Note: entry of the last name, first and optionally the middle name, filters the parameter to select the referring you wish to report on. The last name must match the whole last name, first and middle name can be a starts with.

## Referring Physician Details

Date Range: 01-01-2015 - 11-08-2016

Physician name searched on: arnold, b

Physician selected: Arnold, Bob ID:17872

Description: This report displays the exams requested by a physician within the dates selected based on date of service. The report is grouped by the date the physicians' requested the study.

Request Date	Accession	Patient	Scheduled Date	Status	
10-19-2016					1
	1026702	(4454) - Rodriguez, Rhonda CT Chest W [71260] - Chest	10-19-2016 10:40 AM	Scheduled	
11-03-2016					1
	1028888	(1378) - Hartman, Loren CT Abdomen W [74160] - Abdomen	11-04-2016 1:30 PM	Scheduled	
11-04-2016					5
	1028891	(8211) - Parsons, Hilary XR Chest Min 4 Views [71030] - Chest	11-04-2016 8:47 AM	Signed1	
	1028889	(9695) - Strickland, Rhonda CT Head W [70460] - Head	11-04-2016 9:55 AM	Signed1	
	1028890	(4454) - Rodriguez, Rhonda CT Cervical Spine W [72126] - Cervical spine	11-04-2016 11:00 AM	Scheduled	
	1029890	(z40000004) - O'Dell, Jayne CT Ear Orbit Sella Wo [70480] - Head	11-04-2016 2:20 PM	Signed1	
	1029889	(1378) - Hartman, Loren CT Abdomen W [74160] - Abdomen	11-07-2016 9:00 AM	Scheduled	
11-07-2016					1
	1030889	(9695) - Strickland, Rhonda CT Abdomen W [74160] - Abdomen	11-07-2016 1:05 PM	ExamDone	
<b>Total</b>					<b>8</b>

Referring Physician Details

Confidential and Proprietary  
Page 1 of 1

Date Run: 11-08-2016 8:04 PM  
RRIS-RSNA\FORMS\GuestUser

**Room Closures**

Description: This report displays the amount of time that a room was closed as well as the Closure Reasons.

Parameters:

- Practice, Site(s), Date Range and Room(s) selected

**Room Closures**

**Practice:** Treesdale Radiology

**Site(s):** All sites for the selected practice(s).

**Date Range:** 10-01-2014 - 11-08-2016

**Description:** This report displays the amount of time that a room was closed as well as the Closure Reasons.

Room	Closure Reason	From Date\Time	To Date\Time	Time Closed
<b>MA1CE</b>				<b>2h 0m</b>
	Room Closed for QC	04-18-2016 10:00 AM	04-18-2016 12:00 PM	2h 0m
<b>Total</b>				<b>2h 0m</b>

Room Closures

Confidential and Proprietary  
Page 1 of 1

Date Run: 11-08-2016 8:07 PM  
RRIS-RSNA\FORMS\GuestUser

**Room Utilization**

Description: This report displays minutes open, closed, scheduled and remaining for rooms within the selected Practice, sites and rooms. Outside read locations are excluded from this report

Parameters:

- Practice, Site(s) and Date

**Room Utilization**

**Practice:** Treesdale Radiology

**Site(s):** All sites for the selected practice(s).

**Date:** 11-09-2016

**Description:** This report displays minutes open, closed, scheduled and remaining for rooms within the selected Practice, sites and rooms.

*\*Outside read locations are excluded from this report*

Site	Room	Open	Closed	Scheduled	Remaining	% Used
☐ Cedarville						7.69 %
	MR1CE	780	0	120	660	15.38 %
	MR3CE	780	0	0	780	0.00 %
						<b>7.69 %</b>

Room Utilization

Confidential and Proprietary  
Page 1 of 1

Date Run: 11-08-2016 8:15 PM  
RRIS-RSNA\FORMS\GuestUser

## RTAT by Radiologist

Description: This report shows the Report Turn Around Times by Radiologist based between Pending Dictation to Signed.

Parameters:

- Practice(s), Site(s), Date Range, Procedure Group(s) and Radiologist(s)

## RTAT by Radiologist

**Practice(s):** Treesdale Radiology

**Site(s):** All sites for the selected practice(s).

**Date Range:** 01-01-2016 - 11-08-2016

**Procedure Groups:** Diagnostic CT,DEXA,Evaluation Mgmt,Fluoroscopy,Mammography,MRI and MRA,Nuclear Med,OTHER,EPET Procedures,PET,Diagnostic Imagt,Radiation Oncol,Intervntnl Rad,Intervntnl Surg,Ultrasound Gen,Ultrasound Vasc,Plain Film/Xray/CR

**Excluding exams signed on weekends**

Pending Dictation - Pending Transcription (% of reports completed in < 3 hours)		
Radiologist	Tree	# Signed Exams
<input checked="" type="checkbox"/> Eckstein, Rhonda	100.00 %	1
<input type="checkbox"/> Powell, Stick	100.00 %	1
	<a href="#">1029890</a>	0.02
<input checked="" type="checkbox"/> Saltmarsh (Rad), Hilary	66.67 %	3
<input checked="" type="checkbox"/> Saltmarsh, Hilary	64.29 %	14
<b>Total</b>	<b>68.42 %</b>	<b>19</b>

Pending Transcription - Pending Signature (% of reports completed in <= 8 hours)		
Radiologist	Tree	# Signed Exams
<input checked="" type="checkbox"/> Eckstein, Rhonda	100.00 %	1
<input checked="" type="checkbox"/> Powell, Stick	100.00 %	1
<input type="checkbox"/> Saltmarsh (Rad), Hilary	100.00 %	3
	<a href="#">1025466</a>	0.00
	<a href="#">1025456</a>	0.00
	<a href="#">1025458</a>	0.00
<input checked="" type="checkbox"/> Saltmarsh, Hilary	100.00 %	14
<b>Total</b>	<b>100.00 %</b>	<b>19</b>

Pending Signature - Signed (% of reports completed in <= 8 hours)		
Radiologist	Tree	# Signed Exams
<input type="checkbox"/> Eckstein, Rhonda	100.00 %	1
	<a href="#">1028889</a>	0.00
<input checked="" type="checkbox"/> Powell, Stick	100.00 %	1
<input checked="" type="checkbox"/> Saltmarsh (Rad), Hilary	100.00 %	3
<input checked="" type="checkbox"/> Saltmarsh, Hilary	100.00 %	14



### RTAT By Referring

Description: This report displays the RTAT for physician specified for the selected Practice/Site(s) and period.

\*RTAT is calculated based on Exam Performed to Signed for the business hours between 8:00am through 6:00pm

Parameters:

- Practice, Site(s), Date Range, first and last name of referring

### RTAT By Referring

Practice: Treesdale Radiology

Site(s): All sites for the selected practice(s).

Date Range: 01-01-2016 - 11-08-2016

Description: This report displays the RTAT for physician specified for the selected Practice/Site(s) and period.

\*RTAT is calculated based on Exam Performed to Signed for the business hours between 8:00am through 6:00pm.

Requested By	Scheduled Date	Details	TAT (hrs)*
<input checked="" type="checkbox"/> Dr Bob Arnold		Average for period	0.02
	11-04-2016	Patient: (8211) Hilary Parsons Procedure: XR Chest Min 4 Views [71030] - Chest Accession: <a href="#">1028891</a> Site: Cedarville	0.02
	11-04-2016	Patient: (9695) Rhonda Strickland Procedure: CT Head W [70460] - Head Accession: <a href="#">1028889</a> Site: Cedarville	0.02
	11-04-2016	Patient: (z40000004) Jayne Demo O'Dell Procedure: CT Ear Orbit Sella Wo [70480] - Head Accession: <a href="#">1029890</a> Site: Cedarville	0.02
Total Records			3

RTAT By Referring

Confidential and Proprietary  
Page 1 of 1

Date Run: 11-08-2016 8:21 PM  
RRIS-RSNA\FORMS\GuestUser

**RTAT Detail**

Description: This report shows the Report Turn Around Times for the practice, and sites selected between Pending Dictation to Signed.

Parameters:

- Practice(s), Site(s), Date Range, and Procedure Group(s)

**RTAT Detail**

**Practice(s):** Treesdale Radiology

**Site(s):** All sites for the selected practice(s).

**Date Range:** 10-01-2016 - 11-08-2016

**Procedure Groups:** (DE) DEXA,(CT) Diagnostic CT,(RD) Diagnostic Imagt,(PE) EPET Procedures,(EM) Evaluation Mgmt,(FL) Fluoroscopy,(RS) Intervntnl Rad,(SR) Intervntnl Surg,(MR) MRI and MRA,(NM) Nuclear Med,(OT) OTHER,(PT) PET,(XR) Plain Film/Xray/CR,(RO) Radiation Oncol,(UG) Ultrasound Gen,(UV) Ultrasound Vasc

**All Modalities Except Mammo**

**Excluding exams signed on weekends**

Radiologist: Pending Dictation - Pending Transcription (% of reports completed in < 3 hours)		
Week Ending	Tree	# Signed Exams
11-04-2016	0.00 %	3
<b>Total</b>	<b>0.00 %</b>	<b>3</b>

Transcriptionist: Pending Transcription - Pending Signature (% of reports completed in < 3 hours)		
Week Ending	Tree	# Signed Exams
11-04-2016	100.00 %	3
<b>Total</b>	<b>100.00 %</b>	<b>3</b>

Radiologist: Pending Signature - Signed (% of reports completed in < 2 hours)		
Week Ending	Tree	# Signed Exams
11-04-2016	100.00 %	3
<b>Total</b>	<b>100.00 %</b>	<b>3</b>

**RTAT Detail - Monthly**

Description: This report shows the Report Turn Around Times for the practice(s), site(s) and month based between Pending Dictation to Signed.

Parameters:

- Practice(s), Site(s), Date Range, and Month ending

**RTAT Detail - Monthly**

**Practice(s):** Treesdale Radiology  
**Site(s):** All sites for the selected practice(s).  
**Date Range:** 11-01-2016 - 11-08-2016  
**All Modalities Except Mammo**  
**Excluding exams signed on weekends**

Practice	Accession #	Radiologist: Pending Dictation - Pending Transcription				Radiologist: Pending Signature - Signed					
		Count	< 1.5 hrs	1.5 - 3 hrs	3 - 6 hrs	> 6 hrs	Count	< 1 hr	1 - 2 hrs	2 - 4 hrs	> 4 hrs
☐ Tree		3	100.00 %	0.00 %	0.00 %	0.00 %	3	100.00 %	0.00 %	0.00 %	0.00 %
	1028889	0.02	1				0	1			
	1028891	0.02	1				0	1			
	1029890	0.02	1				0	1			
<b>Total</b>		<b>3</b>	<b>100.00 %</b>	<b>0.00 %</b>	<b>0.00 %</b>	<b>0.00 %</b>	<b>3</b>	<b>100.00 %</b>	<b>0.00 %</b>	<b>0.00 %</b>	<b>0.00 %</b>

RTAT Detail - Monthly

Confidential and Proprietary  
 Page 1 of 1

Date Run: 11-08-2016 8:25 PM  
 RRIS-RSNAVFormsGuestUser

**RTAT Detail Mammo**

Description: This report shows the Mammography specific Report Turn Around Times for the selected practice(s), site(s), date range and procedure groups based on the time between Pending Dictation to Signed.

Parameters:

- Practice(s), Site(s), Date Range, Procedure Group(s)

**RTAT Detail Mammo**

**Practice(s):** Treesdale Radiology  
**Site(s):** All sites for the selected practice(s).  
**Date Range:** 09-01-2016 - 11-08-2016  
**Procedure Group(s):** (MA) Mammography  
**Excluding exams signed on weekends**

Radiologist: Pending Dictation - Pending Transcription (% of reports completed in < 3 hours)		
Week Ending	Tree	# Signed Exams
09-02-2016	0.00 %	1
<a href="#">1025483</a>		60.81
<b>Total</b>	<b>0.00 %</b>	<b>1</b>

Transcriptionist: Pending Transcription - Pending Signature (% of reports completed in < 3 hours)		
Week Ending	Tree	# Signed Exams
09-02-2016	100.00 %	1
<a href="#">1025483</a>		0.00
<b>Total</b>	<b>100.00 %</b>	<b>1</b>

Radiologist: Pending Signature - Signed (% of reports completed in < 2 hours)		
Week Ending	Tree	# Signed Exams
09-02-2016	100.00 %	1
<a href="#">1025483</a>		0.00
<b>Total</b>	<b>100.00 %</b>	<b>1</b>

RTAT Detail Mammo

Confidential and Proprietary  
 Page 1 of 1

Date Run: 11-08-2016 8:30 PM  
 RRIS-RSNAFormsGuestUser

### RTAT RAD Signing

Description: This report shows the report turnaround time from pending signature to signed by radiologist for selected year and practice.

Parameters:

- Practice, Year from and to Month

#### RTAT RAD Signing

Practice: Treedale Radiology

Date Range: 07-01-2016 - 09-30-2016

Description: This report shows the report turnaround time from pending signature to signed by radiologist for selected year and practice.

	September					August					July				
	<1 hr	1-2hr	RTAT < 2hrs	Monthly Total	% signed < 2hrs	<1 hr	1-2hr	RTAT < 2hrs	Monthly Total	% signed < 2hrs	<1 hr	1-2hr	RTAT < 2hrs	Monthly Total	% signed < 2hrs
Saltmarsh, Hilary	2		2	2	100.00 %	1		1	1	100.00 %	1		1	1	100.00 %
<b>Total</b>	<b>2</b>		<b>2</b>	<b>2</b>	<b>100.00 %</b>	<b>1</b>		<b>1</b>	<b>1</b>	<b>100.00 %</b>	<b>1</b>		<b>1</b>	<b>1</b>	<b>100.00 %</b>

RTAT RAD Signing

Confidential and Proprietary  
Page 1 of 1

Date Run: 11-08-2016 8:36 PM  
RRIS-RSNA\FORMS\GuestUser

### RTAT Summary

Description: This summary report shows the report turnaround time from pending signature to signed by week ending for selected date range and practice.

Parameters:

- Date Range and Practice(s)

### RTAT Summary

**Date Range:** 08-01-2016 - 11-08-2016

(Time exam ended to time report signed; % of reports finalized in < 8 hours)

Excluding exams signed on weekends.

*\*Note the grouping is by week ending, but you would need to specify From date to start on a monday and end date to be a friday when picking parameters to get the entire week(s).*

Non-Mammo Exams (Procedure Category is not BRE)		
Week Ending	Tree	# Signed Exams
08-19-2016	100.00 %	1
09-02-2016	0.00 %	1
09-09-2016	100.00 %	1
11-04-2016	100.00 %	3
<a href="#">1029890</a>		0.02
<a href="#">1028889</a>		0.02
<a href="#">1028891</a>		0.02
<b>Total</b>	<b>83.33 %</b>	<b>6</b>

Mammo Exams (Procedure Category is BRE)	
Week Ending	# Signed Exams
<b>Total</b>	<b>0</b>

RTAT Summary

Confidential and Proprietary  
Page 1 of 1

Date Run: 11-08-2016 8:38 PM  
RRIS-RSNA\FORMS\GuestUser

### Scheduled vs Performed

Description: This report shows exams that the procedure was changed from scheduled to performed for the dates and sites selected.

Parameters:

- Practice, Site(s) and Date Range

#### Scheduled vs Performed

Practice: Treesdale Radiology

Site(s): All sites for the selected practice(s).

Date Range: 01-01-2015 - 11-08-2016

Description: This report shows exams that the procedure was changed from scheduled to performed for the dates and sites selected.

Tech	DOS	Acc #	Tech modified to this procedure	Scheduled procedure	Status
Eckstein, Rhonda					1
	11-04-2016	<a href="#">1028889</a>	70460 - CT Head W	70460 - CT Head W	Signed1
<b>Grand Total</b>					<b>1</b>

Scheduled vs Performed



Confidential and Proprietary  
Page 1 of 1

Date Run: 11-08-2016 8:41 PM  
RRIS-RSNA\FORMS\GuestUser

### Special Accommodations

Description: This report displays studies where special accommodations have been identified as required for the selected practice/site(s), date and special accommodation filter.

Parameters:

- Practice, Site(s), Exam Date, ADA flag (American Disability Act flag from special accommodation lookup), and Special Accommodations

#### Special Accommodations

**Practice:** Treesdale Radiology

**Site(s):** All sites for the selected practice(s).

**Date:** 11-07-2016

**Special Accommodation(s):** Hearing Assistance, Interpreter Assistance, Lift Assistance, Vision Assistance, Wheel Chair Assistance

**Description:** This report displays studies where special accommodations have been identified as required for the selected practice/site(s), date and special accommodation filter.

Time Sched	Patient Name	Home Phone	Address	Exam(s)	ADA
<b>Cedarville</b>					
1:05 PM	Strickland, Rhonda	(256) 224-5160	39 Hague St. Santa Ana 64307	Acc#: 1030889 CT Abdomen W [74160] - Abdomen	Y
<b>ADA Flag:Y</b>	<b>Special Accomodation:</b> Hearing Assistance		<b>Other:</b>		
<b>Mapleton</b>					
10:50 AM	Freddieson, Fred	(902) 374-5967	907 East New Parkway Panama City 32406	Acc#: 1027871 CT Head Wo [70450] - Head	Y
<b>ADA Flag:Y</b>	<b>Special Accomodation:</b> Lift Assistance		<b>Other:</b> h		
<b>Pine Grove</b>					
9:00 AM	Hartman, Loren		325 Clarendon Boulevard Riverside 87864	Acc#: 1029889 CT Abdomen W [74160] - Abdomen	Y
<b>ADA Flag:Y</b>	<b>Special Accomodation:</b> Interpreter Assistance		<b>Other:</b> Interpreter Assistance		

Special Accommodations

Confidential and Proprietary  
Page 1 of 1

Date Run: 11-08-2016 8:43 PM  
RRIS-RSNAVFormsGuestUser



### Time to Schedule (seconds)

Description: This report displays the amount of time in seconds representing the start of a new appointment to the point of saving the booking.

Parameters:

- Practice, Site(s), and Date Range

### Time to Schedule (seconds)

**Practice:** Treesdale Radiology

**Site(s):** All sites for the selected practice(s).

**Date Range:** 10-01-2016 - 11-08-2016

**Description:** This report displays the amount of time in seconds representing the start of a new appointment to the point of saving the booking.

Created Date	Scheduled By	Procedures	Acc #	Time to Schedule
10-19-2016	Eckstein, Rhonda (eckstein)		<b>Average</b>	<b>125</b>
		CT Chest W [71260] - Chest	1026702	125
				<b>Average for Date</b>
10-28-2016	Saltmarsh, Hilary (saltmarsh)		<b>Average</b>	<b>175</b>
		CT Head W & Sinus W [CT22] - Head	1027403	294
	US Abdominal Complete [76700] - Abdomen	1027404	55	
			<b>Average for Date</b>	<b>175</b>
11-01-2016	Powell, Stick (stick)		<b>Average</b>	<b>16</b>
		PT/CT Sodium Fluoride PET Bone Scan [PT4]	1027888	17
	PT/CT Sodium Fluoride PET Bone Scan [PT4]	1027889	14	
			<b>Average for Date</b>	<b>16</b>
11-04-2016	Eckstein, Rhonda (eckstein)		<b>Average</b>	<b>137</b>
		CT Head W [70460] - Head	1028889	165
	CT Cervical Spine W [72126] - Cervical spine	1028890	109	
	Powell, Stick (stick)		<b>Average</b>	<b>70</b>
		CT Ear Orbit Sella Wo [70480] - Head	1029890	70
	Saltmarsh, Hilary (saltmarsh)		<b>Average</b>	<b>111</b>
XR Chest Min 4 Views [71030] - Chest	1028891	111		
			<b>Average for Date</b>	<b>114</b>
11-07-2016	Eckstein, Rhonda (eckstein)		<b>Average</b>	<b>278</b>
		CT Abdomen W [74160] - Abdomen	1030889	278
			<b>Average for Date</b>	<b>278</b>
11-08-2016	Eckstein, Rhonda (eckstein)		<b>Average</b>	<b>160</b>
		MR Abdomen W [74182] - Abdomen	1030894	160
		MR Abdomen Wo & Pelvis Wo [MR41] - Abdomen	1030893	160
		MR Brain W/Wo [70553] - Head	1030892	160
		MR Chest Wo [71550] - Chest	1030891	160
			<b>Average for Date</b>	<b>160</b>
			<b>Total Average</b>	<b>134</b>

Time to Schedule (seconds)

Confidential and Proprietary  
Page 1 of 1

Date Run: 11-08-2016 8:44 PM  
RRIS-RSNAFormsGuestUser

### Time to Schedule by User

Description: This report displays the amount of time in seconds representing the start of a new appointment to the point of saving the booking.

Parameters:

- Practice(s), Site(s), Date Range and Scheduled By (user)

### Time To Schedule By User

**Practice:** Treesdale Radiology

**Site(s):** All sites for the selected practice(s).

**Date:** 11-08-2016

**Description:** This report displays the amount of time in seconds representing the start of a new appointment to the point of saving the booking.

Created Date	Scheduled By	Procedures	Acc #	Time to Schedule
11-08-2016	Eckstein, Rhonda (eckstein)		<b>Average</b>	<b>160</b>
		MR Abdomen W [74182] - Abdomen	1030894	160
		MR Abdomen Wo & Pelvis Wo [MR41] - Abdomen	1030893	160
		MR Brain W/Wo [70553] - Head	1030892	160
		MR Chest Wo [71550] - Chest	1030891	160
			<b>Average for Date</b>	<b>160</b>
			<b>Total Average</b>	<b>160</b>

Time To Schedule By User

Confidential and Proprietary  
Page 1 of 1

Date Run: 11-08-2016 8:46 PM  
RRIS-RSNA\FORMS\GuestUser

### Transportation Report

Description: This report displays transportation requirements for the practice, site(s) and date range selected (based on DOS).

Parameters:

- Practice, Site(s) and Date

#### Transportation Report

Practice: Treesdale Radiology

Site(s): All sites for the selected practice(s).

Date Range: 01-01-2016 - 11-08-2016

Description: This report displays transportation requirements for the practice, site(s) and date range selected (based on DOS).

DOS	Time	Patient Name	Home Phone	Address	Exam(s)	Referring	
<b>Cedarville</b>							<b>1</b>
11-07-2016							<b>1</b>
	1:05 PM	Strickland, Rhonda	(256) 224-5160	39 Hague St. Santa Ana 64307	CT Abdomen W [74160] - Abdomen	Arnold, Bob	
<b>Transportation Notes:</b> Taxi Needed							
<b>Mapleton</b>							<b>1</b>
11-07-2016							<b>1</b>
	10:50 AM	Freddieson, Fred	(902) 374-5967	907 East New Parkway Panama City 32406	CT Head Wo [70450] - Head	Flowers, Tim	
<b>Transportation Notes:</b>							
<b>Total Count</b>							<b>2</b>

Transportation Report

Confidential and Proprietary  
Page 1 of 1

Date Run: 11-08-2016 8:47 PM  
RRIS-RSNA\FORMS\GuestUser

### Unfinished Exams

Description: This report lists the exams currently in the statuses selected between the dates specified for the selected practice(s), site(s).

Parameters:

- Practice(s), Site(s), Date Range and Status(s)

### Unfinished Exams

**Practice(s):** Treedale Radiology

**Site(s):** All sites for the selected practice(s).

**Date:** 11-01-2016

**Status(s):** Scheduled, CheckingIn, Arrived, Started

**Description:** This report lists the exams currently in the statuses selected between the dates specified for the selected practice(s), site(s).

Date	Patient	MRN	Accession	Procedure Code	Procedure Description	Appt Time
<b>Cedarville</b>						
<b>11-01-2016</b>						
	Cam 1 Cammerson	y3	<a href="#">1027867</a>	72141	MR Cervical Spine Wo [72141] - Cervical spine	10:30 AM
	<b>Status:</b> Scheduled		<b>Modality:</b> MR 1 at Cedar			
			<b>User:</b> -			
	Ed 1 Effortson	y5	<a href="#">1027869</a>	70551	MR Brain Wo [70551] - Head	11:00 AM
	<b>Status:</b> Scheduled		<b>Modality:</b> MR 1 at Cedar			
			<b>User:</b> -			
	Lou 1 Lammerson	y12	<a href="#">1027878</a>	73080	XR Elbow Min 3 Views [73080] - Elbow - Left	11:35 AM
	<b>Status:</b> Arrived		<b>Modality:</b> XR 1 at Cedar			
			<b>User:</b> -			
	Adam 1 Adamson	y20	<a href="#">1027879</a>	70450	CT Head Wo [70450] - Head	1:15 PM
	<b>Status:</b> Arrived		<b>Modality:</b> CT 1 at Cedar			
			<b>User:</b> stick - Powell, Stick			
	Bart 1 Bartson	y21	<a href="#">1027880</a>	73080	XR Elbow Min 3 Views [73080] - Elbow - Left	1:20 PM
	<b>Status:</b> Scheduled		<b>Modality:</b> XR 1 at Cedar			
			<b>User:</b> -			
	Gran 1 Grannyson	y26	<a href="#">1027885</a>	73080	XR Elbow Min 3 Views [73080] - Elbow - Right	1:25 PM
	<b>Status:</b> Scheduled		<b>Modality:</b> XR 1 at Cedar			
			<b>User:</b> -			
	Dent 1 Dentson	y23	<a href="#">1027882</a>	72148	MR Lumbar Spine Wo [72148] - Lumbar spine	2:00 PM
	<b>Status:</b> Scheduled		<b>Modality:</b> MR 1 at Cedar			
			<b>User:</b> -			

### Unknown Referring

Description: This report displays exams where requested by physician is “Unknown referring” for a practice, site(s) and selected date range (DOS) excluding outside reads. Unlike the unknown referring physician, this report looks for a referring with a last\_name like Refer and first name like Unknown.

Parameters:

- Practice, Site(s), Date Range and Statuses to ignore

#### UNKNOWN REFERRING

Practice: Treesdale

Site(s): Cedarville

Date: 11-08-2016

Description: This report displays exams where requested by physician is like “Unknown” and last name is like “Refer” for a practice and site(s). Both active and inactive physicians are included in this report.

\*Excluding status(s): Cancelled,Discontinued,Aborted,Scheduled,Void

Patient Name	MRN	Procedures	Referring	Accession	Status
<b>CEDARVILLE</b>				<b>Site Total :</b>	<b>2</b>
<b>11-08-2016</b>				<b>Date Total :</b>	<b>2</b>
<b>JOHNSON . JANE</b>	7010	US Renal and Bladder [76770A] - Abdomen	REFER, MD, UNKNOWN	7077	ExamDone
<b>PLAT, SUE</b>	37707	MA Digital Screening Mammo W CAD [G0202] - Breast - B	REFER, MD, UNKNOWN	1867	ExamDone
<b>Grand Total</b>					<b>2</b>

Unknown Referring

Confidential and Proprietary  
Page 1 of 1

Date Run: 11-08-2016 10:04 PM  
RADNET\ecsteinR

**Unread Exams**

Description: This report displays exams in an unread status based on selected practice(s) and site(s).

Parameters:

- Practice(s), Site(s) and Exam Types (modality type)

**Unread Exams**

**Practice(s):** Treesdale Radiology

**Site(s):** All sites for the selected practice(s).

**Date Range:** 11-01-2016 - 11-08-2016

**Exam Types:** Angiography,Bone Density,Special Charges,CT,Endoscopy,Fluoroscopy,Mammography,MRI,Nuclear Medicine,PET,Radiation Therapy,Secondary Capture,Special Procedures,Ultrasound,Xray

**Description:** This report displays exams in an unread status based on selected practice(s) and site(s) and date range selected (DOS).

Performed	Patient	MRN	Acc #	Procedures	Performed By	Modality Type	Site
11-01-2016 10:44 AM	Zulu, Average	z12345	1027790	CT Ankle Arthrogram [CT31] - Ankle joint - Left	Waite, Andrew R	CT	Mapleton
11-01-2016 10:44 AM	Adams, Average	z10000001	1027768	CT Ankle Arthrogram [CT31] - Ankle joint - Left	Waite, Andrew R	CT	Mapleton
11-01-2016 10:44 AM	O'Brian, Hillary	z20000002	1027778	CT Ankle Arthrogram [CT31] - Ankle joint - Left	Waite, Andrew R	CT	Mapleton
11-01-2016 10:44 AM	O'Conner, Average	z30000003	1027788	CT Ankle Arthrogram [CT31] - Ankle joint - Left	Waite, Andrew R	CT	Mapleton
11-01-2016 10:44 AM	O'Dell, Jayne	z40000004	1027798	CT Ankle Arthrogram [CT31] - Ankle joint - Left	Waite, Andrew R	CT	Mapleton
11-01-2016 10:44 AM	O'Neil, Average	z50000005	1027808	CT Ankle Arthrogram [CT31] - Ankle joint - Left	Waite, Andrew R	CT	Mapleton
11-01-2016 10:44 AM	Jones, Sue	z50000006	1027818	CT Ankle Arthrogram [CT31] - Ankle joint - Left	Waite, Andrew R	CT	Mapleton
11-07-2016 10:39 AM	Strickland, Rhonda	9695	1030889	CT Abdomen W [74160] - Abdomen	Eckstein, Rhonda	CT	Cedarville
							<b>Count: 8</b>

Unread Exams

Confidential and Proprietary  
Page 1 of 1

Date Run: 11-08-2016 8:53 PM  
RRIS-RSNA\FORMS\GuestUser

### Wait Room Times

Description: This report displays the waiting room times. Calculation is based on if arrival date is greater than scheduled start date the calculation is based on arrived to started. Otherwise it is based on scheduled date to actual start date.

Parameters:

- Practice, Site(s), Date Range and Procedure group(s)

### Wait Room Times

**Practice:** Treesdale Radiology

**Site(s):** All sites for the selected practice(s).

**Date Range:** 11-01-2016 - 11-08-2016

**Description:** This report displays the waiting room times. Calculation is based on if arrival date is greater than scheduled start date the calculation is based on arrived to started. Otherwise it is based on scheduled date to actual start date.

Site	Procedure Type	Accession #	WRT <= 10 minutes	WRT > 10 minutes	AVG WRT
Cedarville			100.00 %	0.00 %	1.25
	☒ CT		100.00 %	0.00 %	1.33
		<a href="#">1028889</a>		2	
		<a href="#">1030889</a>		2	
		<a href="#">1029890</a>			
	☒ XR		100.00 %	0.00 %	1.00
		<a href="#">1028891</a>		1	
Mapleton			100.00 %	0.00 %	0.00
	☒ CT		100.00 %	0.00 %	0.00
<b>Total</b>			<b>100.00 %</b>	<b>0.00 %</b>	<b>0.45</b>

Wait Room Times

Confidential and Proprietary  
Page 1 of 1

Date Run: 11-08-2016 8:58 PM  
RRIS-RSNA\FORMS\GuestUser

### Wait Room Turnover

Description: This report displays the wait time by interval counts and percentages are grouped by procedure group.

Parameters:

- Practice, Site(s), Date Range and Procedure group(s)

#### Wait Room Turnover

Practice: Treesdale Radiology

Site(s): All sites for the selected practice(s).

Date Range: 11-01-2016 - 11-08-2016

Description: This report displays the wait time by interval counts and percentages are grouped by procedure group.

Site	Interval	CT	XR	Total
<b>Cedarville</b>		<b>3 27.27 %</b>	<b>1 9.09 %</b>	<b>4</b>
	0-5	3 100.00 %	1 100.00 %	4 100.00 %
	5-10	0.00 %	0.00 %	0.00 %
	10-15	0.00 %	0.00 %	0.00 %
	15-20	0.00 %	0.00 %	0.00 %
	20-30	0.00 %	0.00 %	0.00 %
	30-40	0.00 %	0.00 %	0.00 %
	>40	0.00 %	0.00 %	0.00 %
<b>Mapleton</b>		<b>7 63.64 %</b>	<b>0 0.00 %</b>	<b>7</b>
	0-5	7 100.00 %	0.00 %	7 100.00 %
	5-10	0.00 %	0.00 %	0.00 %
	10-15	0.00 %	0.00 %	0.00 %
	15-20	0.00 %	0.00 %	0.00 %
	20-30	0.00 %	0.00 %	0.00 %
	30-40	0.00 %	0.00 %	0.00 %
	>40	0.00 %	0.00 %	0.00 %
<b>Grand Total</b>		<b>10 90.91 %</b>	<b>1 9.09 %</b>	<b>11</b>

Wait Room Turnover

Confidential and Proprietary  
Page 1 of 1

Date Run: 11-08-2016 8:59 PM  
RRIS-RSNA\FORMS\GuestUser



## Management Reports \ Peer Review

### Peer Review Report

Description: There are two reports within this report. The first is an overall on each practice selected, while the second set of report is on each peer review radiologist within the practices and reporting period selected. Each practice and radiologist will be on their own page when the report is run.


Parameters: Date Range, and Practice

#### Peer Review Report - Practice Overall Report

Practice: **Treesdale**  
 Date Range: 01-01-2016 - 11-08-2016  
 # Practices in report: 1

Description: This report shows all the peer review performed and respective scoring categories within the selected practice(s) and date range.

	Total		1		2		3		4		Clinically Significant	
<b>Nuclear Medicine</b>	1	0.09%	1	100.00%	0	0.00%	0	0.00%	0	0.00%	0	0.00%
<b>Dexa</b>	2	0.19%	2	100.00%	0	0.00%	0	0.00%	0	0.00%	0	0.00%
<b>PET/CT</b>	8	0.75%	5	62.50%	3	37.50%	0	0.00%	0	0.00%	0	0.00%
<b>CT</b>	78	7.36%	64	82.05%	7	8.97%	6	7.69%	1	1.28%	0	0.00%
<b>MRI</b>	129	12.17%	111	86.05%	9	6.98%	9	6.98%	0	0.00%	3	2.33%
<b>CR XR</b>	149	14.06%	113	75.84%	19	12.75%	13	8.72%	4	2.68%	13	8.72%
<b>Ultrasound</b>	161	15.19%	156	96.89%	1	0.62%	3	1.86%	1	0.62%	0	0.00%
<b>DR XR</b>	204	19.25%	181	88.73%	9	4.41%	12	5.88%	2	0.98%	6	2.94%
<b>Digital Mammography</b>	328	30.94%	303	92.38%	18	5.49%	4	1.22%	3	0.91%	3	0.91%
<b>Total</b>	1060	100.00%	936	88.30%	66	6.23%	47	4.43%	11	1.04%	25	2.36%

	<b>All included Practices</b>	36.52%	2.58%	1.83%	0.43%	0.98%
	Peer Review Report	Confidential and Proprietary		Date Run: 11-08-2016 9:56 PM		
		Page 1 of 74		RADNET/E: ksteinR		

**Review Summary by Radiologist**

Description: This report shows the number of peer reviews performed by a radiologist by modality type within the practice and date range selected.

Parameters:

- Date Range, Practice

**Review Summary by Radiologist**

Practice *Treesdale*

Date Range: 09-01-2016 - 11-08-2016

# Radiologists in this report: 10

Description: This report shows the number of peer reviews performed by a radiologist by modality type within the practice and date range selected.

Radiologist	CT	DMA	DR	MR	US	Total
1000033435					2	2
1000033439	1			2		3
1000033440					1	1
1000033446		2				2
1000033447		2				2
1000033449				2		2
1000033450					2	2
1000033451				1		1
1000036006		1				1
1000041004			2			2
<b>Total</b>	<b>1</b>	<b>5</b>	<b>2</b>	<b>5</b>	<b>5</b>	<b>18</b>



Confidential and Proprietary  
Page 1 of 1

Date Run: 11-08-2016 10:00 PM  
RADNET\EcksteinR

## Management Reports \ Sales

### Referrals by Physician

Description: This report displays procedure counts based on scheduled exams falling between the dates specified for the chosen practice, sites and modality types. Specifically you can choose the Top N physicians.

*\*\*excluding statuses cancelled, discontinued and voided exams.*

Parameters:

- Practice(s), Site(s), Date Range, and Modality type

### Referrals by Physician

**Practice(s):** Treesdale Radiology


**Site(s):** All sites for the selected practice(s).

**Date Range:** 10-01-2016 - 11-08-2016

**Top N:** 30

**Description:** This report displays procedure counts based on scheduled exams falling between the dates specified for the chosen practice, sites and modality types. Specifically you can choose the Top N physicians.

*\*\*excluding statuses cancelled, discontinued and voided exams.*

Physician	+ -	Procedure Type	Counts
<input type="checkbox"/> Arnold, Bob		Total	8
		CT	7
		XR	1
<input type="checkbox"/> Waite, Andrew		Total	8
<input type="checkbox"/> Edwards, John		Total	5
<input type="checkbox"/> Graves, Henry		Total	3
<input type="checkbox"/> Bihand, Bob		Total	1
<input type="checkbox"/> Brown, Bob		Total	1
<input type="checkbox"/> Dolengo, Tim		Total	1
<input type="checkbox"/> Flowers, Tim		Total	1
<input type="checkbox"/> Armiger, Tom		Total	1
<input type="checkbox"/> Howard, Tom		Total	1
<input type="checkbox"/> Cook, Darryl		Total	1
<input type="checkbox"/> Krinsky, Darryl		Total	1
<input type="checkbox"/> Cho, Clifton		Total	1
<input type="checkbox"/> Kim, Clifton		Total	1
<input type="checkbox"/> Joseph, Henry		Total	1
<input type="checkbox"/> Harding, Luca		Total	1
<input type="checkbox"/> Washington, George		Total	1

### Referrals by Physician Patient List

Description: This report displays procedure counts based on scheduled exams falling between the dates specified for the chosen practice, sites and modality types. It is broken down by distinct patients with the total exams scheduled.

*\*\*excluding statuses cancelled, discontinued and voided exams.*

Parameters:

- Practice(s), Site(s), Date Range and Modality Types

### Referrals by Physician Patient List

**Practice(s):** Treesdale Radiology

**Site(s):** All sites for the selected practice(s).

**Date Range:** 11-01-2016 - 11-08-2016

**Description:** This report displays procedure counts based on scheduled exams falling between the dates specified for the chosen practice, sites and modality types. It is broken down by distinct patients with the total exams scheduled.

*\*\*excluding statuses cancelled, discontinued and voided exams.*

Referring Physician	MRN (distinct)	# Exams	Patient Name	Procedure Type	Procedures	Modality
<input type="checkbox"/> Armiger, Tom G		1	1			
<input type="checkbox"/> Arnold, Bob		5	7			
	1378		Hartman, Loren	CT	CT Abdomen W [74160] - Abdomen	CT 1 at Pine
	8211		Parsons, Hilary	XR	XR Chest Min 4 Views [71030] - Chest	XR 1 at Cedar
	9695		Strickland, Rhonda	CT	CT Head W [70460] - Head	CT 1 at Cedar
	4454		Rodríguez, Rhonda	CT	CT Cervical Spine W [72126] - Cervical spine	CT 1 at Cedar
	z40000004		O'Dell, Jayne	CT	CT Ear Orbit Sella Wo [70480] - Head	CT 1 at Cedar
	1378		Hartman, Loren	CT	CT Abdomen W [74160] - Abdomen	CT 1 at Pine
	9695		Strickland, Rhonda	CT	CT Abdomen W [74160] - Abdomen	CT 1 at Cedar
<input type="checkbox"/> Bihand, Bob		1	1			
<input type="checkbox"/> Cho, Clifton J		1	1			
<input type="checkbox"/> Cook, Darryl J		1	1			
<input type="checkbox"/> Dolengo, Tim		1	1			
<input type="checkbox"/> Edwards, John E		5	5			
<input type="checkbox"/> Flowers, Tim		1	1			
<input type="checkbox"/> Graves, Henry		3	3			
<input type="checkbox"/> Harding, Luca		1	1			
<input type="checkbox"/> Howard, Tom H		1	1			
<input type="checkbox"/> Joseph, Henry		1	1			
<input type="checkbox"/> Kim, Clifton		1	1			
<input type="checkbox"/> Krinsky, Darryl		1	1			
<input type="checkbox"/> Waite, Andrew R		8	8			
<b>Total Patients:</b>		<b>29</b>	<b>34</b>			

### Referring Physician Volumes

Description: This report displays the orders placed by the referring's select for the practice, sites, date range, referring practices and specialties selected.

\*Date chosen will assume the end of the selected month was intended.

Parameters:

- Practice, Site(s), Date Range, Referring Practice(s), Specialties and Referring(s)

### Referring Physician Volumes

**Practice:** Treesdale Radiology

**Site(s):** All sites for the selected practice(s).

**Date Range:** 11-01-2016 - 11-08-2016

**Description:** This report displays the orders placed by the referrals select for the practice,sites,date range, referring practices and specialties selected.

\*Date chosen will assume the end of the selected month was intended.

Referring Physician	Order Status	Volume
<input type="checkbox"/> Arnold, Bob		8
	OrderSigned	8
<input type="checkbox"/> Edwards, John		3
	OrderSigned	3
<input type="checkbox"/> Graves, Henry		3
<input type="checkbox"/> Jackson, Bob		1
<input type="checkbox"/> Bihand, Bob		1
<input type="checkbox"/> Dolengo, Tim		1
<input type="checkbox"/> Flowers, Tim		1
<input type="checkbox"/> Armiger, Tom		1
<input type="checkbox"/> Howard, Tom		1
<input type="checkbox"/> Cook, Darryl		1
<input type="checkbox"/> Krinsky, Darryl		1
<input type="checkbox"/> Cho, Clifton		1
<input type="checkbox"/> Kim, Clifton		1
<input type="checkbox"/> Waite, Andrew		1
<input type="checkbox"/> Joseph, Henry		1
<b>Total</b>		<b>26</b>

Referring Physician Volumes

Confidential and Proprietary  
Page 1 of 1

Date Run: 11-08-2016 9:08 PM  
RRIS-RSNA\FORMS\GuestUser

## Management Reports \ Statistics

### Management Reports \ Statistics \ Editors

#### Edited Jobs Activity Report

Description: This report displays the number of transcribed reports by editor broken down by procedure groups for the selected practice, site and date range.

Parameters:

- Practice, Site(s), Date Range and Editor(s)

#### Edited Jobs Activity Report

**Practice:** Treesdale Radiology

**Site(s):** Cedarville, Mapleton, Pine Grove

**Date Range:** 11-01-2012 - 11-10-2016

**Description:** This report displays the number of transcribed reports by editor broken down by procedure groups for the selected practice, site and date range.

Editor	CT	UG	Total Edits	% of Total
<b>Treesdale Radiology</b>	4	1	5	100.00 %
<i>Cedarville</i>		1	1	20.00 %
Test User, System Admin		1	1	20.00 %
<i>Mapleton</i>	2		2	40.00 %
Noye, Darcy Thane	1		1	20.00 %
Test User, System Admin	1		1	20.00 %
<i>Pine Grove</i>	2		2	40.00 %
Waite, Andrew R	2		2	40.00 %
<b>Grand Total</b>	<b>4</b>	<b>1</b>	<b>5</b>	

Edited Jobs Activity Report

Confidential and Proprietary  
Page 1 of 1

Date Run: 11-10-2016 9:06 PM  
RRIS-RSNA\FORMS\GuestUser

### Edited Jobs By User Detail

Description: This report details the editor activity for editors within the specified practice, sites and dates provided.

Parameters:

- Practice, Site(s), Date Range, and Editor(s)

#### Edited Jobs By User Detail

Practice: Treesdale Radiology

Site(s): All sites for the selected practice(s).

Date Range: 10-30-2012 - 11-10-2016

Description: This report details the editor activity for editors within the specified practice, sites and dates provided.

Edited By	Accession	Edited Date	Status	Dictated By	Dictation Time	Dictated Date	Turnaround
<input type="checkbox"/> Waite, Andrew <span style="float: right;">2</span>							
	<a href="#">1011708</a>	11-24-2012 1:17 PM	Signed1	Waite, Andrew		11-24-2012 1:14 PM	00:03:00
	<a href="#">1011711</a>	11-24-2012 1:17 PM	Signed1	Waite, Andrew		11-24-2012 1:14 PM	00:03:00
<b>Total</b> <span style="float: right;"><b>2</b></span>							

Edited Jobs By User Detail

Confidential and Proprietary  
Page 1 of 1

Date Run: 11-10-2016 11:51 AM  
RRIS-RSNA\FORMS\GuestUser

**Editors Activity**

Description: This report displays the counts by editor for each site and procedure group code for the selected Practice, Site(s) Date Range and procedure group(s).

Parameters: Practice, Site(s), Date Range, Editor(s) and Procedure Group(s)

**Editors Activity**

**Practice:** Treesdale Radiology

**Site(s):** All sites for the selected practice(s).

**Date Range:** 10-30-2012 - 11-10-2016

**Procedure Groups:** DE, CT, RD, PE, EM, FL, RS, SR, MA, MR, NM, OT, PT, XR, RO, UG, UV

**Description:** This report displays the counts by editor for each site and procedure group code for the selected Practice, Site(s) date range and procedure group(s).

	CT	Total
<b>Mapleton</b>	<b>1</b>	<b>1</b>
darcyn	1	1
<b>Pine Grove</b>	<b>3</b>	<b>3</b>
Andrew	2	2
denis	1	1
<b>Total</b>	<b>4</b>	<b>4</b>

Editors Activity

Confidential and Proprietary  
Page 1 of 1

Date Run: 11-10-2016 11:48 AM  
RRIS-RSNAFormsGuestUser



## Management Reports \ Statistics \ IVT

### IVT User Activity

Description: This report displays the counts # of patients IVT activity has been logged by users within the selected sites/practice. Report is tracking any edits user(s) have made to study verification and pre-certification notes.

\*Outside reads are excluded from this report.

Parameters:

- Practice, Site(s), Date Range, User(s) and Room(s)

### IVT User Activity

**Practice:** Treesdale Radiology

**Site(s):** All sites for the selected practice(s).

**Date Range:** 01-01-2014 - 11-08-2016

**Description:** This report displays the counts # of patients IVT activity has been logged by users within the selected sites/practice. Report is tracking any edits user(s) have made to study verification and pre-certification notes.

\*Outside reads are excluded from this report.

Date / Practice / Site - User	MRN	Patient Name	Accession #
<b>11-07-2016</b>	<b>1</b>		
<b>Treesdale Radiology</b>		<b>1</b>	
<input type="checkbox"/> Cedarville		<b>1</b>	
<input type="checkbox"/> eckstein		<b>1</b>	
Notes made on:	9695	Strickland Rhonda	1030889
<b>Total Patients:</b>		<b>1</b>	

IVT User Activity

Confidential and Proprietary  
Page 1 of 1

Date Run: 11-08-2016 9:13 PM  
RRIS-RSNA\FORMS\GuestUser

## Management Reports \ Statistics \ Radiologists

### Dictations by User

Description: This report displays the dictations by user id and displays the total amount of their reads and the percentage of self edits. \*Linked counts represent the additional number of studies associated with the reports.

Parameters:

- Practice, Site(s), Date Range and Radiologists

### Dictations by User

**Practice:** Advanced Radiology

**Site(s):** All sites for the selected practice(s).

**Date Range:** 01-01-2015 - 07-24-2015

**Description:** This report displays the dictations by user id and displays the total amount of their reads and the percentage of self edits. \*Linked counts represents the additional number of studies associated with the reports.

Dictated By User	Dictated Date	# Reads	# Self Edits	% Self Edits	%	Linked Counts
Yeo Jr, Dale E (dale)		7	5	71.43%	100.00 %	2
	☒ 01-21-2015	1	0	0.00%		
	☒ 03-27-2015	1	1	100.00%		
	☒ 05-11-2015	5	4	80.00%		2
<b>Total</b>		<b>7</b>	<b>5</b>			<b>2</b>

### QA Radiologist Review

Description: This report displays the Rejected QA reports and history associated with those studies for exams reviewed within the dates specified.

Parameters:

- Practice(s), Site(s) and Date Range

### QA Radiologist Review

Practice(s): Treesdale

Site(s): All sites for the selected practice(s).

Date Range: 08-01-2014 - 11-10-2016

Description: This report displays the Rejected QA reports and history associated with those studies for exams reviewed within the dates specified.

Dictated By	Modality Type	Accession	Status	Quality Code	Reviewed by	Reviewed Date
						217
	<input checked="" type="checkbox"/> CT					22
	<input checked="" type="checkbox"/> DMA					75
	<input checked="" type="checkbox"/> DR					55
	<input checked="" type="checkbox"/> FL					1
	<input checked="" type="checkbox"/> MR					8
	<input type="checkbox"/> PCT					3
		<b>123456</b>	<b>PETCT Tumor Imaging Skull-Thigh [78815] - Body</b>			
			Signed1	REJECTED		09-02-2014 2:13 PM
		<b>234565</b>	<b>PETCT Tumor Imaging Skull-Thigh [78815] - Body</b>			
			Signed1	REJECTED		09-19-2014 3:14 PM
		<b>432543</b>	<b>PETCT Tumor Imaging Skull-Thigh [78815] - Body</b>			
			Signed1	REJECTED		09-19-2014 5:56 PM
	<input checked="" type="checkbox"/> US					53

### Radiologist Dictation Activity by Date

Description: This report details the number of dictations by radiologist, by dictation date, the % of dictations done and the breakdown by procedure category. \*Linked counts represent the additional number of studies associated with the reports.

Parameters:

- Practice, Site(s), Date Range, Procedure Group(s) and Radiologist(s)

### Radiologist Dictation Activity by Date

**Practice:** Treesdale Radiology

**Site(s):** All sites for the selected practice(s).

**Date Range:** 01-01-2016 - 11-08-2016

**Procedure Group(s):** All procedure groups.

**Description:** This report details the number of dictations by radiologist, by dictation date, the % of dictations done and the breakdown by procedure category. \*Linked counts represents the additional number of studies associated with the reports.

Radiologist	Date	# of Dictated	% of Total	CT	MA	MR	UG	XR	Linked Counts
Eckstein, Rhonda		1	4.76%	1					
	11-04-2016	1	4.76%	1					
Powell, Dr Stick		1	4.76%	1					
	11-04-2016	1	4.76%	1					
Saltmarsh (Rad), Hilary		4	19.05%	2	2				1
	03-16-2016	2	9.52%	2					
	05-17-2016	1	4.76%		1				
	05-21-2016	1	4.76%		1				1
Harding, Clifton		15	71.43%	6	4	2	2	1	
	02-24-2016	1	4.76%		1				
	03-11-2016	1	4.76%		1				
	03-16-2016	1	4.76%	1					
	04-28-2016	2	9.52%				2		
	06-28-2016	1	4.76%	1					
	06-30-2016	3	14.29%	3					
	07-08-2016	1	4.76%	1					
	08-17-2016	1	4.76%		1				
	09-02-2016	1	4.76%		1				
	09-07-2016	1	4.76%			1			
	10-11-2016	1	4.76%			1			
	11-04-2016	1	4.76%					1	
	<b>Total</b>		<b>21</b>	<b>100.00%</b>	<b>10</b>	<b>6</b>	<b>2</b>	<b>2</b>	<b>1</b>

### Radiologist Productivity by User - Mammo

Description: This report displays the selected mammography exams read by Radiologist for a specified period. All Exams are Signed exams.\*Linked counts represent the additional number of studies associated with the reports.

Parameters:

- Practice, Site(s), Date Range, Mammo Procedure(s) and Radiologist(s)

### Radiologist Productivity by User - Mammo

**Practice:** Advanced Radiology

**Site(s):** All sites for selected practice

**Date Range:** [fromDate][dash]«Exp»

**Procedure Codes:** All mammography procedures

**Description:** This report displays the selected mammography exams read by Radiologist for a specified period. All Exams are Signed exams.\*Linked counts represents the additional number of studies associated with the reports.

Procedure	Count	%	Linked Counts
<b>Yeo, Dale (dale)</b>	1	100.00%	
MA DIGITAL MAMMO SCREEN BIL	1	100.00%	
<b>Total</b>	1	100.00%	

### Radiologist Signing Activity

Description: This report displays the signing counts by radiologist broken down by procedure group for the Practice\Site(s) and date range. \*Linked counts represent the additional number of studies associated with the reports.

Parameters:

- Practice, Site(s), Date Range, Procedure Group(s) and Radiologist(s)

### Radiologist Signing Activity

**Practice:** Advanced Radiology

**Site(s):** Eldersburg, Fisher, Lutherville, PowerVille, OutsideMed (outside read), Test Outside Read (outside read)

**Date Range:** 01-01-2015 - 07-24-2015

**Procedure Group(s):** All procedure groups.

**Description:** This report displays the signing counts by radiologist broken down by procedure group for the Practice\Site(s) and date range. \*Linked counts represents the additional number of studies associated with the reports.

Radiologist	# of Signed	% of Total	CT	MA	Linked Counts
Yeo Jr, Jr Dale E	5	100.00%	4	1	2
<b>Total</b>	<b>5</b>	<b>100.00%</b>	<b>4</b>	<b>1</b>	<b>2</b>

## Management Reports \ Statistics \ Registration

### Front Desk Detailed Activity by Hour

Description: This report displays a chart of the counts of arrived patients broken down by hour of the day. Also included is a grid of the details.

Parameters:

- Practice, Site(s), Date and User(s)


### Front Desk Detailed Activity by Hour

Practice: Treesdale

Site(s): All sites for the selected practice(s).

Date: 11-07-2016

Description: This report displays a chart of the counts of arrived patients broken down by hour of the day. Also included is a grid of the details.



Arrived by user id	Arrival Time	Accession #	9:00	10:00	11:00	12:00	13:00	14:00	15:00	16:00	17:00	18:00	19:00	20:00	Total
<input type="checkbox"/> eckstein	Total		1	2	4	1	1	1		4					14
<input type="checkbox"/> saltmarsh	Total							4	4	5	1	9	4	3	30
<input type="checkbox"/> powell	Total		1	3	12	3	10	6	5	1	8	3			52
<b>Total</b>			<b>2</b>	<b>5</b>	<b>16</b>	<b>4</b>	<b>11</b>	<b>11</b>	<b>9</b>	<b>10</b>	<b>9</b>	<b>12</b>	<b>4</b>	<b>3</b>	<b>96</b>

Front Desk Detailed Activity by Hour

Confidential and Proprietary  
Page 1 of 1

Date Run: 11-08-2016 9:45 PM  
RADNET\EcksteinR

### Front Desk Detailed Activity

Description: This report displays the detailed information of checkins, scheduled events, # walkins, # verifications for the selected user within the selected sites/practice.

\*Exams that are created and scheduled for the same day are considered walkins.

\*Outside reads are excluded from this report

Parameters:

- Practice, Site(s), Date Range and User ID

### Front Desk Detailed Activity

**Practice:** Treesdale Radiology

**Site(s):** All sites for the selected practice(s).

**Date Range:** 10-01-2016 - 11-08-2016

**Description:** This report displays the detailed information of checkins, scheduled events, # walkins, # verifications for the selected user within the selected sites/practice.

\*Exams that are created and scheduled for the same day are considered walkins.

\*Outside reads are excluded from this report

	Check in	Scheduling Events			Verified
		Scheduled	Walk in	Total	
<b>Eckstein, Rhonda (eckstein)</b>	2	4	4	8	0
☐ 11-08-2016	0	4	0	4	0
1030894		1		1	
1030893		1		1	
1030892		1		1	
1030891		1		1	
☐ 11-07-2016	1	0	1	1	0
1030889	1		1	1	
☐ 11-04-2016	1	0	2	2	0
1028890			1	1	
1028889	1		1	1	
☐ 10-19-2016	0	0	1	1	0
1026702			1	1	
<b>Total</b>	<b>2</b>	<b>4</b>	<b>4</b>	<b>8</b>	<b>0</b>



## Management Reports \ Statistics \ Reviewers

### QA Review Statistics

Description: This report displays the statistics by reviewer for the practice(s) and date range specified.

Parameters: Practice(s), Site(s) and Date Range

### QA Review Statistics

Practice(s): Treesdale Radiology

Date Range: 01-01-2011 - 11-10-2016

Description: This report displays the statistics by reviewer for the practice(s) and date range specified.

Reviewer	Review Type	Accession	Reviewed Date	
<b>Treesdale Radiology</b>				<b>2</b>
<input type="checkbox"/> Eckstein, Rhonda				2
	<a href="#">TranscribedReportQA</a>			2
		<a href="#">1010777</a>	11/10/2016 12:02:06 PM	
		<a href="#">1010960</a>	11/10/2016 12:01:19 PM	
<b>Total</b>				<b>2</b>

QA Review Statistics

Confidential and Proprietary  
Page 1 of 1

Date Run: 11-10-2016 12:02 PM  
RRIS-RSNAVFormsGuestUser



## Management Reports \ Statistics \ Schedulers

### Scheduler Activity

Description: This report displays the number of exams a scheduler created, or cancelled for selected period. Broken down to the hour of day. As well as displays the number of exams a scheduler created, or cancelled per day.

Parameters:

- Practice, Site(s), Date Range, User(s), Room(s), Show scheduled(Y/N) and Show Cancelled(Y/N)

### Scheduler Activity

**Practice:** Treesdale Radiology

**Site(s):** Cedarville, Mapleton, Pine Grove

**Date Range:** 11-01-2016 - 11-08-2016

**Room(s):**

CT1CE,CT1MP,CT1PG,DE1MP,DE1PG,MA1CE,MA1MP,MA1PG,MA2CE,MA2PG,MA3PG,MR1CE,MR1MP,MR1PG,MR2CE,MR2PG,MR3CE,NM1PG,NM2PG,PT1MP,PT1PG,US1CE,US1MP,US1PG,US2CE,US2MP,US2PG,US3PG,XR1CE,XR1MP,XR1PG,XR2PG

**Description:** This report displays the number of exams a scheduler created, or cancelled for selected period. Broken down to the hour of day. As well as displays the number of exams a scheduler created, or cancelled per day.

Scheduled by		9:00 AM	10:00 AM	11:00 AM	12:00 PM	2:00 PM	4:00 PM	Total
Harding, Clifton	<input checked="" type="checkbox"/> Scheduled							0
	<input type="checkbox"/> Cancelled		1	1				2
	11-01-2016		1	1				2
			1	1				2
Eckstein, Rhonda (eckstein)	<input type="checkbox"/> Scheduled	2	1					3
	11-04-2016	2						2
	11-07-2016		1					1
	<input checked="" type="checkbox"/> Cancelled							0
		2	1					3
Powell, Stick (stick)	<input checked="" type="checkbox"/> Scheduled	2	4	3	9	1	7	26
	<input checked="" type="checkbox"/> Cancelled							0
		2	4	3	9	1	7	26
Saltmarsh, Hilary (saltmarsh)	<input checked="" type="checkbox"/> Scheduled			1				1
	<input checked="" type="checkbox"/> Cancelled							0
				1				1
<b>Total</b>		4	6	5	9	1	7	32

Scheduler Activity

Confidential and Proprietary  
Page 1 of 1

Date Run: 11-08-2016 9:24 PM  
RRIS-RSNA\FORMS\GuestUser

### Scheduling Event Trend

Description: This report displays the counts of created exams broken down by hour for each day within the Date Range.

Parameters:

- Practice, Site(s) and Date Range


### Scheduling Event Trend

Practice: Treesdale

Site(s): Cedarville

Date Range: 11-07-2016 - 11-08-2016

Description: This report displays the counts of created exams broken down by hour for each day within the Date Range.



Date	6	7	8	9	10	11	12	13	14	15	16	20	Total
11-07-2016			20	31	27	30	20	26	43	25	12	1	235
11-08-2016	1	2	16	27	36	31	27	15	38	32	14		239
<b>Total</b>	<b>1</b>	<b>2</b>	<b>36</b>	<b>58</b>	<b>63</b>	<b>61</b>	<b>47</b>	<b>41</b>	<b>81</b>	<b>57</b>	<b>26</b>	<b>1</b>	<b>474</b>

Scheduling Event Trend

Confidential and Proprietary  
Page 1 of 1

Date Run: 11-08-2016 9:35 PM  
RADNET\EcksteinR

## Management Reports \ Statistics \ Technologists

### Technologist Activity by Modality

Description: This report displays the counts of performed, canceled and discontinued exams for each technologist by site and scanner.

Parameters:

- Practice, Site(s), Date Range and Technologist(s)

### Technologist Activity by Modality

**Practice:** Treesdale Radiology

**Site(s):** All sites for the selected practice(s).

**Date Range:** 10-01-2016 - 11-08-2016

**Description:** This report displays the counts of performed, canceled and discontinued exams for each technologist by site and scanner.

Technologist	Site	Scanner	# Performed	# Cancelled	# Discontinued
Aiken, Darcy (darcy)	Cedar	CT1CE	1		
		XR1CE	1		
		Site Totals	2	0	0
<b>Tech Totals</b>			<b>2</b>	<b>0</b>	<b>0</b>
Eckstein, Rhonda (eckstein)	Cedar	CT1CE	2		
		Site Totals	2	0	0
		<b>Tech Totals</b>			<b>2</b>
Powell, Stick (stick)	Maple	PT1MP		2	
		Site Totals	0	2	0
		<b>Tech Totals</b>			<b>0</b>
Waite, Andrew (andrew)	Maple	CT1MP	7		
		Site Totals	7	0	0
		<b>Tech Totals</b>			<b>7</b>
<b>Total</b>			<b>11</b>	<b>2</b>	<b>0</b>

Technologist Activity by Modality

Confidential and Proprietary  
Page 1 of 1

Date Run: 11-08-2016 9:31 PM  
RRIS-RSNA\FORMS\GuestUser

## Management Reports \ Utilization Management

### Coverage by Medical Group

Description: This will allow the user to select a Medical Group, Modality Type, and Coverage Type. In other words, it is a means to run a report against the UM authorization rules (matrices).

Parameters:

- Medical Group, Coverage and Modality Type(s)

### Coverage by Medical Group

**Medical Group:** Vantage Medical Group

**Modality Type(s):** CT,MR

**Coverage Code(s):** All coverages selected.

**Description:** This report shows the authorization rules by medical group.

Coverage	Description
A	Capitated service, authorization required
EX	Excluded
FFS	Fee for service, authorization required
FFS_N	Fee for service, authorization not required
X	Capitated service, direct referral (no authorization required)

CPT Code	Description	Coverage
<b>CAT Scan</b>		
70110	RADIOLOGIC EXAMINATION, MANDIBLE; COMPLETE, MINIMUM OF 4 VIEWS	X
70460	COMPUTED TOMOGRAPHY, HEAD OR BRAIN; WITH CONTRAST MATERIAL(S)	A
70470	COMPUTED TOMOGRAPHY, HEAD OR BRAIN; WITHOUT CONTRAST MATERIAL, FOLLOWED BY CONTRAST MATERIAL(S) AND FURTHER SECTIONS	A
70480	COMPUTED TOMOGRAPHY, ORBIT, SELLA, OR POSTERIOR FOSSA OR OUTER, MIDDLE, OR INNER EAR; WITHOUT CONTRAST MATERIAL	A
70481	COMPUTED TOMOGRAPHY, ORBIT, SELLA, OR POSTERIOR FOSSA OR OUTER, MIDDLE, OR INNER EAR; WITH CONTRAST MATERIAL(S)	A
70482	CT ORBIT SELLA/POST FOSSA/EAR C-/C+	X
70486	COMPUTED TOMOGRAPHY, MAXILLOFACIAL AREA; WITHOUT CONTRAST MATERIAL	A
70487	COMPUTED TOMOGRAPHY, MAXILLOFACIAL AREA; WITH CONTRAST MATERIAL(S)	A
70488	COMPUTED TOMOGRAPHY, MAXILLOFACIAL AREA; WITHOUT CONTRAST MATERIAL, FOLLOWED BY CONTRAST MATERIAL(S) AND FURTHER SECTIONS	A
70490	COMPUTED TOMOGRAPHY, SOFT TISSUE NECK; WITHOUT CONTRAST MATERIAL	A
70491	COMPUTED TOMOGRAPHY, SOFT TISSUE NECK; WITH CONTRAST MATERIAL(S)	A
70492	COMPUTED TOMOGRAPHY, SOFT TISSUE NECK; WITHOUT CONTRAST MATERIAL FOLLOWED BY CONTRAST MATERIAL(S) AND FURTHER SECTIONS	A

### Inbound Document Activity

Description: This report will allow the user to select a date range and receive the following information:

- The number of documents not processed and the average age in hours.
- The number of documents discarded, sorted by reason (e.g. Junk Mail, Duplicate, etc.).
- The number of documents processed by data entry user.
- A breakdown of incoming vs. processed documents by hour of day.


Parameters:

- From Date and To Date

### Inbound Document Activity

Date Range: 08-27-2016 - 11-25-2016

Description: This report shows a breakdown of inbound document activity based on processed date.



Document Status	Count	Percent	Age(Minutes)
<input type="checkbox"/> Discarded	7	30.43 %	29,977
Duplicate	3	42.86 %	29,756
Incomplete	1	14.29 %	40,177
Junk Mail	2	28.57 %	20,105
Multiple Patient	1	14.29 %	40,181
<input type="checkbox"/> Processed	16	69.57 %	13,655
spencer	16	100.00 %	13,655
<b>Total</b>	<b>23</b>		<b>18,622</b>

### Medical Group Activity

Description: This report will allow the user to select a date range and get a listing of all medical group activity that required utilization review and are in a final UM status for the period:

- A breakdown of STAT, Urgent, and Routine.
- Total procedure count, average turnaround hours, and percent compliant.
- Drill down capabilities from UM Status such as Recommended all the way down to Modality Type (CT, MR, etc.).

Parameters:

- From Date and To Date

### Medical Group Activity

Date Range: 01-01-2016 - 11-24-2016

Description: This report shows procedure count by medical group categorized by review priority, average turnaround time, and the percentage of reviews completed with the time allotted.

Medical Group	STAT			Urgent			Routine			Total		
	Count	Hours	Comp	Count	Hours	Comp	Count	Hours	Comp	Count	Hours	Comp
<input checked="" type="checkbox"/> Choice							1	1	100 %	1	1	100 %
<input checked="" type="checkbox"/> Prospect Medical Group							1	135	0 %	1	135	0 %
<input checked="" type="checkbox"/> Regal Medical Group							5	1	100 %	5	1	100 %
<input checked="" type="checkbox"/> Received							2	0	100 %	2	0	100 %
<input checked="" type="checkbox"/> Recommended							3	1	100 %	3	1	100 %
<input checked="" type="checkbox"/> As Added Exam							1	1	100 %	1	1	100 %
<input checked="" type="checkbox"/> As Alternative							1	1	100 %	1	1	100 %
<input checked="" type="checkbox"/> MRI and MRA							1	1	100 %	1	1	100 %
<input checked="" type="checkbox"/> As Requested							1	1	100 %	1	1	100 %
<input checked="" type="checkbox"/> Riverside Physician Network				3	1	100 %				3	1	100 %
<input checked="" type="checkbox"/> Seaview IPA	2	25	100 %	7	79	71 %	23	4	91 %	32	21	88 %
<input checked="" type="checkbox"/> Vantage Medical Group	21	0	100 %	52	35	75 %	217	20	88 %	290	21	87 %
<b>Total</b>	<b>23</b>	<b>2</b>	<b>100 %</b>	<b>62</b>	<b>38</b>	<b>76 %</b>	<b>247</b>	<b>19</b>	<b>88 %</b>	<b>332</b>	<b>21</b>	<b>87 %</b>

Medical Group Activity

Confidential and Proprietary  
Page 1 of 1

Date Run: 11-24-2016 2:53 PM  
SpencerLaptp-HPFormsGuestUser

### Referring Physician Activity

Description: This report will allow the user to select a date range and get a listing of all medical group activity that required utilization review and are in a final UM status for the period:

- A breakdown of STAT, Urgent, and Routine.
- Total procedure count, average turnaround hours, and percent compliant.
- Drill down capabilities from UM Status such as Recommended all the way down to Modality Type (CT, MR, etc.).

Parameters:

- From Date and To Date

### Referring Physician Activity

Date Range: 01-01-2016 - 11-24-2016

Description: This report shows utilization management procedure counts by referring physician categorized by review priority, average turnaround time, and the percentage of reviews completed with the time allotted.

Referring Physician	STAT			Urgent			Routine			Total		
	Count	Hours	Comp	Count	Hours	Comp	Count	Hours	Comp	Count	Hours	Comp
<input type="checkbox"/> Abbott, Alberto3 24490							8	10	88 %	8	10	88 %
<input type="checkbox"/> Not Recommended							2	0	100 %	2	0	100 %
<input type="checkbox"/> Recommended							5	16	80 %	5	16	80 %
<input type="checkbox"/> As Requested							5	16	80 %	5	16	80 %
<input type="checkbox"/> Diagnostic CT							4	20	75 %	4	20	75 %
<input type="checkbox"/> MRI and MRA							1	0	100 %	1	0	100 %
<input type="checkbox"/> Withdrawn							1	0	100 %	1	0	100 %
<input type="checkbox"/> James, Alfred0 88052							1	1	100 %	1	1	100 %
<input type="checkbox"/> Jefferson, Abraham6 45231							3	55	67 %	3	55	67 %
<input type="checkbox"/> Jones, Aimee 04700							2	68	50 %	2	68	50 %
<input type="checkbox"/> Jones, Allan 64317							2	0	100 %	2	0	100 %
<input type="checkbox"/> Jones, Artur 12418							1	0	100 %	1	0	100 %
<input type="checkbox"/> ref, test							1	0	100 %	1	0	100 %
<input type="checkbox"/> Refer, Self							3	0	100 %	3	0	100 %
<input type="checkbox"/> Smith, Alfons	14	4	100 %	47	48	70 %	186	21	89 %	247	25	86 %
<input type="checkbox"/> Smith, Alliso 97245				3	1	100 %	4	1	100 %	7	1	100 %
<input type="checkbox"/> Smith, Amelia0 89165							4	1	100 %	4	1	100 %
<input type="checkbox"/> Smith, Bobbie 46382	9	0	100 %	12	7	92 %	27	8	85 %	48	6	90 %
<input type="checkbox"/> Williams, Alexis 24384							5	25	60 %	5	25	60 %
<b>Total</b>	<b>23</b>	<b>2</b>	<b>100 %</b>	<b>62</b>	<b>38</b>	<b>76 %</b>	<b>247</b>	<b>19</b>	<b>88 %</b>	<b>332</b>	<b>21</b>	<b>87 %</b>

Referring Physician Activity

Confidential and Proprietary  
Page 1 of 1

Date Run: 11-24-2016 2:55 PM  
SpencerLapt-HP\FFormsGuestUser



### Utilization Review User Activity

Description: This report will allow the user to select a date range and get the following information:

- A listing of all the reviewers who finalized an authorization request during the period.
- Average hours to final status and percent of total.
- Ability to drill-down to see the UM Statuses such as Recommended As Requested.

Parameters:

- From Date and To Date

### Utilization Review User Activity

Date Range: 01-01-2016 - 11-24-2016

Description: This report shows utilization review procedure count by reviewer.

Review Finalized By	Proc Count	Avg Hours	% of Total
☐ MacDougall, Spencer	243	30	99.18 %
☒ MD Review	1	1	0.41 %
☐ Not Recommended	50	18	20.58 %
Medical Necessity Not Established	43	17	86.00 %
Not Medically Necessary	7	24	14.00 %
☒ Received	3	20	1.23 %
☒ Recommended	173	27	71.19 %
☒ Reopen	2	110	0.82 %
☒ Returned, Review Completed	3	41	1.23 %
☒ Returned, Review Not Initiated	1	299	0.41 %
☒ Withdrawn	10	108	4.12 %
☒ Mills, Terry	2	69	0.82 %
<b>Total</b>	<b>245</b>	<b>30</b>	

Utilization Review User Activity

Confidential and Proprietary  
Page 1 of 1

Date Run: 11-24-2016 3:10 PM  
SpencerLaptop-HP\FormsGuestUser

### UM Daily Dashboard

Description: This report is limited to one day’s activity and would typically be scheduled and distributed to the UM management team. It contains the following information:

- Number of unprocessed inbound documents in the data entry queue.
- Breakdown of the reviewer queue by priority and status.
- Finalized procedure count by priority and status.
- Finalized procedure count broken down by average turnaround time.
- Finalized procedure seven day trend bar chart.
- Breakdown of procedure count by Modality Type.
- Breakdown of procedure count by medical group.
- Top ten referring physician procedure count.

Parameters:

- Date

<b>Data Entry Queue</b>	25
-------------------------	----

Reviewer Queue	STAT	Urgent	Routine	Total	%
Hold	1	0	1	2	1.60%
MDReview	3	5	9	17	13.60%
NurseReview	0	0	2	2	1.60%
Received	12	8	80	100	80.00%
Reopen	0	0	4	4	3.20%
<b>Total</b>	<b>16</b>	<b>13</b>	<b>96</b>	<b>125</b>	

Final Status	Urgent	Routine	Total	%
NotRecommended	0	2	2	1.60%
Recommended	1	10	11	8.00%
ReturnedNotInitiated	1	0	1	0.80%
<b>Total</b>	<b>2</b>	<b>12</b>	<b>14</b>	

Turnaround Times	Urgent	Routine	Total	%
<= 10 Hours	0	2	2	14.29%
>= 31 Hours	2	10	12	85.71%
<b>Total</b>	<b>2</b>	<b>12</b>	<b>14</b>	
<b>Average</b>	<b>300.50</b>	<b>190.67</b>	<b>206.36</b>	

Modality Type	Total	%
Diagnostic CT	12	85.71%
MRI and MRA	2	14.29%
<b>Total</b>	<b>14</b>	

Medical Group	Total	%
Prospect Medical Group	1	7.14%
Vantage Medical Group	13	92.86%
<b>Total</b>	<b>14</b>	

Top 10 Requesting Physicians	Total
Smith, Alfons	9
Jones, Artur 12418	1
Jones, Aimee 04700	2
Jefferson, Abraham6 45231	1
Abbott, Alberto3 24490	1

