



eRAD RIS

RELEASE ANNOUNCEMENT

Build 4.2024.011

UPDATED JANUARY 30, 2024

TABLE OF CONTENTS

Summary 2

 Release Announcement 2

 New Known Limitations 2

Version Details 3

 Code Stream 3

PUBLICATION HISTORY

Revision	Author	Description
January 30, 2024	Kevin Brooks / Hilary Saltmarsh	▪ Commercial release.
January 31, 2024	Kevin Brooks	▪ Added newly identified Known Limitation BUG #35113.

SUMMARY

Release Announcement

This release of ERAD RIS 4.2024.011 merges the content of ERAD RIS 4.2024.008.1 which addresses a security issue and resolves a scanning issue which prevented attachments from being accessed (and includes ERAD RIS 4.2024.008.0.2 which introduces the ability to bypass Insurance for Mobile Workflow and resolves a URL site group parameter issue).

No additional functionality was added.

Refer to the previously released ERAD RIS CUSTOMER RELEASE NOTES for details of these changes.

New Known Limitations

The following new Known Limitations were identified that affects this ERAD RIS 4.2024.011 release (as well as the previously released HOTFIX AND RAPID RELEASES 3.2018.5.4.0.1 THROUGH 4.2023.076.1):

- **BUG #35113- NO PROVIDERS FOUND WHEN MAKING AN APPOINTMENT OUTSIDE OF LOG IN**
 - **Issue:** With the security fix hotfix in place, searching for a provider no longer returns any results in the guest appointment workflow.
 - **Impact:** When a patient who is not logged in (the guest appointment workflow), the provider search will not return any results, however they must make a provider selection to proceed.
 - **Workaround:** Patients that are logged in do not experience this limitation. Alternatively, they may select the `Provider Not Listed` option.



Carefully review all Known Limitations identified here and in the release notes of all other versions being applied and communicate relevant items to all those affected.

VERSION DETAILS

Code Stream

The following source code branches have been merged into this release:

