



eRAD RIS

RELEASE ANNOUNCEMENT

Rapid Release 4.2023.074.0.2

UPDATED SEPTEMBER 28, 2023

TABLE OF CONTENTS

Summary 2

 Release Announcement2

 Resolved Issues2

 New Known Limitations2

Version Details 3

 Code Stream3

PUBLICATION HISTORY

Revision	Author	Description
September 28, 2023	Kevin Brooks / Hilary Saltmarsh	▪ Commercial release.

SUMMARY

Release Announcement

This rapid release of eRAD RIS 4.2023.074.0.2 resolves issues related to C-CDA generation and transmission, but has no customer-facing changes.

Resolved Issues

This release resolves the following issues:

Category	Redmine #	Subject
Interfaces	34556	Resolves Interfaces issue with the Wedge Darena Plugin.
CURES	33328	Resolved CURES issue where C-CDA included visit-level data multiple times.

New Known Limitations

The following new Known Limitations were identified with this release:

- EXPECTED BEHAVIOR BY DESIGN (FROM #34556) - C-CDA TRANSMISSION FAILURE MONITORING
 - **Issue:** When the Wedge Darena Plugin fails to transmit a C-CDA, the failure is noted but there is no automated retry workflow.
 - **Impact:** When an exception occurs on the initial C-CDA submission to Darena, the text of the C-C-CDA is stored in the last_exception column of the c_action_queue row.
 - **Workaround:** Manual monitoring for errors is required. If a failure is identified, details should be forwarded to Darena for review.

VERSION DETAILS

Code Stream

The following source code branches have been merged into this release:

