



# eRAD RIS

## RELEASE ANNOUNCEMENT

### Build 4.2024.059

UPDATED MAY 9, 2024

#### TABLE OF CONTENTS

**Summary** ..... 2

    Release Announcement..... 2

    New Features..... 2

**Feature Details** ..... 3

    System Config..... 3

    Feature #34737 – Ability to set Web Viewer as the Default Image Viewer..... 3

**Version Details** ..... 4

    Code Stream..... 4

#### PUBLICATION HISTORY

Revision	Author	Description
May 9, 2024	Kevin Brooks, Michelle Mahabir / Hilary Saltmarsh	▪ Commercial release.

# SUMMARY

## Release Announcement

This release of *ERAD RIS 4.2024.059* introduces a System Config enhancement that allows the Web Viewer to be set as the Default Image Viewer for new Referring Portal accounts and merges the content of release *ERAD RIS 4.2024.056.0.1* (adds a service-only feature with new hosting options for Online Registration application but has no customer-facing changes) into this version.

Refer to the previously released *ERAD RIS CUSTOMER RELEASE NOTES* for details of these changes.

## New Features

This release introduces the following features and enhancements:

Category	Redmine #	Subject	Description
System Config	34747	Ability to set Web Viewer as the Default Image Viewer	This enhancement to System Config adds the ability to set Web Viewer as a Default Image Viewer for Referring Portal accounts.

SORTED BY CATEGORY AND REDMINE

Refer to the **FEATURE DETAILS** section below for configuration and usage information.

# FEATURE DETAILS

## System Config

### Feature #34737 – Ability to set Web Viewer as the Default Image Viewer

#### Summary

This enhancement to System Config adds the ability to set Web Viewer as the Default Image Viewer for Referring Portal accounts.


#### Background

Previously, the Full Viewer was the only option available as a default viewer when creating new Referring Portal accounts.

#### Feature Description

With this change, a new `RPDefaultViewer` appconfig setting adds the capability to set the user's default viewer to either the Full Viewer or the Web Viewer.

This appconfig can be overridden at the practice level, however, in instances where the practice level configuration is null, the system will default to the appconfig with the value of FullViewer or configured value.



The screenshot shows a navigation bar with tabs: My Profile, Settings, Security, Staff Settings, and Image Viewer. Below the tabs, the 'Image Viewer' section is active. It contains a form with a label 'Default Image Viewer' and a dropdown menu currently set to 'Web Viewer'. Below the dropdown are two buttons: 'SAVE' (dark blue) and 'CANCEL' (light blue). The entire form area is highlighted with a red border.

DEFAULT VIEWER FIELD IN ADMIN PORTAL

#### Configuration Instructions

No System Administrator actions are necessary to enable this feature; however, optional configuration is available:

#### Portals

##### Changes to applicationsettings.config

- Set `RPDefaultViewer` configuration to desired default image viewer.
- Setting can be overridden at the Practice level.

# VERSION DETAILS

## Code Stream

The following source code branches have been merged into this release:

