



eRAD RIS

RELEASE ANNOUNCEMENT

Build 4.2023.051

UPDATED JULY 10, 2023

TABLE OF CONTENTS

Summary 2

- Release Announcement2
- New Features2
- Resolved Issues2

Feature Details 3

- Lookup Tables3
 - Feature #31727 – Support sending scanned documents to PACS.....3
- Worklists3
 - Feature #33991 – Add Alerts Column to the Clinical Task WL.....3

PUBLICATION HISTORY

Revision	Author	Description
July 10, 2023	Kevin Brooks, Michelle Mahabir / Hilary Saltmarsh	<ul style="list-style-type: none"> ▪ Commercial release.

SUMMARY

Release Announcement

This release of ERAD RIS 4.2023.051 resolves issues relating to Insurance Eligibility, DB, and System Configurations and introduces two workflow related enhancements.

New Features

This release introduces the following features and enhancements:

Category	Redmine #	Subject	Description
Lookup Tables	31727	Support sending scanned documents to PACS	This enhancement adds ability to identify scanned documents that can be sent to PACS in conjunction with a third party application.
Worklists	33991	Add Alerts column to the Clinical Task Worklist	This enhancement adds the Alerts column to the Clinical Task worklist, consistent with other worklists.

SORTED BY CATEGORY AND REDMINE

Resolved Issues

This release resolves the following issues:

Category	Redmine #	Subject
DB	33193	This is a service-only change with no customer-facing functionality. Resolved DB issue where three missing indexes were identified.
Insurance Eligibility	33239	This is a service-only change with no customer-facing functionality. Resolved Insurance Eligibility issue that was presenting deadlocks due to missing code in the InsuranceEligibilityService.cs.
System Config	32676	Resolved System Config issue that was incorrectly prompting the user with the Zip Code/County Selection and Primary Chair selection when in View Study mode.

SORTED BY CATEGORY AND REDMINE

FEATURE DETAILS

Lookup Tables

Feature #31727 – Support sending scanned documents to PACS

Summary

This enhancement to Lookup Tables adds the ability to identify scanned documents that can be sent to PACS in conjunction with a third party application such as ScanLink.

Feature Description

With this change, a new column has been added to the `ScanType` Lookup Table to assist the system with sending scanned documents to PACS as DICOM objects.

This can be used in conjunction with the ScanLink application and will allow control over which documents are permitted to be sent.

Configuration Instructions

System Administrators must complete the following actions to enable this feature:

RIS Client

Changes to ScanType Lookup Table Settings

- Enable this feature by setting `SendToPACSFlag` to `Y`.
 - Allowed values consist of `Y` and `N`
 - Default value will be set to `N`

Worklists

Feature #33991 – Add Alerts Column to the Clinical Task WL

Feature Description

The `Alert` column has been added to the `Clinical Task WL`, consistent with other worklists. This enhancement will assist with current and future workflows.

Alerts	Flags	Status	Procedures
NotNull:			
		Scheduled	US Abdominal Ltd Or Fu [76705]
		Report Drafted	XR CERVICAL SPINE MIN 4 VIEWS [7
		Ordered	CT Ear Orbit Sella W [70481] - Head - f
		Exam Done	XR SHOULDER MIN 2 VIEWS [73030]
		Signed	MRA Spine W/Wo [72159]
		Signed	CT Head W/Wo & Sinus W/Wo [CT15] -
		Signed	CT Chest W/o & Abdomen W/o [CT1] - C
		Ordered	CT Cervical Spine W [721261] - Cervica

Patient-level alert:
Iodine Allergy
Updated: 08-13-2018 by Amy MacPherson

Configuration Instructions

No System Administrator actions are necessary to enable this feature.