



eRAD RIS

RELEASE ANNOUNCEMENT

Build 4.2023.017

UPDATED APRIL 7, 2023

PUBLICATION HISTORY

Revision	Author	Description
April 7, 2023	Kevin Brooks, Michelle Mahabir / Hilary Saltmarsh	<ul style="list-style-type: none">Commercial release.

SUMMARY

Release Announcement

This release of ERAD RIS 4.2023.017 resolves a read-only field concern for digital form mapping.

New Features

No new features were introduced with this release.

Resolved Issues

This release resolves the following issues:

Category	Redmine #	Subject
Digital Forms	32759	Resolved a digital form issue that was causing answer fields to be read-only for questions mapped to portal scheduling questions.

New Known Limitations

While no new known limitations were identified with this release, the following expected behaviors were noted:

- EXPECTED BEHAVIOR BY DESIGN
 - Previously, when an answer on the Connect Portal tab was modified and saved, the answer would be reflected on the (read-only) digital form. With this change, values now may be modified in either location, however changes made in one location will not be reflected in the other.