



eRAD RIS

RELEASE ANNOUNCEMENT (SERVICE TEAM EDITION)

Rapid Release 4.2023.106.0.2

UPDATED DECEMBER 27, 2023

TABLE OF CONTENTS

Summary 2

 Release Announcement2

 New Features2

Feature Details 3

 Patient WF: Make Appointment3

 Feature #34832 - Hide site address and skip prescription workflow step in
 MobileMammo workflow3

Version Details 5

 Code Stream5

PUBLICATION HISTORY

| Revision | Author | Description |
|-------------------|------------------------------------|---|
| December 27, 2023 | Kevin Brooks / Hilary Saltmarsh | <ul style="list-style-type: none"> ▪ Commercial release. |

SUMMARY

Release Announcement

This rapid release of ERAD RIS 4.2023.106.0.2 updates the Mobile Mammo workflow to optionally skip the prescription information workflow step and optionally hide site address information for selected sites.

New Features

This release introduces the following features and enhancements:

| Category | Redmine # | Subject | Description |
|------------------------------|-----------|---|---|
| Patient WF: Make Appointment | 34832 | Hide site address and skip prescription workflow step in MobileMammo workflow | This enhancement to Patient WF: Make Appointment updates the Mobile Mammo workflow to optionally skip the prescription information workflow step and optionally hide site address information for selected sites. |

FEATURE DETAILS

Patient WF: Make Appointment

Feature #34832 - Hide site address and skip prescription workflow step in MobileMammo workflow

Summary

This enhancement to Patient WF: Make Appointment updates the Mobile Mammo workflow to optionally skip the prescription information workflow step and optionally hide site address information for selected sites.


Feature Description

With this change, the page presented to users scheduling from the the Mobile Mammo workflow will now skip the prescription information workflow step by default.



Note that this functionality is optional but enabled by default.

Additionally, site address information may now be suppressed for selected sites.



Service Team Upgrade Instructions



Additional technical details and deploy instructions are available in the SERVICE TEAM EDITION of this document.

The Service Team must complete the following actions to deploy this feature:

Web

- To **disable** the skip prescription information workflow step, check `Show Prescription Step` via the admin URL Generation Tool, located at:
`{Portal URL}/Admin/adminsystem/MobileScheduleURL`

GENERATE RESTRICTED SCHEDULE URL

Site_Group

Show Prescription Step

StartDate *

Start time *

EndDate *

End time *

SUBMIT

Web Services

Changes to applicationsettings.config

- To suppress site address information, update the new `SitesForMobileMammoToHideAddress` with a list of sites in CSV format, e.g., `LU,EL,MA`

Configuration Instructions

Service Team assistance is required to enable this feature.

VERSION DETAILS

Code Stream

The following source code branches have been merged into this release:

