

# **Server Update**

**For eRAD RIS  
Version 1.0  
Build 38  
Update 38.006**

## **Table of Contents**

1. Purpose.....	3
Defects resolved.....	3
2. Package Contents .....	3
3. Intended Audience .....	4
4. Who is affected .....	4
5. eRAD RIS Release Version Numbers .....	5
6. Installing .....	6

## **RADNET, Inc. – eRAD RIS Server Update**

### **1. Purpose**

This server update has been produced to fix bugs required for production setting, or to complete an enhancement requested for production.

This server update includes the following:

#### **Defects resolved**

Redmine #	Subject	Priority
2869	Reporting - Interactive dictation getting lost and set to 0 bytes	Urgent
2876	Transcription Dashboard not reporting self-editing radiologists properly	Normal

### **2. Package Contents**

Files included for this server update:

1. GUI ZIP Updates folder contains
  - a. AudioLib.dll
  - b. AudioLib.pdb
  - c. Ris.exe
  - d. Ris.pdb
2. SQL folder contains
  - a. B38-006.sql

### **3. Intended Audience**

This document is created by the RIS Development team for the RadNet RIS Implementation team.

### **4. Who is affected**

Previous versions of build 38.

Server updates may not be cumulative; therefore this server updated must be applied directly to the previous server update.

## 5. eRAD RIS Release Version Numbers

Build	Patch	UI Version	Core Version	WS Version	DB Version	Notes
37	-	1.2.37.12843	1.2.37.12782	1.2.37.12843	1.2.37.12843	
37	1	1.2.37.13040	1.2.37.13040	1.2.37.13039	1.2.37.13003	
37	2	1.2.37.13069	1.2.37.13040	1.2.37.13069	1.2.37.13003	
37	3	1.2.37.13069	1.2.37.13040	1.2.37.13069	1.2.37.13003	Just WebAPI updated, no version changes
37	4	1.2.37.13171	1.2.37.13040	1.2.37.13165	1.2.37.13171	
37	5	1.2.37.13347	1.2.37.13040	1.2.37.12979	1.2.37.12979	Bug was discovered with version numbers, this is why it appears we went backwards in the version number but we really did not.
37	6	1.2.37.13665	1.2.37.13040	1.2.37.13580	1.2.37.12979	
37	7	1.2.37.13807	1.2.37.13040	1.2.37.13799	1.2.37.13807	
37	8	1.2.37.13807	1.2.37.13040	1.2.37.13827	1.2.37.13807	
37	9	1.2.37.13807	1.2.37.13040	1.2.37.13839	1.2.37.13807	
37	10	1.2.37.13895	1.2.37.13040	1.2.37.13895	1.2.37.13895	
37	11	1.2.37.14161	1.2.37.13040	1.2.37.14161	1.2.37.14161	
37	12	1.2.37.14267	1.2.37.13040	1.2.37.14161	1.2.37.14267	
37	13	1.2.37.14326	1.2.37.13040	1.2.37.14326	1.2.37.14326	
38	-	1.2.38.14486	1.2.38.14486	1.2.38.14486	1.2.38.14486	
38	1	1.2.38.14748	1.2.38.14748	1.2.38.14748	1.2.38.14748	
38	2	1.2.38.15182	1.2.38.14748	1.2.38.14748	1.2.38.14748	
38	3	1.2.38.15251	1.2.38.14748	1.2.38.14748	1.2.38.14748	
38	4	1.2.38.15327	1.2.38.14748	1.2.38.15290	1.2.38.14748	
38	5	1.2.38.15513	1.2.38.14748	1.2.38.15513	1.2.38.15451	
38	6	1.2.38.15618	1.2.38.14748	1.2.38.15513	1.2.38.15618	

## **6. Installing**

- 1) Copy the files in the GUI ZIP Updates folder and replace those in the rRIS\_1.2.38.15513.zip folder. Rename the zip file to be **rRIS\_1.2.38.15618.zip**
- 2) Run the B38\_006.sql file against the eRAD RIS database instance.