

Server Update

for eRAD RIS
Version 1.0
Build 36
Update 36.005

Table of Contents

1. Purpose.....	3
2. Intended Audience	4
3. Who is affected	4
4. Installing	4

1. Purpose

This server update has been produced to fix bugs required for production setting, or to complete an enhancement requested for production.

This server update fixes 2 issues.

- 1.) Redmine Bug 1673 – setBuillingStatus NA action causing billing post to fail.
When completing a TechOnly exam that goes to billing, an error was occurring because it was being set at not applicable. It has been changed back to be set as UI_StudyCompleted

Files for this fix are included in the following folders.

- 1.) Folder GUI ZIP Updates contains the following
 - a. ris.exe
- 2.) Folder RIS Service Updates contains the following
 - a. WCFService.dll
 - b. WCFService.pdb

2. Intended Audience

This document is created by the RIS Development team for the RadNet RIS Implementation team.

3. Who is affected

All previous versions of build 36.

NOTE: Previous server updates and hotfixes must be applied before this update.

4. Installing

- 1.) Copy and replace the files **from** the GUI ZIP Updates folder **to** the rRIS_1.1.36.12653.zip file. Rename the .zip to rRIS_1.1.36.12953.zip
(NOTE: the .zip file name must be different for the changes to take effect on the client. Changing it to .12953 will keep it in sync with the RIS UI Version)
- 2.) Stop the following services as needed: copy the files from the RIS Service Updates folder to rRISService folder and restart the service(s).
These directories may include:
 - a. rRISService (core service)
 - b. CDS rRISService
 - c. Wedge rRISService
 - d. DDE rRISService

(NOTE: your installation may not have all services listed)

