



eRAD RIS

RELEASE ANNOUNCEMENT

Build 4.2024.076

UPDATED JULY 8, 2024

TABLE OF CONTENTS

Summary 2

 Release Announcement2

 New Features2

 Resolved Issues2

Feature Details 3

 Extra Data3

 Service-Only Feature #35559 – Add Site Extra Info values to outbound Wedge message3

Version Details 4

 Code Stream4

PUBLICATION HISTORY

Revision	Author	Description
July 8, 2024	Kevin Brooks, Michelle Mahabir / Hilary Saltmarsh	▪ Commercial release.

SUMMARY

Release Announcement

This release of ERAD RIS 4.2024.076 introduces a service-only update to outbound Wedge messages, resolves filtering performance within Materialized Worklist and merges the content of ERAD RIS 4.2024.072.1 (which resolves an error loading orders in the IVT workflow and an error accessing the Historical Carrier dialog).

New Features

This release introduces the following features and enhancements:

Category	Redmine #	Subject	Description
Extra Data	35559	Add Site Extra Info values to outbound Wedge message	<p>This is a service-only change with no customer-facing functionality.</p> <p>This enhancement to Extra Data adds Site Extra Info values to outbound Wedge messages to remain in compliance with FDA regulations.</p>

Refer to the FEATURE DETAILS section below for configuration and usage information.

Resolved Issues

This release resolves the following issues:

Category	Redmine #	Subject
Performance	35469	Resolved Performance issues when utilizing filters in Materialized Worklists.

FEATURE DETAILS

Extra Data

Service-Only Feature #35559 – Add Site Extra Info values to outbound Wedge message

Summary

This enhancement to Extra Data adds Site Extra Info values to outbound Wedge messages to remain in compliance with FDA regulations.

There are no customer-facing changes.

Feature Description

With this change, sites that have multiple extra info codes will now include all in the message individually.

In the outbound message resolver, when expanding site_code and the columns below it, a sub node for reach l_site_extra_info attached to the site is now included, e.g.:

```
<c_study>
  <scheduled_modality_code>
  ...
  <site_code>
    <l_site>
      <site_code>LU</site_code>
    </l_site>
    <l_site_extra_info>
      ...
    </l_site_extra_info>
    <l_site_extra_info>
      ...
    </l_site_extra_info>
  </site_code>
</c_study>
```

Service Team Upgrade Instructions

No Service Team actions are necessary to deploy this feature.

Configuration Instructions

No System Administrator actions are necessary to enable this feature.

VERSION DETAILS

Code Stream

The following source code branches have been merged into this release:

