



eRAD RIS

# RELEASE ANNOUNCEMENT

Build 4.2023.060

UPDATED AUGUST 4, 2023

## TABLE OF CONTENTS

<b>Summary</b> .....	<b>2</b>
Release Announcement .....	2
Resolved Issues .....	2
New Known Limitations .....	2
<b>Version Details</b> .....	<b>3</b>
Code Stream .....	3

## PUBLICATION HISTORY

Revision	Author	Description
August 4, 2023	Kevin Brooks / Hilary Saltmarsh	▪ Commercial release.
August 4, 2023	Kevin Brooks	▪ Clarification to Known Limitation

# SUMMARY

## Release Announcement

*This release of eRAD RIS 4.2023.060 resolves several Insurance Management issues.*

## Resolved Issues

This release resolves the following issues:

Category	Redmine #	Subject
Insurance Management	34310	Resolved Insurance Management issue where patient insurance information was lost when scheduling an order.
Insurance Management	34319	Resolved Insurance Management issue where viewing or associating an insurance card intermittently returned an error.

## New Known Limitations

The following new Known Limitations were identified with this release:

- **BUG #34337 - AFTER RE-PROCESS ELIGIBILITY THE WL WARNING FLAG ICON DOES NOT DISPLAY AND IMPACTS NEXT ATTEMPT TO RE-PROCESS ELIGIBILITY**
  - **Issue:** When re-processing eligibility for carriers that do not go thru Availity, the worklist Warning Flag icon does not display and prevents further attempts to re-process eligibility. This issue does not happen with Availity based carriers.
  - **Impact:** RIS will indicate "There are no items matching the specified criteria." and the shield will no longer appear on the worklist. **Note:** This issue does not happen with Availity based carriers.
  - **Workaround:** By opening the order in View/Edit and re-saving, the shield will reappear and reprocessing will work as normal.

# VERSION DETAILS

## Code Stream

The following source code branches have been merged into this release:

