



eRAD RIS

# CUSTOMER RELEASE NOTES

Hotfix Build v3.2021.5.24.2

UPDATED JULY 26, 2021

# Table of Contents

<b>Summary</b> .....	<b>3</b>
Intended Audience .....	3
Who Is Affected .....	3
<b>New Settings</b> .....	<b>4</b>
<b>New Features</b> .....	<b>5</b>
<b>Resolved Issues and Known Limitations</b> .....	<b>6</b>
Resolved Issues .....	6
New Known Limitations .....	6
<b>Version Details</b> .....	<b>7</b>
Package Contents .....	7
Code Stream .....	7
eRAD RIS Release Version Numbers .....	8

## Publication History

Revision	Author	Description
July, 26, 2021	Kevin Brooks / Hilary Saltmarsh	<ul style="list-style-type: none"> <li>▪ Commercial release.</li> </ul>

# SUMMARY

## Intended Audience

The intended audience for this CUSTOMER RELEASE NOTES document is the RIS Administration team for all eRAD RIS customers.

This document describes the purpose, configuration, and operation of new features made available with this release, identifies issues resolved in the release, and highlights any unresolved known limitations.

This information should be used by all customers to determine how these changes affect their organization's workflow, and to plan their organization's upgrade strategy for eRAD RIS.



*There is no separate SERVICE RELEASE NOTES edition of this document.*

## Who Is Affected

This hotfix release of eRAD RIS v3.2021.5.24.2 resolves several urgent issues, including multiple Scheduling issues. This version is recommended to be applied on eRAD RIS v3.2018.5.24.1 installations.



*Plsased carefully review these release notes even if your system will not be upgraded immediately, to identify and communicate any issues that may affect your organization.*

# NEW SETTINGS

No changes in this release.

# NEW FEATURES

No features in this release.

# RESOLVED ISSUES AND KNOWN LIMITATIONS

## Resolved Issues

This release resolves the following issues:

Redmine #	Subject
30262	Resolved database issue where the Performed Modality Code was incorrectly updated for Scheduled exams.
30308	Resolved Scheduling issue where invalid phone numbers were being deleted and then returning an error.
30319	Resolved Patient Portal issue where blank phone numbers save as a question mark.
30322	Resolved Digital Forms issue where saving study-level forms without a study returned an error.
30345	Resolved EMR issue where a missing placer_order_number resulted in a Scheduling Error.
30358	Resolved EMR issue where choosing a Procedure Plan for an EMR order resulted in a Scheduling Error.
30368	Resolved DB performance issue where database indexing was causing a backlog on the c_action_queue.















## New Known Limitations

No new known limitations were identified with this release.

# VERSION DETAILS


## Package Contents

The release package includes the following folders:

 _ReleaseNotes	7/23/2021 4:41 PM
 Client Application	7/23/2021 2:55 PM
 DB	7/23/2021 2:55 PM
 External WebAPI	7/23/2021 2:55 PM
 Identity Service	7/23/2021 2:55 PM
 PACS Citrix Bridge	7/23/2021 2:56 PM
 RIS Service	7/23/2021 2:55 PM
 Service Tools	7/23/2021 2:57 PM
 Web Digital Forms	7/23/2021 2:55 PM
 Web Patient Connect	7/23/2021 2:56 PM
 Web Referring Connect	7/23/2021 2:56 PM
 Web UM Connect	7/23/2021 2:56 PM
 Build_2021.5.24.2.zip	7/23/2021 2:58 PM
 RISServerMasterCert.pfx	3/31/2016 1:38 PM

## Code Stream

The following source code branches have been merged into this release:



~\* Previous years patches released in 2021  
 ~ Previous years patches released in 2020  
 \*\*\* = previous years patches released in 2019

LEGEND:

- Light Green = Previously Released software
- Gray = Internal version, non-release version
- Bright Green = Current Release

## eRAD RIS Release Version Numbers

The following table details the version identifiers for components in this release:



Build	Patch	UI Version	Core Version	WS Version	DB Version	Digital Forms	Patient Portal	UM Portal	Provider Portal	Notes
2018.5	6	3.18.5.6(3GB)	3.18.5.6	3.18.5.6	3.18.5.6.0.02571320	3.18.5.6.0	3.18.5.6.0.1050	3.18.5.6.0.1050	3.18.5.6.0.1050	Full version release
2018.5	6.1	3.18.5.6.1(3GB)	3.18.5.6.1	3.18.5.6.1	3.18.5.6.1.02601339	3.18.5.6.1	3.18.5.6.1.1070	3.18.5.6.1.1070	3.18.5.6.1.1070	GUI, Web Services, DB, Patient/Provider/UM Portals and Digital Forms
2018.5	6.2	3.18.5.6.2(3GB)	3.18.5.6.2							GUI
2018.5	6.3	3.18.5.6.3(3GB)	3.18.5.6.3	3.18.5.6.3	3.21.5.10.002823110					GUI and Webservices
2020.12.21	-	3.20.12.21.0(3GB)	3.18.5.6	3.20.12.21.0	3.20.12.21.0.02608693	3.20.12.21.0	3.20.12.21.0.897	3.20.12.21.0.897	3.20.12.21.0.897	Full version release
2021.1.4	-	3.21.1.4.0(3GB)	3.21.1.4				3.21.1.4.0.1075	3.21.1.4.0.1075	3.21.1.4.0.1075	GUI and Patient/Provider/UM portals
2021.1.18	-	3.21.1.18.0(3GB)	3.21.1.18	3.21.1.18.0	3.21.1.18.002652334	3.21.1.18.0	3.21.1.18.0.1081	3.21.1.18.0.1081	3.21.1.18.0.1081	GUI, Web Services, DB, Patient/Provider/UM portals and Digital Forms
2021.2.1	-	3.21.2.1.0(3GB)	3.21.2.1	3.21.2.1.0	3.21.2.1.002672074	3.21.2.1.0	3.21.2.1.0.1090	3.21.2.1.0.1090	3.21.2.1.0.1090	GUI, Web Services, DB, Patient/Provider/UM portals and Digital Forms
2021.2.15	-	3.21.2.15.0(3GB)	3.21.2.15	3.21.2.15.0	3.21.2.15.002698266	3.21.2.15.0	3.21.2.15.0.1111	3.21.2.15.0.1111	3.21.2.15.0.1111	GUI, Web Services, DB, Patient/Provider/UM portals and Digital Forms
2021.2.15	1	3.21.2.15.1(3GB)	3.21.2.15.1	3.21.2.15.1		3.21.2.15.1	3.21.2.15.1.1113	3.21.2.15.1.1113	3.21.2.15.1.1113	GUI, Web Services, Patient/Provider/UM portals and Digital Forms
2021.2.15	2	3.21.2.15.2(3GB)	3.21.2.15.2	3.21.2.15.2		3.21.2.15.2	3.21.2.15.2.1122	3.21.2.15.2.1122	3.21.2.15.2.1122	GUI, Web Services, Patient/Provider/UM portals and Digital Forms
2021.2.15	3	3.21.2.15.3(3GB)	3.21.2.15.3							GUI
2021.2.15	4	3.21.2.15.4(3GB)	3.21.2.15.4	3.21.2.15.4	3.21.2.15.402732838	3.21.2.15.4	3.21.2.15.4.1133	3.21.2.15.4.1133	3.21.2.15.4.1133	GUI, Web Services, DB, Patient/Provider/UM portals and Digital Forms
2021.2.15	5	3.21.2.15.5(3GB)	3.21.2.15.5	3.21.2.15.5	3.21.2.15.502738557					GUI, Web Services and DB
2021.2.15	6	3.21.2.15.6(3GB)	3.21.2.15.6	3.21.2.15.6						GUI and Web Services
2021.3.1	-	3.21.3.1(3GB)	3.21.3.1	3.21.3.1	3.21.3.1.002712308	3.21.3.1	3.21.3.1.0.1119	3.21.3.1.0.1119	3.21.3.1.0.1119	GUI, Web Services, DB, Patient/Provider/UM portals and Digital Forms
2021.3.15	-	3.21.3.15(3GB)	3.21.3.15	3.21.3.15	3.21.3.15.002747230	3.21.3.15	3.21.3.15.0.1144	3.21.3.15.0.1144	3.21.3.15.0.1144	GUI, Web Services, DB, Patient/Provider/UM portals and Digital Forms
2021.4.12	-	3.21.4.12(3GB)	3.21.4.12	3.21.4.12	3.21.4.12.002778929	3.21.4.12	3.21.4.12.0.1167	3.21.4.12.0.1167	3.21.4.12.0.1167	GUI, Web Services, DB, Patient/Provider/UM portals and Digital Forms
2021.5.10	-	3.21.5.10(3GB)	3.21.5.10	3.21.5.10	3.21.5.10.002823110	3.21.5.10	3.21.5.10.0.1189	3.21.5.10.0.1189	3.21.5.10.0.1189	GUI, Web Services, DB, Patient/Provider/UM portals and Digital Forms
2021.5.24	-	3.21.5.24(3GB)	3.21.5.24	3.21.5.24	3.21.5.24.002842120	3.21.5.24	3.21.5.24.0.1205	3.21.5.24.0.1205	3.21.5.24.0.1205	GUI, Web Services, DB, Patient/Provider/UM portals and Digital Forms
2021.5.24	1	3.21.5.24.1(3GB)	3.21.5.24.1	3.21.5.24.1	3.21.5.24.102880695	3.21.5.24.1	3.21.5.24.1.1224	3.21.5.24.1.1224	3.21.5.24.1.1224	GUI, Web Services, DB, Patient/Provider/UM portals and Digital Forms
2021.5.24	2	3.21.5.24.2(3GB)	3.21.5.24.2	3.21.5.24.2	3.21.5.24.202923248	3.21.5.24.2	3.21.5.24.0.1244	3.21.5.24.0.1244	3.21.5.24.0.1244	GUI, Web Services, DB, Patient/Provider/UM portals and Digital Forms
2021.6.7	-	3.21.6.7(3GB)	3.21.6.7	3.21.6.7	3.21.6.7.002898418	3.21.6.7	3.21.6.7.0.1234	3.21.6.7.0.1234	3.21.6.7.0.1234	GUI, Web Services, DB, Patient/Provider/UM portals and Digital Forms