

Server Update

For eRAD RIS
Version 1.0
Build 40
Update 40.008

Table of Contents

| | |
|--|----|
| 1. New Features | 3 |
| Reporting - Order Notes Added To Exam Details Pane | 3 |
| Reporting – Alter and Save Default Layouts | 4 |
| 2. Defects Resolved | 6 |
| 2. Package Contents | 7 |
| 3. Intended Audience | 8 |
| 4. Who is affected | 8 |
| 5. eRAD RIS Release Version Numbers | 9 |
| 6. Installing | 10 |

1. New Features

Reporting - Order Notes Added To Exam Details Pane

Order Notes were missing from reporting window. Previously the users would have to access the Order Notes from the data nugget or the View / Edit window on the Patient Folder WL.




Figure 1.1 – Order Notes in Exam Details pane

In Server release 1.40.008, the Order notes have been added to the Exam Details data pane in the reporting window. Users must be aware that in order to view the Order notes they must dock the Exam Details data pane and save layout for it to display in future reporting window sessions.

Reporting – Alter and Save Default Layouts

In the reporting windows view menu users had the option to choose from of a list of default layouts. These layouts came “canned” with eRAD RIS and if changes were desired to the layout, they could not be saved.




Figure 1.2 – View Menu - Load Layout

It was requested that these layouts could be altered and changes saved. We have added a new menu and sub-menu to the View menu. Under the report editor screens View menu, administrators will find a new menu item. **Administer Layouts...** this will lead to a **Save As** option for each of the standard layouts. The administrator simply designs an open screen as they desire then use this menu item to save to the database as one of the default layouts for all users. Availability of this set of menu items is based on a **new access string “View.Administration”**. Users must have full access to the mentioned access string in order to use this feature.

RADNET, Inc. – eRAD RIS Server Update




Figure 1.3 – View Menu – Administer Layouts

2. Defects Resolved

| Redmine # | Subject | Priority |
|-----------|---|----------|
| 3696 | Attachments - Wrong Patient Key Data | Urgent |
| 3422 | On the Order tab, "Visit location for" field value is incorrect for the physician that is highlighted | High |
| 3614 | Cancelled orders not appearing on Patient Folder | High |
| 3693 | Scheduling availability templates rules - times can appear as open in schedule search | High |
| 3713 | Open/Wide Board Modalities are difficult to identify from search results | High |
| 3715 | Insert macro/template voice commands not working. | High |
| 3718 | Error applying default view to patient folder when patient folder is pinned. | High |
| 3720 | Prompt for more pages does cause issues if multiple documents are scanned before the prompt comes up | High |
| 3737 | the binary string truncated error | High |
| 3636 | IVT WL showing cancelled studies as scheduled | Normal |
| 3672 | External Report checkbox does not stay checked | Normal |
| 3673 | Reporting - reports of double dictation text intermittently | Normal |
| 3678 | Billing messages to correctly specify the correct site code when the study is an outside read | Normal |
| 3694 | booking locks not released in b40 when cancelling studies in some cases | Normal |
| 3695 | booking locks not released when moving a study to a new order | Normal |
| 3702 | override templates are not painting in calendar control | Normal |

2. Package Contents

Files included for this server update:

- 1) GUI ZIP Updates folder contains
 - a) rRIS_1.3.40.17364.zip
- 2) RIS Service Updates folder contains
 - a) CDSJobManager.dll
 - b) CDSJobManager.pdb
 - c) CommonLib.dll
 - d) CommonLib.pdb
 - e) DDEngine.dll
 - f) DDEngine.pdb
 - g) rRISServices.exe
 - h) rRISServices.pdb
 - i) rRISServices.vshost.exe
 - j) ScheduledTaskManager.dll
 - k) ScheduledTaskManager.pdb
 - l) SignatureApprovalService.exe
 - m) SignatureApprovalService.pdb
 - n) WCFService.dll
 - o) WCFService.pdb
 - p) WebApi.dll
 - q) WebApi.pdb
 - r) Wedge.exe
 - s) Wedge.pdb
- 3) SQL folder contains:
 - a) B40_008.sql

3. Intended Audience

This document is created by the RIS Development team for the RadNet RIS Implementation team.

4. Who is affected

Build 40 installs. This server update must be applied directly to the build 40 core install. Server updates are not cumulative. Any previous server updates must be applied before this one.

5. eRAD RIS Release Version Numbers

| Build | Patch | UI Version | Core Version | WS Version | DB Version | Notes |
|-------|-------|--------------|--------------|--------------|--------------|---|
| 37 | - | 1.2.37.12843 | 1.2.37.12782 | 1.2.37.12843 | 1.2.37.12843 | Full version release |
| 37 | 1 | 1.2.37.13040 | 1.2.37.13040 | 1.2.37.13039 | 1.2.37.13003 | Core install, GUI.zip, web service and database updates |
| 37 | 2 | 1.2.37.13069 | 1.2.37.13040 | 1.2.37.13069 | 1.2.37.13003 | GUI.zip and Web service updates |
| 37 | 3 | 1.2.37.13069 | 1.2.37.13040 | 1.2.37.13069 | 1.2.37.13003 | WebAPI updated, no version changes |
| 37 | 4 | 1.2.37.13171 | 1.2.37.13040 | 1.2.37.13165 | 1.2.37.13171 | GUI.zip, web service, and database updates |
| 37 | 5 | 1.2.37.13347 | 1.2.37.13040 | 1.2.37.12979 | 1.2.37.12979 | Bug was discovered with version numbers, this is why it appears we went backwards in the version number but we did not. |
| 37 | 6 | 1.2.37.13665 | 1.2.37.13040 | 1.2.37.13580 | 1.2.37.12979 | GUI.zip and Web service updates |
| 37 | 7 | 1.2.37.13807 | 1.2.37.13040 | 1.2.37.13799 | 1.2.37.13807 | GUI.zip and Web service updates. Also MIRTH mapping changes |
| 37 | 8 | 1.2.37.13807 | 1.2.37.13040 | 1.2.37.13827 | 1.2.37.13807 | Web service update |
| 37 | 9 | 1.2.37.13807 | 1.2.37.13040 | 1.2.37.13839 | 1.2.37.13807 | Web service update |
| 37 | 10 | 1.2.37.13895 | 1.2.37.13040 | 1.2.37.13895 | 1.2.37.13895 | GUI.zip, Web service and database updates |
| 37 | 11 | 1.2.37.14161 | 1.2.37.13040 | 1.2.37.14161 | 1.2.37.14161 | GUI.zip, Web service, database updates and SQL Management reporting |
| 37 | 12 | 1.2.37.14267 | 1.2.37.13040 | 1.2.37.14161 | 1.2.37.14267 | GUI.zip, Web service and database updates |
| 37 | 13 | 1.2.37.14326 | 1.2.37.13040 | 1.2.37.14326 | 1.2.37.14326 | GUI.zip, Web service and database updates |
| 38 | - | 1.2.38.14486 | 1.2.38.14486 | 1.2.38.14486 | 1.2.38.14486 | Full version release |
| 38 | 1 | 1.2.38.14748 | 1.2.38.14748 | 1.2.38.14748 | 1.2.38.14748 | Core install required |
| 38 | 2 | 1.2.38.15182 | 1.2.38.14748 | 1.2.38.14748 | 1.2.38.14748 | GUI.zip update only |
| 38 | 3 | 1.2.38.15251 | 1.2.38.14748 | 1.2.38.14748 | 1.2.38.14748 | GUI.zip update only |
| 38 | 4 | 1.2.38.15327 | 1.2.38.14748 | 1.2.38.15290 | 1.2.38.14748 | GUI.zip and Web service updates |
| 38 | 5 | 1.2.38.15513 | 1.2.38.14748 | 1.2.38.15513 | 1.2.38.15451 | GUI.zip, Web service and database updates |
| 38 | 6 | 1.2.38.15618 | 1.2.38.14748 | 1.2.38.15513 | 1.2.38.15618 | GUI.zip and database updates |
| 38 | 7 | 1.2.38.15768 | 1.2.38.14748 | 1.2.38.15768 | 1.2.38.15618 | GUI.zip and Web Service updates |
| 39 | - | 1.3.39.15389 | 1.3.39.15389 | 1.3.39.15377 | 1.3.39.15355 | Full version release |
| 39 | 1 | 1.3.39.16026 | 1.3.39.15943 | 1.3.39.15943 | 1.3.39.15957 | Requires Core Install. GUI.zip, Web Service, and Database updates |
| 40 | - | 1.3.40.16628 | 1.3.40.16628 | 1.3.40.16628 | 1.3.40.16628 | Full version release |
| 40 | 1 | 1.3.40.16849 | 1.3.40.16659 | 1.3.40.16744 | 1.3.40.16727 | GUI.zip, Web service, database updates and SQL Management reporting |
| 40 | 2 | 1.3.40.16868 | 1.3.40.16659 | 1.3.40.16868 | 1.3.40.16727 | GUI.zip and Web Service |

RADNET, Inc. – eRAD RIS Server Update

| | | | | | | |
|----|---|--------------|--------------|--------------|--------------|--|
| 40 | 3 | 1.3.40.16868 | 1.3.40.16659 | 1.3.40.16868 | 1.3.40.16727 | DB updated only |
| 40 | 4 | 1.3.40.16967 | 1.3.40.16659 | 1.3.40.16932 | 1.3.40.16961 | GUI.zip, Web service, database updates |
| 40 | 5 | 1.3.40.17087 | 1.3.40.16659 | 1.3.40.17069 | 1.3.40.17087 | GUI.zip, Web service, database updates |
| 40 | 6 | 1.3.40.17146 | 1.3.40.16659 | 1.3.40.17146 | 1.3.40.17146 | GUI.zip, Web service, database updates |
| 40 | 7 | 1.3.40.17268 | 1.3.40.16659 | 1.3.40.17268 | 1.3.40.17259 | GUI.zip, Web service, database updates |
| 40 | 8 | 1.3.40.17364 | 1.3.40.16659 | 1.3.40.17362 | 1.3.40.17343 | GUI.zip, Web service, database updates |

6. Installing

- 1) Copy and replace the current rRIS...zip file with the rRIS_1.3.40.17364.zip file provided with this release.
Be sure to make a backup of the ris.exe.config file contained within the current .zip folder.

- 2) Stop each service listed below. Copy the files from the RIS Service Updates folder to the directory for each rRISService. Be sure to make backup of all files, but especially the rRISServices.exe.config and rRISServices.vshost.exe.config files.
rRISService directories include:

- a) *rRISService*
- b) *CDS rRISService*
- c) *Wedge rRISService*
- d) *DDE rRISService*

Restart each service.

- 3) Run the B40_008.sql file against the eRAD RIS database instance.