

Server Update

**For eRAD RIS
Version 1.0
Build 38
Update 38.003**

Table of Contents

1. Purpose.....	3
Defects resolved.....	3
2. Package Contents	4
3. Intended Audience	5
4. Who is affected	5
5. eRAD RIS Release Version Numbers	5
6. Installing	6

RADNET, Inc. – eRAD RIS Server Update

1. Purpose

This server update has been produced to fix bugs required for production setting, or to complete an enhancement requested for production.

This server update includes the following:

Defects resolved

Redmine #	Subject	Priority
2649	Reporting - Field's type attribute datatype cast exception. Translation – Inserting old, pre-existing macros will cause error.	High

2. Package Contents

Files included for this server update:

1. GUI ZIP Updates folder contains
 - a. AudioLib.dll
 - b. AudioLib.pdb
 - c. CDSJobManagerConsoleClient.exe
 - d. CDSJobManagerConsoleClient.pdb
 - e. CommonLib.dll
 - f. CommonLib.pdb
 - g. Diagram Paint.pdb
 - h. DynamicFormEditor.dll
 - i. DynamicFormEditor.pdb
 - j. MessageHelper.pdb
 - k. Microsoft.Practices.EnterpriseLibrary.Caching.pdb
 - l. Microsoft.Practices.EnterpriseLibrary.Common.pdb
 - m. Microsoft.Practices.EnterpriseLibrary.Logging.pdb
 - n. Microsoft.Practices.ServiceLocation.pdb
 - o. Microsoft.Practices.Unity.Interception.pdb
 - p. Microsoft.Practices.Unity.pdb
 - q. ris.exe
 - r. ris.pdb
 - s. rupdate.pdb
 - t. UILib.dll
 - u. UILib.pdb
 - v. WCFClient.exe
 - w. WCFClient.pdb
 - x. WSNewCropClient.exe
 - y. WSNewCropClient.pdb

3. Intended Audience

This document is created by the RIS Development team for the RadNet RIS Implementation team.

4. Who is affected

Previous versions of build 38.

Server updates may not be cumulative; therefore this server update must be applied directly to the previous server update.

5. eRAD RIS Release Version Numbers

Build	Patch	UI Version	Core Version	WS Version	DB Version	Notes
37	-	1.2.37.12843	1.2.37.12782	1.2.37.12843	1.2.37.12843	
37	1	1.2.37.13040	1.2.37.13040	1.2.37.13039	1.2.37.13003	
37	2	1.2.37.13069	1.2.37.13040	1.2.37.13069	1.2.37.13003	
37	3	1.2.37.13069	1.2.37.13040	1.2.37.13069	1.2.37.13003	Just WebAPI updated, no version changes
37	4	1.2.37.13171	1.2.37.13040	1.2.37.13165	1.2.37.13171	
						Bug was discovered with version numbers, this is why it appears we went backwards in the version number but we really did not.
37	5	1.2.37.13347	1.2.37.13040	1.2.37.12979	1.2.37.12979	
37	6	1.2.37.13665	1.2.37.13040	1.2.37.13580	1.2.37.12979	
37	7	1.2.37.13807	1.2.37.13040	1.2.37.13799	1.2.37.13807	
37	8	1.2.37.13807	1.2.37.13040	1.2.37.13827	1.2.37.13807	
37	9	1.2.37.13807	1.2.37.13040	1.2.37.13839	1.2.37.13807	
37	10	1.2.37.13895	1.2.37.13040	1.2.37.13895	1.2.37.13895	
37	11	1.2.37.14161	1.2.37.13040	1.2.37.14161	1.2.37.14161	
37	12	1.2.37.14285	1.2.37.13040	1.2.37.14300	1.2.37.14285	
37	13	1.2.37.14326	1.2.37.13040	1.2.37.14326	1.2.37.14326	
38	-	1.2.38.14486	1.2.38.14486	1.2.38.14486	1.2.38.14486	
38	1	1.2.38.14748	1.2.38.14748	1.2.38.14748	1.2.38.14748	
38	2	1.2.38.15182	1.2.38.14748	1.2.38.14748	1.2.38.14748	
38	3	1.2.38.15251	1.2.38.14748	1.2.38.14748	1.2.38.14748	

6. Installing

- 1) Copy the files in the GUI ZIP Updates folder and replace those in the rRIS_1.2.38.15182.zip folder. Rename the zip file to be **rRIS_1.2.38. 15251.zip**

NOTE: A core install is not required for this update. Simply exit the application and reopen to get the latest changes.