



eRAD RIS

RELEASE ANNOUNCEMENT

Build 4.2024.075

UPDATED JUNE 28, 2024

TABLE OF CONTENTS

Summary 2

 Release Announcement..... 2

 New Features..... 2

Feature Details 3

 Insurance Eligibility..... 3

 Feature #33872 - Addition of Unique Key Identifier within Warning Messages.. 3

Version Details 4

 Code Stream..... 4

PUBLICATION HISTORY

Revision	Author	Description
June 28, 2024	Kevin Brooks, Michelle Mahabir / Hilary Saltmarsh	<ul style="list-style-type: none"> Commercial release.

SUMMARY

Release Announcement

This release of ERAD RIS 4.2024.075 introduces an Insurance Eligibility enhancement that allows for efficient detection of triggered Eligibility Special Conditions warnings.

New Features

This release introduces the following features and enhancements:

Category	Redmine #	Subject	Description
Insurance Eligibility	33872	Addition of a Unique Key Identifier within Warning Messages	This enhancement to Insurance Eligibility provides the ability to identify which Eligibility Special Condition warning is being triggered and displayed to the end user.

SORTED BY CATEGORY AND REDMINE

Refer to the [FEATURE DETAILS](#) section below for configuration and usage information.

FEATURE DETAILS

Insurance Eligibility

Feature #33872 - Addition of Unique Key Identifier within Warning Messages

Summary

This enhancement to Insurance Eligibility provides the ability to identify which Eligibility Special Condition warning is being triggered and being displayed to the end user.


Background

Previously, the process of identifying the logic that was triggered involved analyzing multiple line items within the table to detect which logic line was causing the warning to be displayed.

Feature Description

To better identify the triggered logic line, a **Key** column has been added to the Eligibility Special Conditions table that will contain the unique key identifier from the warning.

This identifier will also appear within the HTML view of the source code of the Eligibility Details. To review this source code, **Right-Click** anywhere in the Eligibility details to open the context menu and select **View source** to open in Notepad. Then, press **CTRL + F** within Notepad and search the source code for `<!--` to find the unique key identifier.



Configuration Instructions

System Administrators must complete the following actions to enable this feature:

RIS Client

Changes to Eligibility Special Conditions Lookup Table

A **Key** column has been added to assist with identifying the unique key identifier within the warning message.

- Review the **Warning Message** text of all current (and future) items in the Eligibility Special Conditions table and ensure the warning includes the `<!--` prefix, e.g.:

```
<!-- ' + cast(eligibility_special_condition_key as varchar) + '-->
```



Caution - Ensure all items added in the future also include this prefix.

VERSION DETAILS

Code Stream

The following source code branches have been merged into this release:

