

Server Update

for eRAD RIS

Version 1.0

Build 37

Update 37.003

Table of Contents

1. Purpose.....	3
2. Intended Audience	3
3. Who is affected	3
4. Installing	4

1. Purpose

This server update has been produced to fix bugs required for production setting, or to complete an enhancement requested for production.

This server update includes the following:

1. Logged as Redmine bug 1676 - Users who are domain users cannot login when they do not have a resource type.

Files included are:

- a) WebApi.dll
- b) WebApi.pdb

2. Intended Audience

This document is created by the RIS Development team for the RadNet RIS Implementation team.

3. Who is affected

All previous versions of build 37.

4. Installing

- 1) Stop each service listed below. Copy the files from the RIS Service Updates folder to the directory for each rRISService.

rRISService directories include:

- a) rRISService*
- b) CDS rRISService*
- c) Wedge rRISService*
- d) DDE rRISService*

Restart each service