



## eRAD RIS

# RELEASE ANNOUNCEMENT

### Build 4.2023.033

UPDATED MAY 30, 2023

#### TABLE OF CONTENTS

**Summary** ..... 2

    Summary Release Announcement..... 2

    New Features..... 2

    Resolved Issues..... 2

**Feature Details** ..... 3

    Provider Portal Administration..... 3

        Feature #31661 – Add User Type Confirmation Step When Creating New User in Referring Admin..... 3

#### PUBLICATION HISTORY

Revision	Author	Description
May 30, 2023	Kevin Brooks, Michelle Mahabir / Hilary Saltmarsh	▪ Commercial release.

# SUMMARY

## Summary Release Announcement

*This release of eRAD RIS 4.2023.033 resolves an issue creating Provider Portal users, and adds the ability to select the desired portal user account type during account creation.*

## New Features

This release introduces the following features and enhancements:

Category	Redmine #	Subject	Description
Provider Portal	31661	Add User Type Confirmation Step When Creating New User in Referring Admin	This enhancement to Provider Portal User Management adds a confirmation step when creating a new user.

Refer to the Feature Details section below for configuration and usage information.

## Resolved Issues

This release resolves the following issues:

Category	Redmine #	Subject
Provider Portal	33750	Resolves a Provider Portal Administration issue where creating a new referring user that was previously affiliated with another referring user returns an error.

# FEATURE DETAILS

## Provider Portal Administration

### Feature #31661 – Add User Type Confirmation Step When Creating New User in Referring Admin


#### Summary

This enhancement to Provider Admin User Management adds a workflow step to confirm the user type when creating a new user within the Referring Admin Portal.

#### Feature Description


Previously, users were erroneously creating Providers as Staff Users.

With this change, a new dialog has been added after the Create New User button is selected to clarify their intent:



If **STAFF** is selected, the user is directed to the existing **Create New User** screen.

However, if **PROVIDER** is selected, the user is prompted to search for an existing provider account before a new account can be created:



#### Known Limitations

The following significant limitations have been identified and should be communicated to affected users:

- **BUG #33846**

Within the UM Admin Portal, attempting to create an account using the Add button returns an error.



*Carefully review all Known Limitations identified here and in the release notes of all other versions being applied and communicate relevant items to all those affected.*

## Configuration Instructions

No System Administrator actions are necessary to enable this feature.