



eRAD RIS

CUSTOMER RELEASE NOTES

Hotfix Build v3.2018.3.4

UPDATED APRIL 14, 2021




Table of Contents

Summary	2
Intended Audience	2
Who Is Affected	2
New Features	3
Scheduling	3
Feature #26293 - Ability to directly edit or remove a closure from the Appointment Book	3
Resolved Issues and Known Limitations	4
Resolved Issues	4
New Known Limitations	4
Version Details	5
Package Contents	5
Code Stream	5
eRAD RIS Release Version Numbers	6

Publication History

Revision	Author	Description
April 14, 2021	Kevin Brooks / Hilary Saltmarsh	<ul style="list-style-type: none"> ▪ Documentation for commercial release.

SUMMARY

Intended Audience

The intended audience for this CUSTOMER RELEASE NOTES document is the RIS Administration team for all eRAD RIS customers.

This document describes the purpose, configuration, and operation of new features made available with this release, identifies issues resolved in the release, and highlights any unresolved known limitations.

This information should be used by all customers to determine how these changes affect their organization's workflow, and to plan their organization's upgrade strategy for eRAD RIS.



There is no separate SERVICE RELEASE NOTES edition of this document.

Who Is Affected

This hotfix release of eRAD RIS v3.2018.3.4 includes a single feature enhancement and one resolved issue.

This version contains new features and bug fixes that are recommended to be applied on eRAD RIS v3.2018.3.3 installations.

NEW FEATURES

Scheduling

Feature #26293 - Ability to directly edit or remove a closure from the Appointment Book

Summary

This enhancement to Scheduling adds the ability to edit or remove a closure or restriction directly from the Appointment Book.

Previously, in version 18.3.2 (#25391) we added the ability to create a closure (or restriction) directly from the Appointment Book.

With this change users can also remove a closure when no longer needed.

Users with Full access to `Config.AppointmentBook.ModalityClosure` will now see `Edit Closure` (or `Edit Restriction`) and `Remove Closure` (or `Remove Restriction`) context menus when right-clicking on an already closed or restricted timeslot in the Appointment Book.

RESOLVED ISSUES AND KNOWN LIMITATIONS

Resolved Issues

This release resolves the following issues:

Redmine #	Subject
26293	Resolved Scheduling issue where Appointment Search results were off by one minute.

New Known Limitations

No new known limitations were identified with this release.

VERSION DETAILS


Package Contents

The release package includes the following folders:

 @Hotfixes	4/20/2012 8:51 AM
 _Documentation	3/12/2020 2:28 PM
 _ReleaseNotes	3/26/2012 11:02 AM
 Client Application	3/12/2020 2:25 PM
 DB	3/12/2020 2:26 PM
 External WebAPI	3/12/2020 2:26 PM
 Identity Service	3/12/2020 2:26 PM
 Management Reports	3/12/2020 2:26 PM
 Mirth Config	3/12/2020 2:26 PM
 MWLServerConfig	3/12/2020 2:26 PM
 PACS Citrix Bridge	3/12/2020 2:26 PM
 RADAR Integration	3/11/2020 3:38 PM
 RIS Service	3/12/2020 2:26 PM
 Service Tools	3/12/2020 2:27 PM
 Web Digital Forms	3/12/2020 2:27 PM
 Web Patient Connect	3/12/2020 2:27 PM
 Web Referring Connect	3/12/2020 2:27 PM
 Web UM Connect	3/12/2020 2:28 PM
 WebHelp	3/12/2020 2:28 PM
 Build_2018.3.4_955.7z	3/11/2020 3:42 PM
 RISServerMasterCert.pfx	3/31/2016 1:38 PM

Code Stream

The following source code branches have been merged into this release:



LEGEND:

- Light Green = Previously Released software
- Gray = Internal version, non-release version
- Bright Green = Current Release

eRAD RIS Release Version Numbers

The following table details the version identifiers for components in this release:

Build	Patch	UI Version	Core Version	WS Version	DB Version	Digital Forms	Patient Portal	UM Portal	Provider Portal	Notes
2018.1	-	3.18.1.0(3GB)	3.18.1.0	3.18.1.0	3.18.1.0.01228009	3.18.1.0	3.18.1.0.1229289	3.18.1.0.1229290	3.18.1.0.1229290	Full version release
2018.1	1	3.18.1.1(3GB)	3.18.1.0	3.18.1.1	3.18.1.1.01288419	3.18.1.1	3.18.1.1.1299765	3.18.1.1.1299766	3.18.1.1.1299765	GUI, Web Service, DB, Digital Forms, Patient, Provider and UM portals
2018.1	2	3.18.1.2(3GB)	3.18.1.0	3.18.1.2	3.18.1.2.01341771	3.18.1.2	3.18.1.2.1355946	3.18.1.2.1355947	3.18.1.2.1355947	GUI, Web Service, DB, Digital Forms, Patient, Provider and UM portals
2018.1	2.1	3.18.1.2.1(3GB)	3.18.1.0	3.18.1.2.1	3.18.1.2.01341771	3.18.1.2	3.18.1.2.1355946	3.18.1.2.1355947	3.18.1.2.1355947	GUI and Web Service
2018.1	3	3.18.1.3(3GB)	3.18.1.0	3.18.1.3	3.18.3.0.01471392	3.18.1.3	3.18.1.3.1466761	3.18.1.3.1466762	3.18.1.3.1466762	GUI, Web Service, DB, Management Reports, Digital Forms, Patient, Provider and UM portals
2018.1	4	3.18.1.4(3GB)	3.18.1.0	3.18.1.4	3.18.1.4.01518863	3.18.1.4	3.18.1.4.1527377	3.18.1.4.1527379	3.18.1.4.1527378	GUI, Web Service, DB, Digital Forms, Patient, Provider and UM portals
2018.1	5	3.18.1.5(3GB)	3.18.1.0	3.18.1.5	3.18.1.5.01547447	3.18.1.5	3.18.1.5.1561651	3.18.1.5.1561652	3.18.1.5.1561653	GUI, Web Service, DB, Digital Forms, Patient, Provider and UM portals
2018.1	5.1	3.18.1.5.1(3GB)	3.18.1.0	3.18.1.5.1	3.18.1.5.1.01613980	3.18.1.5.1	3.18.1.5.1.743	3.18.1.5.1.743	3.18.1.5.1.743	GUI, Web Service, DB, Digital Forms, Patient, Provider and UM portals
2018.1	5.2	3.18.1.5.1(3GB)	3.18.1.0	3.18.1.5.2	3.18.1.5.1.01613980	3.18.1.5.1	3.18.1.5.1.743	3.18.1.5.1.743	3.18.1.5.1.743	Web Services only
2018.1	6	3.18.1.6(3GB)	3.18.1.0	3.18.1.6	3.18.1.6.0.01694594	3.18.1.6	3.18.1.6.0.776	3.18.1.6.0.776	3.18.1.6.0.776	GUI, Web Service, DB, Digital Forms, Patient, Provider and UM portals
2018.1	6.1	3.18.1.6(3GB)	3.18.1.0	3.18.1.6	3.18.1.6.1.01720366	3.18.1.6	3.18.1.6.0.776	3.18.1.6.0.776	3.18.1.6.0.776	DB updates
2018.2	-	3.18.2(3GB)	3.18.2.0	3.18.2.0	3.18.2.0.01412126	3.18.2.0	3.18.2.0.1416370	3.18.2.0.1416371	3.18.2.0.1416371	Full Version Release. Including Patient, Provider and UM Portals
2018.2	1	3.18.2.1(3GB)	3.18.2.0	3.18.2.1	3.18.2.1.01490049	3.18.2.1	3.18.2.1.1509822	3.18.2.1.1509823	3.18.2.1.1509822	GUI, Web Service, DB, Management Reports, Digital Forms, Patient, Provider and UM portals
2018.2	2	3.18.2.2(3GB)	3.18.2.0	3.18.2.2	3.18.2.2.01563299	3.18.2.2	3.18.2.2.1583286	3.18.2.2.1583287	3.18.2.2.1583286	GUI, Web Service, DB, Digital Forms, Patient, Provider and UM portals
2018.2	3	3.18.2.3(3GB)	3.18.2.0	3.18.2.3	3.18.2.2.01563299	3.18.2.3	3.18.2.3.1592002	3.18.2.3.1592003	3.18.2.3.1592003	GUI, Web Service, Digital Forms, Patient, Provider and UM portals
2018.2	4	3.18.2.4(3GB)	3.18.2.0	3.18.2.4	3.18.2.4.01610643	3.18.2.4	3.18.2.3.1592002	3.18.2.3.1592003	3.18.2.3.1592003	GUI, Web Service, DB and Digital Forms
2018.3	-	3.18.3(3GB)	3.18.3.0	3.18.3.0	3.18.3.0.01547822	3.18.3	N/A	N/A	N/A	Full Version Release. Including Patient, Provider and UM Portals
2018.3	1	3.18.3.1(3GB)	3.18.3.1	3.18.3.1	3.18.3.1.01693170	3.18.3.1	N/A	N/A	N/A	Full Version Release.
2018.3	2	3.18.3.2.0(3GB)	3.18.3.2	3.18.30200.0	3.18.3.2.0.01905964	3.18.3.2	N/A	N/A	N/A	GUI, Web Service, DB updates
2018.3	3	3.18.3.3.0(3GB)	3.18.3.3	3.18.30300.0	3.18.3.3.0.02152589	3.18.3.3	N/A	N/A	N/A	GUI, Web Service, DB updates
2018.3	4	3.18.3.4.0(3GB)	3.18.3.4	3.18.30400.0						GUI and Web Service