



eRAD RIS

# **CUSTOMER & SERVICE RELEASE NOTES**

Hotfix Build v3.2022.1.17.2

UPDATED FEBRUARY 18, 2022



# Table of Contents

|  |          |
|--|----------|
| <b>Summary</b> .....                               | <b>2</b> |
| Intended Audience .....                            | 2        |
| Who Is Affected .....                              | 2        |
| <b>New Settings</b> .....                          | <b>3</b> |
| <b>New Features</b> .....                          | <b>4</b> |
| <b>Resolved Issues and Known Limitations</b> ..... | <b>5</b> |
| Resolved Issues .....                              | 5        |
| New Known Limitations .....                        | 5        |
| <b>Version Details</b> .....                       | <b>6</b> |
| Package Contents .....                             | 6        |
| Code Stream .....                                  | 7        |
| eRAD RIS Release Version Numbers .....             | 8        |

## Publication History

| Revision          | Author                             | Description   |
|-------------------|------------------------------------|---|
| February 18, 2022 | Kevin Brooks<br>/ Hilary Saltmarsh | <ul style="list-style-type: none"> <li>▪ Commercial release.</li> </ul> |

# SUMMARY

## Intended Audience

The intended audience for this CUSTOMER RELEASE NOTES document is the RIS Administration team for all eRAD RIS customers.

This document describes the purpose, configuration, and operation of new features made available with this release, identifies issues resolved in the release, and highlights any unresolved known limitations.

This information should be used by all customers to determine how these changes affect their organization's workflow, and to plan their organization's upgrade strategy for eRAD RIS.



*There is no separate SERVICE RELEASE NOTES edition of this document.*

## Who Is Affected

This hotfix release of eRAD RIS v3.2022.1.17.2 resolves an urgent Registration memory handling issue.

This version is recommended to be applied on eRAD RIS v3.2022.1.17 or eRAD RIS v3.2022.1.17.1 installations.



*Please carefully review these release notes even if your system will not be upgraded immediately, to identify and communicate any issues that may affect your organization.*

# NEW SETTINGS

No settings were added or updated with this release.

# NEW FEATURES

No new features were introduced with this release.

# RESOLVED ISSUES AND KNOWN LIMITATIONS

## Resolved Issues

This release resolves the following issues:

| Redmine # | Subject  |
|-----------|--|
| 31488     | Resolved memory leak issue in the Registration screen. |

## New Known Limitations

No new known limitations were identified with this release.

# VERSION DETAILS


## Package Contents

The release package includes the following folders:

|   |                   |
|---|-------------------|
|  _ReleaseNotes           | 2/18/2022 8:54 AM |
|  Client Application      | 2/15/2022 5:03 PM |
|  Build_2022.1.17.2.zip   | 2/15/2022 5:07 PM |
|  RISServerMasterCert.pfx | 3/31/2016 1:38 PM |

## Code Stream

The following source code branches have been merged into this release:



LEGEND:

- Light Green = Previously Released software
- Gray = Internal version, non-release version
- Bright Green = Current Release

## eRAD RIS Release Version Numbers

The following table details the version identifiers for components in this release:

| Build      | Patch | UI Version        | Core Version | WS Version   | DB Version           | Digital Forms | Patient Portal    | UM Portal         | Provider Portal   | Notes  |
|------------|-------|-------------------|--------------|--------------|----------------------|---------------|-------------------|-------------------|-------------------|--|
| 2021.2.15  | -     | 3.21.2.15.0(3GB)  | 3.21.2.15    | 3.21.2.15.0  | 3.21.2.15.002698266  | 3.21.2.15.0   | 3.21.2.15.0.1111  | 3.21.2.15.0.1111  | 3.21.2.15.0.1111  | GUI, Web Services, DB, Patient/Provider/UM portals and Digital Forms |
| 2021.2.15  | 1     | 3.21.2.15.1(3GB)  | 3.21.2.15.1  | 3.21.2.15.1  |                      | 3.21.2.15.1   | 3.21.2.15.1.1113  | 3.21.2.15.1.1113  | 3.21.2.15.1.1113  | GUI, Web Services, Patient/Provider/UM portals and Digital Forms     |
| 2021.2.15  | 2     | 3.21.2.15.2(3GB)  | 3.21.2.15.2  | 3.21.2.15.2  |                      | 3.21.2.15.2   | 3.21.2.15.2.1122  | 3.21.2.15.2.1122  | 3.21.2.15.2.1122  | GUI, Web Services, Patient/Provider/UM portals and Digital Forms     |
| 2021.2.15  | 3     | 3.21.2.15.3(3GB)  | 3.21.2.15.3  |              |                      |               |                   |                   |                   | GUI  |
| 2021.2.15  | 4     | 3.21.2.15.4(3GB)  | 3.21.2.15.4  | 3.21.2.15.4  | 3.21.2.15.402732838  | 3.21.2.15.4   | 3.21.2.15.4.1133  | 3.21.2.15.4.1133  | 3.21.2.15.4.1133  | GUI, Web Services, DB, Patient/Provider/UM portals and Digital Forms |
| 2021.2.15  | 5     | 3.21.2.15.5(3GB)  | 3.21.2.15.5  | 3.21.2.15.5  | 3.21.2.15.502738557  |               |                   |                   |                   | GUI, Web Services and DB   |
| 2021.2.15  | 6     | 3.21.2.15.6(3GB)  | 3.21.2.15.6  | 3.21.2.15.6  |                      |               |                   |                   |                   | GUI and Web Services   |
| 2021.3.1   | -     | 3.21.3.1(3GB)     | 3.21.3.1     | 3.21.3.1     | 3.21.3.1.002712308   | 3.21.3.1      | 3.21.3.1.0.1119   | 3.21.3.1.0.1119   | 3.21.3.1.0.1119   | GUI, Web Services, DB, Patient/Provider/UM portals and Digital Forms |
| 2021.3.15  | -     | 3.21.3.15(3GB)    | 3.21.3.15    | 3.21.3.15    | 3.21.3.15.002747230  | 3.21.3.15     | 3.21.3.15.0.1144  | 3.21.3.15.0.1144  | 3.21.3.15.0.1144  | GUI, Web Services, DB, Patient/Provider/UM portals and Digital Forms |
| 2021.4.12  | -     | 3.21.4.12(3GB)    | 3.21.4.12    | 3.21.4.12    | 3.21.4.12.002778929  | 3.21.4.12     | 3.21.4.12.0.1167  | 3.21.4.12.0.1167  | 3.21.4.12.0.1167  | GUI, Web Services, DB, Patient/Provider/UM portals and Digital Forms |
| 2021.5.10  | -     | 3.21.5.10(3GB)    | 3.21.5.10    | 3.21.5.10    | 3.21.5.10.002823110  | 3.21.5.10     | 3.21.5.10.0.1189  | 3.21.5.10.0.1189  | 3.21.5.10.0.1189  | GUI, Web Services, DB, Patient/Provider/UM portals and Digital Forms |
| 2021.5.24  | -     | 3.21.5.24(3GB)    | 3.21.5.24    | 3.21.5.24    | 3.21.5.24.002842120  | 3.21.5.24     | 3.21.5.24.0.1205  | 3.21.5.24.0.1205  | 3.21.5.24.0.1205  | GUI, Web Services, DB, Patient/Provider/UM portals and Digital Forms |
| 2021.5.24  | 1     | 3.21.5.24.1(3GB)  | 3.21.5.24.1  | 3.21.5.24.1  | 3.21.5.24.102880695  | 3.21.5.24.1   | 3.21.5.24.1.1224  | 3.21.5.24.1.1224  | 3.21.5.24.1.1224  | GUI, Web Services, DB, Patient/Provider/UM portals and Digital Forms |
| 2021.5.24  | 2     | 3.21.5.24.2(3GB)  | 3.21.5.24.2  | 3.21.5.24.2  | 3.21.5.24.20293248   | 3.21.5.24.2   | 3.21.5.24.0.1244  | 3.21.5.24.0.1244  | 3.21.5.24.0.1244  | GUI, Web Services, DB, Patient/Provider/UM portals and Digital Forms |
| 2021.5.24  | 3     | 3.21.5.24.3(3GB)  | 3.21.5.24.3  | 3.21.5.24.3  | 3.21.5.24.302932869  | 3.21.5.24.3   | 3.21.5.24.0.1252  | 3.21.5.24.0.1252  | 3.21.5.24.0.1252  | GUI, Web Services, DB, Patient/Provider/UM portals and Digital Forms |
| 2021.5.24  | 4     | 3.21.5.24.4(3GB)  | 3.21.5.24.4  | 3.21.5.24.4  | 3.21.5.24.402961496  | 3.21.5.24.4   | 3.21.5.24.0.1260  | 3.21.5.24.0.1260  | 3.21.5.24.0.1260  | GUI, Web Services, DB, Patient/Provider/UM portals and Digital Forms |
| 2021.6.7   | -     | 3.21.6.7(3GB)     | 3.21.6.7     | 3.21.6.7     | 3.21.6.7.002898418   | 3.21.6.7      | 3.21.6.7.0.1234   | 3.21.6.7.0.1234   | 3.21.6.7.0.1234   | GUI, Web Services, DB, Patient/Provider/UM portals and Digital Forms |
| 2021.6.7   | 1     | 3.21.6.7.1(3GB)   | 3.21.6.7.1   | 3.21.6.7.1   | 3.21.6.7.102961651   | 3.21.6.7.1    | 3.21.6.7.1.1262   | 3.21.6.7.1.1262   | 3.21.6.7.1.1262   | GUI, Web Services, DB, Patient/Provider/UM portals and Digital Forms |
| 2021.7.5   | -     | 3.21.7.5(3GB)     | 3.21.7.5     | 3.21.7.5     | 3.21.07.5.002970391  | 3.21.7.5      | 3.21.7.5.0.1266   | 3.21.7.5.0.1266   | 3.21.7.5.0.1266   | GUI, Web Services, DB, Patient/Provider/UM portals and Digital Forms |
| 2021.9.27  | -     | 3.21.9.27(3GB)    | 3.21.9.27    | 3.21.9.27    | 3.21.9.27.003025038  | 3.21.9.27     | 3.21.9.27.0.1280  | 3.21.9.27.0.1280  | 3.21.9.27.0.1280  | GUI, Web Services, DB, Patient/Provider/UM portals and Digital Forms |
| 2021.10.25 | -     | 3.21.10.25(3GB)   | 3.21.10.25   | 3.21.10.25   | 3.21.10.25.003071251 | 3.21.10.25    | 3.21.10.25.0.1317 | 3.21.10.25.0.1317 | 3.21.10.25.0.1317 | GUI, Web Services, DB, Patient/Provider/UM portals and Digital Forms |
| 2021.10.25 | 1     | 3.21.10.25.1(3GB) |              | 3.21.10.25.1 |                      | 3.21.10.25.1  |                   |                   |                   | Portals  |
| 2021.11.22 | -     | 3.21.11.22(3GB)   | 3.21.11.22   | 3.21.11.22   | 3.21.11.22.003102803 | 3.21.11.22    | 3.21.11.22.0.1346 | 3.21.11.22.0.1346 | 3.21.11.22.0.1346 | GUI, Web Services, DB, Patient/Provider/UM portals and Digital Forms |
| 2021.12.6  | -     | 3.21.12.6(3GB)    | 3.21.12.6    | 3.21.12.6    | 3.21.12.6.003112661  | 3.21.12.6     | 3.21.12.6.0.1352  | 3.21.12.6.0.1352  | 3.21.12.6.0.1352  | GUI, Web Services, DB, Patient/Provider/UM portals and Digital Forms |
| 2021.12.20 | -     | 3.21.12.20(3GB)   | 3.21.12.20   | 3.21.12.20   | 3.21.12.20.003133170 | 3.21.12.20    | 3.21.12.20.0.1360 | 3.21.12.20.0.1360 | 3.21.12.20.0.1360 | GUI, Web Services, DB, Patient/Provider/UM portals and Digital Forms |
| 2021.12.20 | 1     | 3.21.12.20.1(3GB) | 3.21.12.20.1 | 3.21.12.20.1 | 3.21.12.20.103170667 | 3.21.12.20.1  | 3.21.12.20.1.1366 | 3.21.12.20.1.1366 | 3.21.12.20.1.1366 | GUI, Web Services, DB, Patient/Provider/UM portals and Digital Forms |
| 2021.12.20 | 2     | 3.21.12.20.2(3GB) | 3.21.12.20.2 |              |                      |               |                   |                   |                   | GUI, Web Services, DB, Patient/Provider/UM portals and Digital Forms |
| 2022.1.17  | -     | 3.22.1.17(3GB)    | 3.22.1.17    | 3.22.1.17    | 3.22.1.17.003185028  | 3.22.1.17     | 3.22.1.17.0.1382  | 3.22.1.17.0.1382  | 3.22.1.17.0.1382  | GUI, Web Services, DB, Patient/Provider/UM portals and Digital Forms |
| 2022.1.17  | 1     | 3.22.1.17.1(3GB)  | 3.22.1.17.1  | 3.22.1.17.1  |                      | 3.22.1.17.1   |                   |                   |                   | GUI and Web Services   |
| 2022.1.17  | 2     | 3.22.1.17.1(3GB)  | 3.22.1.17.1  |              |                      |               |                   |                   |                   | GUI only   |