



eRAD RIS

CUSTOMER RELEASE NOTES

Hotfix Build 3.2022.8.29.2

UPDATED NOVEMBER 16, 2022

Table of Contents

- Summary 3
 - Intended Audience..... 3
 - Release Content..... 3
 - Who Is Affected..... 3
- New Settings 4
- New Features 5
- Resolved Issues and Known Limitations 6
 - Resolved Issues..... 6
 - New Known Limitations..... 6
- Version Details 7
 - Package Contents..... 7
 - Code Stream..... 8
 - eRAD RIS Release Version Numbers..... 9

Publication History

Revision	Author	Description
November 16, 2022	Kevin Brooks / Hilary Saltmarsh	▪ Commercial release.

SUMMARY

Intended Audience

The intended audience for this CUSTOMER RELEASE NOTES document is the RIS Administration team for all eRAD RIS customers.

This document describes the purpose, configuration, and operation of new features made available with this release, identifies issues resolved in the release, and highlights any unresolved known limitations.

This information should be used by all customers to determine how these changes affect their organization's workflow, and to plan their organization's upgrade strategy for eRAD RIS.



There is no separate SERVICE RELEASE NOTES edition of this document.

Release Content

This hotfix release of eRAD RIS 3.2022.8.29.2 resolves an urgent Radiology Reporting issue, as well as scanning save issue and worklist performance issues.



Please carefully review these release notes even if your system will not be upgraded immediately, to identify and communicate any issues that may affect your organization.

Who Is Affected

This version is recommended to be applied on eRAD RIS v3.2022.8.29.x installations.



Note that the functionality of this release has not yet been merged forward into any other release.
Please carefully review the code stream diagram of any future releases before deploying to ensure that the functionality of this release will not be lost when upgrading.

NEW SETTINGS

No settings were added or updated with this release.

NEW FEATURES

No new features were introduced with this release.

RESOLVED ISSUES AND KNOWN LIMITATIONS

Resolved Issues

This release resolves the following issues:

Redmine #	Subject
33022	Resolved Scanning issue where attachments failed to save.
33024	Resolved Database performance issue when updating worklist flags. Note the related KNOWN LIMITATION - BUG #33134.
33117	Resolved Radiology Reporting issue preventing signing of studies.

New Known Limitations

The following new Known Limitations were identified with this release:

- **BUG #33134- NOTES FLAG MAY NOT BE REMOVED AS EXPECTED**
 - **Issue:** Removing patient/visit notes from a patient might not turn off the associated flag.
 - **Impact:** The note flag will incorrectly remain on worklists.
 - **Workaround:** No known workaround. The issue is related to whitespace characters being saved while RIS Services only checks for null vs not null.



Please carefully review these Known Limitations even if your system will not be upgraded immediately, to identify and communicate any issues that may affect your organization.

VERSION DETAILS


Package Contents

The release package includes the following folders:

 _Documentation	11/9/2022 7:09 PM
 _ReleaseNotes	11/16/2022 11:01 AM
 Client Application	11/9/2022 7:09 PM
 DB	11/9/2022 7:09 PM
 External WebAPI	11/9/2022 7:09 PM
 Identity Service	11/9/2022 7:09 PM
 Management Reports	11/9/2022 7:09 PM
 RIS Service	11/9/2022 7:09 PM
 Service Tools	11/9/2022 7:09 PM
 Web Digital Forms	11/9/2022 7:09 PM
 Web Patient Connect	11/9/2022 7:09 PM
 Web Referring Connect	11/9/2022 7:09 PM
 Web UM Connect	11/9/2022 7:09 PM
 Build_2022.8.29.2.zip	11/9/2022 7:10 PM
 RISServerMasterCert.pfx	3/31/2016 1:38 PM

Code Stream

The following source code branches have been merged into this release:



LEGEND:

Light Green = Previously Released software

Gray = Internal version, non-release version

Bright Green = Current Release



Note that the functionality of this release has not yet been merged forward into any other release.

Please carefully review the code stream diagram of any future releases before deploying to ensure that the functionality of this release will not be lost when upgrading.

eRAD RIS Release Version Numbers

The following table details the version identifiers for components in this release:

Build	Patch	UI Version	Core Version	WS Version	DB Version	Digital Forms	Patient Portal	UM Portal	Provider Portal	Notes
2022.1.17	-	3.22.1.17(3GB)	3.22.1.17	3.22.1.17	3.22.1.17.003185028	3.22.1.17	3.22.1.17.0.1382	3.22.1.17.0.1382	3.22.1.17.0.1382	GUI, Web Services, DB, Portals and Digital Forms
2022.1.17	1	3.22.1.17.1(3GB)	3.22.1.17.1	3.22.1.17.1		3.22.1.17.1				GUI and Web Services
2022.1.17	2	3.22.1.17.2(3GB)	3.22.1.17.2							GUI only
2022.1.17	3	3.22.1.17.3(3GB)	3.22.1.17.3	3.22.1.17.3						GUI and Web Services only
2022.1.17	4	3.22.1.17.4(3GB)	3.22.1.17.4							GUI only
2022.1.17	4.1	3.22.1.17.4.1(3GB)	3.22.1.17.4.1							GUI only
2022.1.17	4.2	3.22.1.17.4.2(3GB)	3.22.1.17.4.2	3.22.1.17.4.2	3.22.1.17.4.203367866					GUI only, Web Services and DB
2022.1.7	4.3						3.22.1.17.4.3	3.22.1.17.4.3	3.22.1.17.4.3	Portals only, first rapid release
2022.1.7	4.4	3.22.1.17.4.4(3GB)	3.22.1.17.4.4	3.22.1.17.4.4			3.22.1.17.4.4	3.22.1.17.4.4	3.22.1.17.4.4	GUI, Web Services and Portals
2022.1.7	4.5	3.22.1.17.4.5(3GB)	3.22.1.17.4.5							GUI
2022.1.17	5	3.22.1.17.5(3GB)	3.22.1.17.5	3.22.1.17.5						GUI and Web Service
2022.2.14	-	3.22.2.14(3GB)	3.22.2.14	3.22.2.14	3.22.2.14.003205179	3.22.2.14	3.22.2.14.0.1394	3.22.2.14.0.1394	3.22.2.14.0.1394	GUI, Web Services, DB, Portals and Digital Forms
2022.2.28	-	3.22.2.28(3GB)	3.22.2.28	3.22.2.28	3.22.2.28.003233569	3.22.2.28	3.22.2.28.0.1403	3.22.2.28.0.1403	3.22.2.28.0.1403	GUI, Web Services, DB, Portals and Digital Forms
2022.3.14	-	3.22.3.14(3GB)	3.22.3.14	3.22.3.14	3.22.3.14.003266996	3.22.3.14	3.22.3.14.0.1410	3.22.3.14.0.1410	3.22.3.14.0.1410	GUI, Web Services, DB, Portals and Digital Forms
2022.3.28	-	3.22.3.28(3GB)	3.22.3.28	3.22.3.28	3.22.3.28.003295497	3.22.3.28	3.22.3.28.0.1428	3.22.3.28.0.1428	3.22.3.28.0.1428	GUI, Web Services, DB, Portals and Digital Forms
2022.4.11	-	3.22.4.11(3GB)	3.22.4.11	3.22.4.11	3.22.4.11.003326062	3.22.4.11	3.22.4.11.0.1436	3.22.4.11.0.1436	3.22.4.11.0.1436	GUI, Web Services, DB, Portals and Digital Forms
2022.4.25	-	3.22.4.25(3GB)	3.22.4.25	3.22.4.25	3.22.4.25.003341901	3.22.4.25	3.22.4.25.0.1441	3.22.4.25.0.1441	3.22.4.25.0.1441	GUI, Web Services, DB, Portals and Digital Forms
2022.5.9	-	3.22.5.9(3GB)	3.22.5.9	3.22.5.9	3.22.5.9.003343130	3.22.5.9	3.22.5.9.0.1450	3.22.5.9.0.1450	3.22.5.9.0.1450	GUI, Web Services, DB, Portals and Digital Forms
2022.5.9	1	3.22.5.9.1(3GB)	3.22.5.9.1							GUI Only
2022.5.9	2	3.22.5.9.2(3GB)	3.22.5.9.2	3.22.5.9.2			3.22.5.9.2	3.22.5.9.2	3.22.5.9.2	GUI Web Services and Portals
2022.6.6	-	3.22.6.6(3GB)	3.22.6.6	3.22.6.6	3.22.6.6.003381845	3.22.6.6	3.22.6.6.0.1469	3.22.6.6.0.1469	3.22.6.6.0.1469	GUI, Web Services, DB, Portals and Digital Forms
2022.6.20	-	3.22.6.20(3GB)	3.22.6.20	3.22.6.20	3.22.6.20.003432703	3.22.6.20	3.22.6.20.0	3.22.6.20.0	3.22.6.20.0	GUI, Web Services, DB, Portals and Digital Forms
2022.8.1	-	3.22.8.1(3GB)	3.22.8.1	3.22.8.1	3.22.8.1.003455261	3.22.8.1	3.22.8.1	3.22.8.1	3.22.8.1	GUI, Web Services, DB, Portals and Digital Forms
2022.8.1	1	3.22.8.1.1(3GB)	3.22.8.1.1	3.22.8.1.1						GUI, Web Services,
2022.8.1	2	3.22.8.1.2(3GB)	3.22.8.1.2	3.22.8.1.2	3.22.8.1.203553632					GUI, Web Services and DB
2022.8.29	-	3.22.8.29(3GB)	3.22.8.29	3.22.8.29	3.22.8.29.003564808	3.22.8.29.0	3.22.8.29.0	3.22.8.29.0	3.22.8.29.0	GUI, Web Services, DB, Portals and Digital Forms
2022.8.29	1	3.22.8.29.1(3GB)	3.22.8.29.1	3.22.8.29.1						GUI, Web Services,
2022.8.29	2	3.22.8.29.2(3GB)	3.22.8.29.2	3.22.8.29.2	3.22.8.29.203607010	3.22.8.29.2	3.22.8.29.0	3.22.8.29.0	3.22.8.29.0	GUI, Web Services, DB, Portals and Digital Forms