



eRAD RIS

RELEASE ANNOUNCEMENT

Build 4.2023.023

UPDATED APRIL 27, 2023

PUBLICATION HISTORY

Revision	Author	Description
April 27, 2023	Kevin Brooks, Michelle Mahabir / Hilary Saltmarsh	▪ Commercial release.

SUMMARY

Release Announcement

This release of ERAD RIS 4.2023.023 resolves a Lay Letter issue that was causing RADAR SMS messages to be sent without a site/practice code.

New Features

No new features were introduced with this release.

Resolved Issues

This release resolves the following issues:

SORTED BY CATEGORY AND REDMINE

Category	Redmine #	Subject
Mammography	33886	Resolved Lay Letter issue that was causing RADAR SMS messages to be sent without a site/practice code.