



eRAD RIS

RELEASE ANNOUNCEMENT

Build 4.2023.049

UPDATED JULY 5, 2023

TABLE OF CONTENTS

Summary 2

 Release Announcement 2

 Resolved Issues 2

Version Details 3

 Code Stream 3

PUBLICATION HISTORY

Revision	Author	Description
July 5, 2023	Kevin Brooks, Michelle Mahabir / Hilary Saltmarsh	<ul style="list-style-type: none"> ▪ Commercial release.

SUMMARY

Release Announcement

This release of eRAD RIS 4.2023.049 resolves issues with attachments when using the Radiology Reporting screen.

Resolved Issues

This release resolves the following issues:

Category	Redmine #	Subject
Radiology Reporting	33903	Resolved Radiology Reporting issue where maximizing an attachment in the Dictation screen returned an error.
Radiology Reporting	33899	Resolved Radiology Reporting issue where modifying a scanned document type returned an error.
Radiology Reporting	32813	Resolved Radiology Reporting issue where loading thumbnail images returned an error during the dictation workflow.

SORTED BY CATEGORY AND REDMINE

VERSION DETAILS

Code Stream

The following source code branches have been merged into this release:

