



## eRAD RIS

# HOTFIX RELEASE ANNOUNCEMENT

### Hotfix 4.2024.022.1.6

UPDATED JULY 11, 2024

#### TABLE OF CONTENTS

**Summary** ..... 2

    Release Announcement .....2

    Resolved Issues .....2

    New Known Limitations .....2

**Version Details** ..... 3

    Code Stream .....3

#### PUBLICATION HISTORY

Revision	Author	Description
July 11, 2024	Kevin Brooks, Michelle Mahabir / Hilary Saltmarsh	<ul style="list-style-type: none"> <li>▪ Commercial release.</li> </ul>

# SUMMARY

## Release Announcement

This hotfix release of ERAD RIS 4.2024.022.1.6 resolves a PDF-to-image conversion issue.

## Resolved Issues

This release resolves the following issues:

Category	Redmine #	Subject
Scanning	35649	Resolved Scanning issue where documents received via Wedge were being truncated by PDF-to-image conversion.

## New Known Limitations


The following new Known Limitations were identified with this release:

- PDF FILES SENT VIA RADAR SECUREPIC MESSAGES OF MAY NOT FULLY DISPLAY.
  - **Issue:** RADAR SecurePIC is not properly displaying PDF files due to conversion problems.
  - **Impact:** For some attachments, the bottom part of the page displays as blacked out.
  - **Workaround:** No known workaround.

# VERSION DETAILS

## Code Stream

The following source code branches have been merged into this release:



**WARNING:** The functionality of this release has not yet been merged forward into any other release.

Please carefully review the code stream diagram of any future releases before deploying to ensure that the functionality of this release will not be lost when upgrading.