



eRAD RIS

RELEASE ANNOUNCEMENT

Build 4.2024.036

UPDATED MARCH 19, 2024

TABLE OF CONTENTS

Summary 2

 Release Announcement 2

 Resolved Issues 2

Version Details 3

 Code Stream 3

PUBLICATION HISTORY

Revision	Author	Description
March 18, 2024	Kevin Brooks, Michelle Mahabir / Hilary Saltmarsh	<ul style="list-style-type: none"> Commercial release.

SUMMARY

Release Announcement

This release of eRAD RIS 4.2024.036 resolves an HL7 EMR transmission issue.

Resolved Issues

This release resolves the following issues:

Category	Redmine #	Subject
HL7	32644	Resolved HL7 issue where an exception would occur upon the transmission of identical EMR orders.

VERSION DETAILS

Code Stream

The following source code branches have been merged into this release:

