



## eRAD RIS

# RELEASE ANNOUNCEMENT

## Rapid Release 4.2023.041.0.1

UPDATED APRIL 3, 2024

### TABLE OF CONTENTS

**Summary** ..... 2

    Release Announcement ..... 2

    Resolved Issues ..... 2

**Version Details** ..... 3

    Code Stream ..... 3

### PUBLICATION HISTORY

Revision	Author	Description
April 3, 2024	Kevin Brooks / Hilary Saltmarsh	<ul style="list-style-type: none"> <li>Commercial release.</li> </ul>

# SUMMARY

## Release Announcement

*This rapid release of ERAD RIS 4.2023.041.0.1 updates Provider Portal to resolve an issue with Single Sign-On (SSO) with Epic EMR systems.*

## Resolved Issues


This release resolves the following issues:

Category	Redmine #	Subject
35395	35395	Resolved character encoding issue when including Accession number in SSO links

# VERSION DETAILS

## Code Stream

The following source code branches have been merged into this release:





*WARNING: The functionality of this release has not yet been merged forward into any other release.*

*Please carefully review the code stream diagram of any future releases before deploying to ensure that the functionality of this release will not be lost when upgrading.*