



# eRAD RIS

## RELEASE ANNOUNCEMENT

### Build 4.2024.105

UPDATED OCTOBER 21, 2024

#### TABLE OF CONTENTS

**Summary** ..... 2

    Release Announcement..... 2

    Resolved Issues..... 2

**Version Details** ..... 3

    Code Stream..... 3

#### PUBLICATION HISTORY

Revision	Author	Description
October 21, 2024	Kevin Brooks, Michelle Mahabir / Hilary Saltmarsh	<ul style="list-style-type: none"> <li>Commercial release.</li> </ul>

# SUMMARY

## Release Announcement

*This release of ERAD RIS 4.2024.105 resolves an issue displaying Key Images in the Provider Portal via V8 and V9 PACS.*

## Resolved Issues

This release resolves the following issues:

Category	Redmine #	Subject
Provider Images	33549	Resolved Provider Images issue where Key Images were not displaying correctly from the V8 or V9 PACS.

SORTED BY CATEGORY AND REDMINE

# VERSION DETAILS

## Code Stream

The following source code branches have been merged into this release:

