

Server Update

for eRAD RIS

Version 1.0

Build 37

Update 37.008

Table of Contents

1. Purpose.....	3
Defects resolved.....	3
2. Package Contents	4
3. Intended Audience	4
4. Who is affected	4
5. Installing	4

RADNET, Inc. – eRAD RIS Server Update

1. Purpose

This server update has been produced to fix bugs required for production setting, or to complete an enhancement requested for production.

This server update includes the following:

Defects resolved

Redmine #	Subject	Priority
1857	New Appointment – Search for time slot takes 20 seconds at ARS	High

2. Package Contents

Files included for this server update:

1. RIS Service Update folder contains:
 - a. WCFSservice.dll
 - b. WCFSservice.pdb

3. Intended Audience

This document is created by the RIS Development team for the RadNet RIS Implementation team.

4. Who is affected

Previous versions of build 37.

5. Installing

- 1) Stop each service listed below. Copy the files from the RIS Service Updates folder to the directory for each rRISService.
rRISService directories include:
 - a) *rRISService*
 - b) *CDS rRISService*
 - c) *Wedge rRISService*
 - d) *DDE rRISService*

Restart each service