



eRAD RIS

RELEASE ANNOUNCEMENT

Build 4.2023.008

UPDATED MARCH 13, 2023

PUBLICATION HISTORY

Revision	Author	Description
March 13, 2023	Kevin Brooks / Hilary Saltmarsh	▪ Commercial release.

SUMMARY

Release Announcement

This release of eRAD RIS 4.2023.008 introduces a feature that adjusts output of generated C-CDA documents.

New Features

This release introduces the following features and enhancements:

Category	Redmine #	Subject	Description
Patient C-CDA	33643	Suppress sections in C-CDA documents	This change to Patient C-CDA modifies the output of generated C-CDA documents to suppress several sections within the document.

Refer to the Feature Details section below for configuration and usage information.

Resolved Issues

No additional issues were resolved with this release.

New Known Limitations

While no new known limitations were identified with this release, the following expected behaviors were noted:

- EXPECTED BEHAVIOR BY DESIGN (FROM #33643) - Until such a time that alternative data is decided to display in them, the Procedure Note, History, and Physical Narrative Note will be removed from the Clinical Notes section of C-CDA documents generated from the RIS Client or Portals. The following new Known Limitations were identified with this release:

FEATURE DETAILS

Patient C-CDA

Feature #33643 - Suppress sections in C-CDA documents

Summary

This change to Patient C-CDA modifies the output of generated C-CDA documents to suppress several sections within the document.

Background

Previously, as part of a CURES update the content of the C-CDA document template was updated. This included adding a new clinical notes section. It has been identified that some information should not display in this section.

Feature Description

Until such a time that alternative data is decided to display in them, the content and headings for `Procedure Note`, `History`, and `Physical Narrative Note` will be removed from the `Clinical Notes` section of C-CDA documents generated from the RIS Client or Portals.

Configuration Instructions

No System Administrator actions are necessary to enable this feature.