



eRAD RIS

RELEASE ANNOUNCEMENT

Build 4.2024.113

UPDATED NOVEMBER 7, 2024

TABLE OF CONTENTS

- Summary** 2
 - Release Announcement..... 2
 - New Features..... 2
 - Resolved Issues..... 2
- Feature Details** 3
 - Installer..... 3
 - Feature #36048 - Installer update to correct file permission changes..... 3
- Version Details** 4
 - Code Stream..... 4

PUBLICATION HISTORY

Revision	Author	Description
November 7, 2024	Kevin Brooks, Michelle Mahabir / Hilary Saltmarsh	<ul style="list-style-type: none"> ▪ Commercial release.

SUMMARY

Release Announcement

This release of ERAD RIS 4.2024.113 resolves a document cropping issue within the scanning process and merges the content of eRAD RIS 4.2024.111.0.1 (which resolves a Mobile Mammography Patient Appointments page redirection issue).

This release also resolves the previously released ERAD RIS SERVICE BULLETIN - INSTALLER FILE PERMISSION CHANGES.

New Features

This release introduces the following features and enhancements:

Category	Redmine #	Subject	Description
Installer	36048	Installer update to correct file permission changes	Resolves the previously released ERAD RIS SERVICE BULLETIN - INSTALLER FILE PERMISSION CHANGES.

Resolved Issues

This release resolves the following issues:

Category	Redmine #	Subject
Scanning	35745	Resolved Scanning issue which caused documents to be incorrectly cropped during the scanning process.

FEATURE DETAILS

Installer

Feature #36048 - Installer update to correct file permission changes

Summary

This enhancement to the RIS backend installer addresses the permissions issue identified in the ERAD RIS SERVICE BULLETIN - INSTALLER FILE PERMISSION CHANGES. A helper script now recursively updates files to remove “Blocked” attribute.

VERSION DETAILS

Code Stream

The following source code branches have been merged into this release:

