



eRAD RIS

RELEASE ANNOUNCEMENT

Build 4.2024.118

UPDATED NOVEMBER 20, 2024

TABLE OF CONTENTS

Summary 2

 Release Announcement2

 Resolved Issues2

Version Details 3

 Code Stream3

PUBLICATION HISTORY

Revision	Author	Description
November 20, 2024	Kevin Brooks, Michelle Mahabir / Hilary Saltmarsh	<ul style="list-style-type: none"> ▪ Commercial release.

SUMMARY

Release Announcement

This release of ERAD RIS 4.2024.118 resolves CD Burner issue and merges the content of ERAD RIS 4.2024.114.1 (which resolved Billing related issues within the IVT workflow.).

Resolved Issues

This release resolves the following issues:

Category	Redmine #	Subject
CD Burn	34033	Resolved CD Burner issue which caused the Epson Disc Producer to fail if the RequestID from PACS and eRAD RIS was longer than 40 characters.

VERSION DETAILS

Code Stream

The following source code branches have been merged into this release:

