



eRAD RIS

RELEASE ANNOUNCEMENT

Build 4.2024.122

UPDATED DECEMBER 13, 2024

TABLE OF CONTENTS

Summary 2

 Release Announcement2

 Resolved Issues2

 New Known Limitations2

Version Details 3

 Code Stream3

PUBLICATION HISTORY

Revision	Author	Description
December 13, 2024	Kevin Brooks, Michelle Mahabir / Hilary Saltmarsh	<ul style="list-style-type: none"> Commercial release.
December 14, 2024	Kevin Brooks, Michelle Mahabir / Hilary Saltmarsh	<ul style="list-style-type: none"> Corrected merge versions.

SUMMARY

Release Announcement

This release of *ERAD RIS 4.2024.122* resolves an IVT/Precert Wedge message issue, merges the content of hotfix *ERAD RIS 4.2024.083.1* (which resolves a Provider Portal issue where the button to switch viewer preferences was not visible) and merges the content of Rapid Release *ERAD RIS 4.2024.118.0.1* (which resolves two scheduling issues within the Patient Portal).

Resolved Issues

This release resolves the following issues:

Category	Redmine #	Subject
IVT/Precert	36149	Resolved IVT/Precert issue where necessary columns were omitted from Assigned IVT Wedge message.

New Known Limitations

The following new Known Limitations were identified with this release:

- EXPECTED BEHAVIOR BY DESIGN – UPGRADE SCRIPT ISSUE AFFECTING V083.1 SITES
 - **Issue:** The upgrade script fails when upgrading from v4.2024.083.1 to v4.2024.122 due to a missing db upgrade path.
 - **Impact:** Only sites upgrading from v4.2024.083.1 are impacted.
 - **Workaround:** There are two options for sites upgrading from v4.2024.083.1:
 1. Defer upgrading to a post-v4.2024.122 release which corrects this script issue, and no action will be required before upgrading.
 2. Utilize the following database upgrade script before starting an upgrade from v4.2024.083.1 to remove one line from the `s_version` table:


```
DELETE FROM s_version WHERE label = '4.24.083.1.004683760'
```

VERSION DETAILS

Code Stream

The following source code branches have been merged into this release:

