



eRAD RIS

RELEASE ANNOUNCEMENT

Build 4.2023.014

UPDATED FEBRUARY 5, 2024

TABLE OF CONTENTS

Summary 2

 Release Announcement2

 New Features2

Feature Details 3

 Lookup Tables3

 Feature #33745 – Add "Allow Routine Mammo Reminders" column to Gender LU table3

 System Config3

 Feature #33746 –Create Manual Reminders by Body Part3

Version Details 5

 Code Stream5

PUBLICATION HISTORY

Revision	Author	Description
February 5, 2024	Kevin Brooks, Michelle Mahabir / Hilary Saltmarsh	<ul style="list-style-type: none"> Commercial release.

SUMMARY

Release Announcement

This release of eRAD RIS 4.2023.014 introduces two features that assist with Lung Cancer Tracking.

New Features

This release introduces the following features and enhancements:

Category	Redmine #	Subject	Description
Lookup Tables	33745	Add "Allow Routine Mammo Reminders" column to Gender Lookup Table	This enhancement to Lookup Tables adds additional parameters to support configuration of the transmission of routine reminders by Gender Code.
System Config	33746	Create Manual Reminders by Body Part	This enhancement to System Config introduces a new configuration setting that allows the creation of manual reminders by body part.

SORTED BY CATEGORY AND REDMINE

Refer to the FEATURE DETAILS section below for configuration and usage information.

FEATURE DETAILS

Lookup Tables

Feature #33745 – Add "Allow Routine Mammo Reminders" column to Gender LU table

Summary

This enhancement to Lookup Tables adds additional parameters to support configuration of the transmission of routine reminders by Gender Code.

Background

Previously, routine reminders were restricted to `F` or `Female` gender code. However, the introduction of LDCT/Lung-RADS tracking led to the necessity of configuring reminders using additional gender codes to support the use of the routine reminder functionality.

Feature Description

An `Allow Routine Mammo Reminders` column has been added to the `Gender` lookup table. When set to `Y`, routine reminders will be created for configured Gender Code.

Configuration Instructions

No System Administrator actions are necessary to enable this feature; however, optional configuration is available.

RIS Client

Changes to Gender Lookup Table Settings

- Optionally, set `Allow Routine Mammo Reminders` value to `Y` for desired Gender Code
 - Default value is `Y` for all Gender Codes with the exception of `MALE` or `M`

System Config

Feature #33746 – Create Manual Reminders by Body Part

Summary

This enhancement to System Config introduces a new configuration setting that allows the creation of manual reminders by body part.

Background

The introduction of LDCT/Lung-RADS tracking, requires the creation of manual reminders by body part. However, the limitation of the current system config, which accepted only one value, hindered the creation of manual reminders for other body parts.

Feature Description

Manual reminders are currently created based on the configured values retrieved from the `BodyPart` lookup table. A new `ReminderBodyPartCode` system configuration has been added that will allow the creation and transmission of manual reminders for multiple body parts.

The `Manual Reminder` option will appear on the right click context menu for all configured `Body Part` codes.

Configuration Instructions

No System Administrator actions are necessary to enable this feature.

RIS Client

Changes to SystemConfig Lookup Table

- Optionally, adjust `ReminderBodyPartCode` with desired body part codes for your installation.
 - Values must be comma separated

The following related settings were added or updated:

Setting	Default	Purpose
ReminderBodyPartCode	Blank	Comma separated list of Body Part code values utilized to initiate creation of reminders. Updated by #33746

VERSION DETAILS

Code Stream

The following source code branches have been merged into this release:

