



eRAD RIS

RELEASE ANNOUNCEMENT

Build 4.2025.033

UPDATED JUNE 5, 2025

TABLE OF CONTENTS

Summary	2
Release Announcement	2
Resolved Issues	2
Version Details	3
Code Stream	3

PUBLICATION HISTORY

Revision	Author	Description
June 5, 2025	Michelle Mahabir/ Kevin Brooks	▪ Commercial release.

SUMMARY

Release Announcement

This release of eRAD RIS 4.2025.033 resolves numerous worklists issues, many relating to PreCert Status and codes, and a Contact Log issue.

Resolved Issues

This release resolves the following issues:

Category	Redmine #	Subject
IVT/Precert	36331	Resolved IVT/Precert issue where the Precert status was incorrectly set to Required if the expiry date was in the past.
IVT/Precert	37113	Resolved IVT/Precert issue where clicking the IVT window returns an error when no row is selected.
Scheduling	36842	Resolved an issue where the "Show All" checkbox in the Contact Log was not functioning as intended, preventing users from viewing all contact log entries across different views.
Worklists	36507	Resolved Worklists issue that prevented the Precert Required column on the OTS and IVT worklist from being filtered correctly.
Worklists	36531	Resolved Worklists issue where patient phone numbers were not displayed on the OTS worklist unless marked as primary in the Patient tab.
Worklists	36880	Resolved Worklists issue where deactivated Precert codes appeared as blank on the worklists, making it difficult for users to update or clear the selection.

SORTED BY CATEGORY AND REDMINE

VERSION DETAILS

Code Stream

The following source code branches have been merged into this release:

